

Jan 29, 2025

Imaging axonal calcium dynamics in *ex vivo* mouse brain slices

 In 1 collection

DOI

dx.doi.org/10.17504/protocols.io.81wgbx9eolpk/v1



Yan-Feng Zhang^{1,2}, Stephanie J Cragg^{3,4}

¹Department of Clinical and Biomedical Sciences, University of Exeter, Exeter, United Kingdom;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, USA;

³Department of Physiology, Anatomy and Genetics, University of Oxford, Oxford OX1 3PT, UKAligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, USA;

⁴Oxford Parkinson's Disease Centre, University of Oxford, Oxford, United Kingdom

Team Cragg



Cláudia C. Mendes

University of Oxford

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.81wgbx9eolpk/v1

Protocol Citation: Yan-Feng Zhang, Stephanie J Cragg 2025. Imaging axonal calcium dynamics in *ex vivo* mouse brain slices. [protocols.io https://dx.doi.org/10.17504/protocols.io.81wgbx9eolpk/v1](https://dx.doi.org/10.17504/protocols.io.81wgbx9eolpk/v1)

Manuscript citation:

Zhang et al. (2024) [An axonal brake on striatal dopamine output by cholinergic interneurons](#), bioRxiv

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: November 27, 2023

Last Modified: January 29, 2025



Protocol Integer ID: 91481

Keywords: calcium, ex vivo, calcium, mouse, dopamine, striatum

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-020370

Abstract

This protocol describes how to image calcium dynamics in striatal dopaminergic axons in *ex vivo* mouse brain slices. We imaged calcium transients in response to single and trains (4 pulses, 100 Hz) of electrical stimulus pulses using genetically encoded calcium indicator GCaMP6f expressed in DAT-Cre:Ai95D mice.

Materials

Equipment:

- Olympus BX51WI microscope equipped with a OptoLED Lite system (CAIRN Research);
- Prime Scientific CMOS (sCMOS) Camera (Teledyne Photometrics);
- x40/0.8 NA water-objective (Olympus UK)
- ITC-18 A/D board (Instrutech)

Software:

- Micro-Manager v1.4
- Matlab vR2019b
- Fiji v1.5
- Igor Pro 6 (WaveMetrics)

Before start

This protocol was performed in heterozygous DAT-Cre:Ai95D (4-7 weeks) mice. These mice were bred from homozygous DAT-Cre mice (B6.SJL-Slc6^{a3tm1.1(cre)Bkmn}/J, JAX stock number 006660) crossed with homozygous Ai95D mice (B6;129S-Gt(ROSA)26Sor^{tm95.1(CAG-GCaMP6f)Hze}/J, JAX stock number 028865).

Note

Ai95(RCL-GCaMP6f)-D (also called Ai95D; RRID:IMSR_JAX:028865) mice are a Cre-dependent, fluorescent, calcium-indicator tool strain. Ai95D has a floxed-STOP cassette preventing transcription of the GCaMP6 fast variant calcium indicator (GCaMP6f; a detector of single neuronal action potentials with fast response kinetics). After Cre exposure, bright EGFP fluorescence is observed following calcium binding (such as neuronal activation).

This Ai95D allele is on a C57BL/6J genetic background.

We prepare *ex vivo* mouse brain slices by performing **steps 1 to 11** from **Protocol: Fast-scan cyclic voltammetry to assess dopamine release in ex vivo mouse brain slices**.

Image Acquisition

- 1 Using a x40/0.8 NA water-objective (Olympus UK), position the stimulating electrode on the surface of the brain slice and centre it in the field of view.
- 2 Change the exposure time to reach a frame rate of around 16.6 Hz every 2.5 min using Micro-Manager v1.4. 16.6 Hz frame rate every 2.5 min using Micro-Manager 1.4.
- 3 Apply electrical stimulus pulses singly and in trains (4 pulses, 100 Hz) using custom-written procedures in Igor Pro 6 (WaveMetrics) and an ITC-18 A/D board (Instrutech).

Note

The order of single and train stimulations was alternated and equally distributed and data were collected in duplicate before and after a change in extracellular experimental condition.

- 4 Record changes in fluorescence intensity using custom-written procedures in Igor Pro 6 (WaveMetrics) and an ITC-18 A/D board (Instrutech).

Image Analysis

- 5 **The following steps were performed in MATLAB vR2019b and Fiji v1.5.**

Extract fluorescence intensity from the region of interest $25\ \mu\text{m} \times 25\ \mu\text{m}$ which was $50\ \mu\text{m}$ away from the electrical stimulating electrode tip.

- 6 After background subtraction, bleach-correct the Ca^{2+} transients by fitting an exponential curve function through both the baseline (2 s prior to stimulation) and the last 1 s in a 7.2 s recording window.
- 7 Expressed data as $\Delta F/F$ where F is the fitted curve.