



Sep 03, 2020

SUPER rapid kit

Hyeon Jin Kim¹, Seong T. Hong²¹SNJ Pharma Inc; ²Jeonbuk National University

1

Works for me

dx.doi.org/10.17504/protocols.io.bkqvkvw6

Hyeon Kim

ABSTRACT

We developed a supersensitive capillary reaction kit intended for the direct and qualitative detection of SARS-CoV-2 virus in saliva and nasal swabs. It enables sensitive and specific detection of the SARS-CoV-2 spike antigens in respiratory specimens more accurately than most PCR method within 5 minutes. Our rapid capillary kit will help not only fight but also end COVID-19 by providing the first sensitive self-diagnosis at home as well as POC clinical diagnosis at the health care facilities.

DOI

dx.doi.org/10.17504/protocols.io.bkqvkvw6

PROTOCOL CITATION

Hyeon Jin Kim, Seong T. Hong 2020. SUPER rapid kit. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bkqvkvw6>



LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Sep 03, 2020

LAST MODIFIED

Sep 03, 2020

PROTOCOL INTEGER ID

41461

- 1 Prepare the specimen by collecting saliva, sputum, or nasal swab for testing
- 2 Dispense 2~5 drops of the specimen into the sampling tube (provided) filled with a general transfer pipette
 - 2.1 To prepare a specimen solution with a nasal swab, insert the testing swab into the sampling tube, roll along the side of the tube 5 times, and discard the used swab to biohazard waste.
- 3 Mix the content by pipetting up & down gently

- 4 Dispense 5 drops of mixed samples into the sample well of the cassette (provided).
- 5 Observe the appearance of the blue line for 1 min and record the result.