



2 ▼

Feb 15, 2022

# First Strand Synthesis with Reverse Transcriptase (M0253) V.2

[New England Biolabs<sup>1</sup>](#)<sup>1</sup>New England Biolabs

1

[dx.doi.org/10.17504/protocols.io.bddyi27w](https://dx.doi.org/10.17504/protocols.io.bddyi27w)**New England Biolabs (NEB)**Tech. support phone: **+1(800)632-7799** email: **info@neb.com****New England Biolabs**  
New England Biolabs

This protocol is for First Strand Synthesis with M-MuLV Reverse Transcriptase (M0253).

DOI

[dx.doi.org/10.17504/protocols.io.bddyi27w](https://dx.doi.org/10.17504/protocols.io.bddyi27w)<https://www.neb.com/protocols/0001/01/01/first-strand-synthesis-protocol-with-reverse-transcriptase>

New England Biolabs 2022. First Strand Synthesis with Reverse Transcriptase (M0253). **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.bddyi27w>[New England Biolabs](#)

M-MuLV reverse transcriptase, first strand reverse transcription, transcription

\_\_\_\_\_ protocol ,

Mar 08, 2020

Feb 15, 2022

33944

## MATERIALS

 **Magnesium Chloride Fisher**

**Scientific Catalog #AC223210010**

 **DTT (Dithiothreitol) (> 99% pure) Protease free Gold**

**Biotechnology Catalog #DTT**

 **Tris-HCl Life**

**Technologies Catalog #AM9855**

 **Potassium chloride Sigma**

**Aldrich Catalog #P9333**

 **M-MuLV Reverse Transcriptase New England**

**Biolabs Catalog #M0253**

### Materials needed:

#### 10X RT buffer:

A	B
Tris-HCl (pH 8.3 @ 25°C)	500 mM
KCl	750 mM
MgCl <sub>2</sub>	30 mM
DTT	100 mM

**dNTP mix** (2.5 mM each in water titrated by Tris-HCl to pH 7.0)

**Oligo-dT primer** (40 µM)

**Random nonamers** (40 µM)

**RNase Inhibitor** (10 U/µL)

**M-MuLV Reverse Transcriptase** (200 units/µL)

Please refer to the Safety Data Sheets (SDS) for health and environmental hazards.

Prepare the following solutions:

#### 10X RT buffer:

A	B
Tris-HCl (pH 8.3 @ 25°C)	500 mM
KCl	750 mM
MgCl <sub>2</sub>	30 mM
DTT	100 mM

**dNTP mix** (2.5 mM each in water titrated by Tris-HCl to pH 7.0)

1 In a sterile microfuge tube add the following:

Reagent	Volume
RNA solution	0.5-2 µg (total RNA) or 50-100 ng (PolyA-selected RNA)
Primer (Oligo-dT primer at 40 µM or Random nonamer N9 at 40 µM)	2 µL
dNTP mix (2.5 mM each in water titrated by Tris-HCl to pH 7.0)	4 µL
nuclease-free H <sub>2</sub> O	to final volume of 16 µL

2 Heat for 3-5 minutes at **65 °C** - **80 °C** .

3 Spin briefly and place promptly **On ice** .

4



Add:

**2 µL 10X RT buffer**

**1 µL RNase inhibitor**

**1 µL M-MuLV Reverse Transcriptase**

Final volume: **20 µL**

5



Incubate at **42 °C** for **01:00:00** .

6 Inactivate enzyme at **90 °C** for **00:10:00** .

7 Store products at **-20 °C** or proceed to next step(s).

