



DEC 27, 2022

## OPEN ACCESS

### DOI:

[dx.doi.org/10.17504/protocols.io.e6nvwjjewlmk/v1](https://doi.org/10.17504/protocols.io.e6nvwjjewlmk/v1)

**Protocol Citation:** Sarah M Prostack, Edgar M Medina, Erik Kalinka, Lillian Fritz-Laylin 2022. Protocol 7: Picking colonies of transformed *Spizellomyces punctatus* (Sp). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.e6nvwjjewlmk/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**Created:** Nov 03, 2022

**Last Modified:** Dec 27, 2022

**PROTOCOL integer ID:** 72238

# Protocol 7: Picking colonies of transformed *Spizellomyces punctatus* (Sp)

In 1 collection

Sarah M Prostack<sup>1</sup>, Edgar M Medina<sup>1</sup>, Erik Kalinka<sup>1</sup>, Lillian Fritz-Laylin<sup>1</sup>

<sup>1</sup>University of Massachusetts at Amherst



fritzlaylinlab.umass

## ABSTRACT





If transformation was successful, antibiotic-resistant *Spizellomyces* colonies should appear after 4-6 days of growth on selection media. These colonies should be rough, white, and opaque. When viewed under a microscope, there should be clustering of sporangia and some movement of zoospores. Picking these colonies and culturing them until there are enough cells for cryopreservation and to use in downstream molecular analysis is the last protocol in the entire transformation process. All steps here, with the exception of centrifugation steps, must be carried out in a sterile environment, either in a laminar flow hood or in the sterile area around an open flame.

## ATTACHMENTS

[Spizellomyces transformation steps.pdf](#)

## MATERIALS TEXT

### Materials

- Active plate of colonies of Sp transformants (see *Protocol 6: Selecting for Spizellomyces punctatus* transformants)
- DS Solution, sterile (see recipe)
- 18G needle, sterile  18G needle **BD Biosciences Catalog #305196**
- 1.5 mL microcentrifuge tubes, sterile such as  Fisherbrand™ Low-Retention Microcentrifuge Tubes **Fisher Scientific Catalog #0681-331**
- K1 agar plates (1.5% w/v) with selection antimicrobials, sterile (see recipe)
- 20-200 µL micropipette such as  Eppendorf Research Plus single channel pipette 2-20 uL yellow operating button for use with 20 uL **pipette.com Catalog #3123000039**
- Filter tips for the micropipettes, sterile such as  TIPONE® FILTER TIPS **USA Scientific Catalog #1122-1830**
- Laminar flow hood and/or open flame, for maintaining sterility
- 70% (v/v) ethanol for maintaining sterility (if using laminar flow hood)

## SAFETY WARNINGS



Proper handling of 18G needles must be observed to avoid personal injury. All sharps should be disposed of into a proper sharps container.

## Steps

- 1 Divide one K1 plate with selection antimicrobials into four sections using a marker on the bottom of the plate.

### Note

One plate can be used to plate up to 4 colonies for one plasmid transformation. This reduces materials used in the initial amplification stage.

- 2 Aliquot  50 µL of DS into one microcentrifuge tube per colony to be picked.

- 3 Using an 18G needle, gently lift the colony of interest from the agar.



### Note

- Be careful not to poke too deeply into the agar. The goal is to lift as little agar as possible from the plate, while picking up most of the colony.
- To prevent cross contamination between colonies, either flame the needle between picking colonies, or use a new, sterile needle for each colony.


- 4 Resuspend the colony into the appropriate tube pre-filled with DS.

### Note

Swirl the needle gently, but with enough force to dislodge the colony from the needle.


- 5 Pipette gently to mix and break up the pellet.



6 Pipette  50  $\mu\text{L}$  of the resuspended colony onto one quadrant of the K1 plate prepared in step 1.



7 Repeat steps 3-6 for each colony to be picked.

8 Let the cells grow for 2-3 days before rehydrating the quadrants with  100  $\mu\text{L}$  of DS and transferring each colony to its own quadrant of a new K1 selection plate.

#### Note

Only a small portion of the quadrant is needed for colony expansion.

9 After 2-4 rounds of colony expansion there is enough inoculum to transfer each isolated colony into its own plate.

10 Continuing subculturing each colony until enough culture exists to freeze and perform downstream analyses.