



Mar 23,
2020

STRIPE-seq library construction v.3

Robert Policastro¹, Gabe Zentner¹

¹Indiana University

1 Works for me [dx.doi.org/10.17504/protocols.io.bdtri6m6](https://doi.org/10.17504/protocols.io.bdtri6m6)



Gabe Zentner
Indiana University



ABSTRACT

Accurate mapping of transcription start sites (TSSs) is key for understanding transcriptional regulation; however, current protocols for genome-wide TSS profiling are laborious and expensive. We present Survey of TRanscription Initiation at Promoter Elements with high-throughput sequencing (STRIPE-seq), a simple, rapid, and cost-effective protocol for sequencing capped RNA 5' ends from as little as 50 ng total RNA. Including depletion of uncapped RNA and bead cleanups, a STRIPE-seq library can be constructed in approximately 5 hours.

MATERIALS

NAME	CATALOG #	VENDOR
Terminator 5-Phosphate-Dependent Exonuclease	TER51020	Lucigen
RNAClean XP	A63987	Beckman Coulter
5M Betain	AAJ77507UCR	Thermo Fisher Scientific
KAPA HiFi HotStart ReadyMix	KK2601	Roche
Sorbitol	DSS23080-500	Dot Scientific
Trehalose	0210309705	MP Biomedicals
dNTPs 10 µM each	97063-232	VWR Scientific
SuperScript II Reverse Transcriptase	18064014	Thermo Fisher Scientific
RNA ScreenTape	5067-5576	Agilent Technologies
High Sensitivity D5000 ScreenTape	5067-5592	Agilent Technologies

BEFORE STARTING

Prepare 3.3 M sorbitol/0.66 M trehalose solution as per Batut and Gingeras (PMID 24510412).

1. Add **2 ml RNase-free H₂O** to a 50 mL tube.
2. Add **8.02 g trehalose** to the tube.
3. Add **3 ml RNase-free H₂O**.
4. Add **17.8 g sorbitol** to the tube.
5. Add **5.5 ml RNase-free H₂O**
6. Bring volume to 30 mL with **0 ml RNase-free H₂O**
7. Transfer to an RNase-free glass bottle and autoclave at 121 °C for 30 min.

Store **1.5 ml** aliquots at **Room temperature** protected from light.

Prepare Total RNA

- 1 Check RNA quality and concentration on an Agilent TapeStation using a High-Sensitivity RNA ScreenTape.

15m



You should have at least 50 to 200 ng of total RNA at a concentration of at least 30 to 125 ng/μl. Your total RNA should also not be highly degraded, as measured by the quality of the rRNA peaks.



TapeStation

Agilent G2991AA [Link](#)

Terminator Exonuclease (TEX) Digestion of Uncapped RNA

- 2 **Prepare TEX Reaction.** TEX preferentially degrades uncapped RNA, thus reducing the amount of rRNA and degraded mRNA fragments in the sample.



TEX is magnesium-dependent, so ensure that the RNA storage buffer does not contain EDTA.

- 2.1 Create TEX master mix. Prepare a sufficient volume for the number of reactions to be performed + 1 to account for volume^{3m} loss during pipetting.
 1. **0.2 μl Terminator Exonuclease** .
 2. **0.2 μl Terminator Exonuclease Reaction Buffer A** .Vortex to mix and spin down.
- 2.2 Prepare TEX reactions in 0.2 mL PCR tubes.
 1. **0.4 μl TEX Master Mix**
 2. Up to **1.6 μl Total RNA** .
 3. Nuclease free water to **2 μl** total reaction volume.Vortex to mix and spin down.

- 3 Incubate the TEX reactions in thermal cycler. 1h
1. **30 °C** for **01:00:00** .
 2. **4 °C Hold** .





This is a good time to prepare the Reverse Transcription Oligo (RTO) annealing and Template Switching Reverse Transcription (TSRT) reaction mixtures from steps 4.1 and 5.1.

Template Switching Reverse Transcription

- 4 **Anneal reverse transcription oligo (RTO) to RNA.** STRIPE-seq primes reverse transcription via a random pentamer adjacent to the full length TrueSeq R2 adapter (including the barcode) in the RTO.
- 4.1 Prepare one RTO annealing mix per sample in 0.2 mL PCR tubes. 5m
1. **1.5 µl Sorbitol/Trehalose Solution** .
 2. **1 µl Reverse Transcription Oligo (RTO)** [] **10 Micromolar (µM)** . Each sample should have its own unique barcode.
 3. **0.5 µl dNTPs** [] **10 Millimolar (mM)** Each .
- Vortex to mix and spin down.
- 4.2 Add **2 µl TEX Reaction** (from step 3) to **3 µl RTO Annealing Mixture** (from step 4.1). Vortex to mix and spin down. 3m
- 4.3 Incubate RTO annealing mixture in thermal cycler. 7m
1. **65 °C** **00:05:00** .
 2. **4 °C** **00:02:00** .
 3. **4 °C Hold** .
- 5 **Prepare template switching reverse transcription (TSRT) reactions.** The process of TSRT enriches for the 5' ends of capped RNA in the final library.
- 5.1 Prepare TSRT reaction master mix (per sample). 5m
1. **2 µl Betaine** [] **5 Molarity (M)** .
 2. **2 µl 5X SuperScript II First Strand Buffer** .
 3. **0.5 µl DTT** [] **0.1 Molarity (M)** .
 4. **0.5 µl SuperScript II Reverse Transcriptase** .
- Vortex to mix and spin down.



Add reverse transcriptase to master mix just prior to adding to samples.


- 5.2 Add  **5 µl TSRT Master Mix** (from step 5.1) into the  **5 µl RTO Annealing Reaction** from step 4.3. Vortex to mix and spin down. 3m

6 TSRT.

- 6.1 First half of TSRT reaction. 25m
1.  **25 °C**  **00:10:00** .
 2.  **42 °C**  **00:05:00** .



Move on to step 6.2 immediately after the end of step 6.1.

- 6.2 Add TSO. Keep the samples in the thermal cycler while adding the TSO. 3m
1.  **0.25 µl TSO** [**M**] **400 Micromolar (µM)** .
 2. Quickly vortex to mix, spin down, and immediately place tubes back in thermal cycler.



Move on to step 6.3 immediately after end of step 6.2.

- 6.3 Second half of TSRT reaction. 30m
1.  **00:25:00**  **42 °C** .
 2.  **00:10:00**  **70 °C** .
 3.  **4 °C Hold** .



This is a good time to prepare the library PCR master mix in step 8.1.

- 7 Cleanup of TSRT product. 20m
1. Transfer the TSRT product from step 6.3 into 0.5 mL tube.
 2. Pipette **8 µl RNAClean XP Beads** up and down 10 times into **10 µl TSRT Reaction** from step 6.3.
 3. Incubate for **00:05:00** at **Room temperature**.
 4. Place tubes on magnetic rack and incubate for **00:05:00** at **Room temperature**.
 5. Carefully aspirate supernatant, leaving ~ **2 µl** in tube to avoid sucking up beads.
 6. While tube is still on rack, wash beads with **175 µl 70% Ethanol**, and immediately discard wash without incubation.
 7. Air dry beads for **00:05:00** at **Room temperature**.
 8. Resuspend beads in **12 µl Nuclease Free Water**, and incubate on magnetic rack for **00:01:00** at **Room temperature**.
 9. Transfer **11 µl Supernatant** into new 0.2 mL PCR tubes.

Library PCR

8 Prepare library PCR reaction.

- 8.1 Create library PCR master mix (per sample). 5m
1. **12.5 µl 2X KAPA HiFi HotStart ReadyMix**.
 2. **0.75 µl Forward Library Oligo (FLO)** **10 Micromolar (µM)**.
 3. **0.75 µl Reverse Library Oligo (RLO)** **10 Micromolar (µM)**.
- Vortex to mix and spin down.

- 8.2 Add **14 µl Library PCR Master Mix** (from step 8.1) into **11 µl Cleaned TSRT Product** (from step 7). Vortex to mix and spin down. 2m

9 Run library PCR reaction. 45m

Initial Denaturation:

- **95 °C 00:03:00**

16-20 cycles:

- **98 °C 00:00:20**
- **63 °C 00:00:15**
- **72 °C 00:00:45**

Final Extension:

- **72 °C 00:02:00**
- **4 °C Hold**

- 10 **Size selection of final library.** SPRI bead size selection is used to remove fragments that are outside the ideal size for Illumina sequencing.

20m

- 10.1 Removal of small fragments.
1. Transfer library PCR product from step 9 into 0.5 mL tube.
 2. Pipette **16.3 µl RNAClean XP Beads** up and down 10 times into **25 µl Library PCR Product** from step 9.
 3. Incubate for **00:05:00** at **Room temperature**.
 4. Place tubes on magnetic rack and incubate for **00:05:00** at **Room temperature**.
 5. Carefully aspirate supernatant, leaving ~ **2 µl** in tube to avoid sucking up beads.
 6. While tube is still on rack, wash beads with **175 µl 70% Ethanol** and immediately discard wash without incubation.
 7. Air dry beads for **00:05:00** at **Room temperature**.
 8. Resuspend beads in **17 µl Nuclease Free Water** and incubate on magnetic rack for **00:01:00** at **Room temperature**.
 9. Transfer **15 µl Supernatant** to new 0.5 mL tube.
 10. **Optional:** Reserve **1 µl Remaining Supernatant** from beads if you would like to see library size distribution after removing small fragments.

40m

- 10.2 Removal of large fragments.
1. Pipette **8.3 µl RNAClean XP Beads** up and down 10 times into **15 µl Cleaned Product** from step 10.1. Make sure to vortex the beads again prior to use.
 2. Incubate for **00:10:00** at **Room temperature**.
 3. Place tubes on magnetic rack and incubate for **00:10:00** at **Room temperature**.
 4. Transfer **22 µl Supernatant** to new tube.
 5. Pipette **22 µl RNAClean XP Beads** up and down 10 times into **22 µl Supernatant** from previous step.
 6. Incubate for **00:05:00** at **Room temperature**.
 7. Place tubes on magnetic rack and incubate for **00:05:00** at **Room temperature**.
 8. Carefully aspirate supernatant, leaving ~ **2 µl** in tube to avoid sucking up beads.
 9. While tube is still on rack, wash beads with **175 µl 70% Ethanol**, and immediately discard wash without incubation.
 10. Air dry beads for **00:05:00** at **Room temperature**.
 11. Resuspend beads in **16 µl Nuclease Free Water**, and incubate on magnetic rack for **00:01:00** at **Room temperature**.
 12. Transfer **15 µl Supernatant** to new tube.

Library Quality Control

- 11 Run final libraries on the Agilent TapeStation using a High Sensitivity D5000 ScreenTape.

15m



Final libraries should be distributed between 250 to 750 bp with a total library amount of 25 to 100 ng.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited