

Sep 06, 2022

# In vitro synthesis of PE2 (nCas9-MMLV RT fusion) polyA mRNA using T7 RNA polymerase

Donald Rio<sup>1</sup>

<sup>1</sup>Dept. of Molecular and Cell Biology, University of California, Berkeley  
Donald Rio: don\_rio@berkeley.edu

1 Works for me

Share

[dx.doi.org/10.17504/protocols.io.b3fmqjk6](https://dx.doi.org/10.17504/protocols.io.b3fmqjk6)

Devin E Snyder

## ABSTRACT

We describe the preparation of synthetic capped and polyadenylated messenger RNA (mRNA) encoding a Cas9 nickase-reverse transcriptase fusion protein (PE2) for use in genome editing with a method called prime editing.

## ATTACHMENTS

[336-739.docx](#)

## DOI

[dx.doi.org/10.17504/protocols.io.b3fmqjk6](https://dx.doi.org/10.17504/protocols.io.b3fmqjk6)

## PROTOCOL CITATION

Donald Rio 2022. In vitro synthesis of PE2 (nCas9-MMLV RT fusion) polyA mRNA using T7 RNA polymerase. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.b3fmqjk6>



## FUNDERS ACKNOWLEDGEMENT

Aligning Science Across Parkinson's  
Grant ID: ASAP-000486

## KEYWORDS

Linearization of plasmid, In vitro transcription, Purification of IVT product, polyA tailing, ASAPCRN

## LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited


CREATED

Jan 04, 2022

LAST MODIFIED

Sep 06, 2022

OWNERSHIP HISTORY

Jan 04, 2022  renuka.s

Jan 10, 2022  Devin E Snyder

PROTOCOL INTEGER ID

56525

MATERIALS TEXT

### Materials

 [HiScribe T7 ARCA mRNA Kit \(with Tailing\) - 20 rxns New England](#)

**Biolabs Catalog #E2060S**

 [Monarch RNA Cleanup Kit \(50 µg\) - 100 preps New England](#)

**Biolabs Catalog #T2040L**

 [PmeI - 2,500 units New England](#)

**Biolabs Catalog #R0560L**

 [DNase I \(RNase-free\) - 5,000 units New England](#)

**Biolabs Catalog #M0303L**

 [E.coli Poly \(A\) Polymerase - 500 units New England](#)

**Biolabs Catalog #M0276L**

- CutSmart buffer
- 5M NaCl
- 100% ethanol
- 10mM Tris-HCl
- 1mM EDTA
- TE buffer
- New England Biolabs Monarch RNA cleanup column
- Nanodrop1000c spectrophotometer

## Protocol Overview





- 1 Preparation of linearized plasmid DNA encoding the PE2 protein for use as a template for in vitro transcription by T7 RNA polymerase.


- 2 Synthesis of mRNA by T7 RNA polymerase, using a 5' cap nucleotide to initiate the mRNA and using E. coli poly A polymerase to add a synthetic poly A tail after transcription.
- 3 Column purification of the final synthetic, capped and polyadenylated mRNA.


#### Part 1: Linearization of pCMV-PE2 plasmid


4 




4h

Digest  **100 µg** plasmid with 200 units Pme I restriction endonuclease in a  **1 mL** reaction using CutSmart buffer for  **04:00:00** . at  **37 °C** .

5 

Purify the cleaved DNA by phenol-chloroform extraction, collect the upper (aqueous) phase and transfer to fresh  **1.5 mL** Eppendorf tubes.

6 

Add  **5 Molarity (M)** NaCl to a final concentration of  **0.1 Molarity (M)** (  **20 µL** ) and 3 volumes of 100% ethanol.

7 

4h

Leave at  **-80 °C** for  **02:00:00** to  **Overnight** to precipitate.


8 

10m

Spin  **16000 x g, 4°C, 00:10:00** .

9 

30m

Remove the ethanol and air dry for  **00:30:00** .

10



Resuspend at **500 undetermined** in TE (**10 millimolar (mM)** Tris-HCl, **7.5**, **1 millimolar (mM)** EDTA) buffer. Store the DNA sample at **-20 °C**.

## Part 2: In vitro transcription and polyA tailing

11

In order get a high yield of the 6500nt long PE2 mRNA by T7 RNA polymerase transcription, set up 8 x **20 µL** in vitro transcription reactions using **1 µg** of cleaved template DNA in each reaction and the New England Biolabs HiScribe T7 ARCA kit with tailing (E2060S), according to the manufacturer's instructions.

12



2h

Incubate at **37 °C** for **02:00:00** in an incubator (not a temp block).

13



For each **20 µL** reaction, add **2 µL** DNase I (stock is **2 U/µL**).

14



15m

Incubate at **37 °C** for **00:15:00**.

15



Dilute the reaction to **50 µL** with **20 µL** RNase-free water, **5 µL** 10X polyA polymerase buffer and **5 µL** E. coli polyA polymerase.

16





30m

Incubate at **37 °C** for **00:30:00**.

### Part 3: Purification of IVT product

5m

- 17 Purify the polyadenylated mRNA on  **50 µg** New England Biolabs Monarch RNA cleanup columns (2 of  **20 µL** reactions per column) (T2040L).

- 18 

Dilute each  **20 µL** reaction with  **100 µL** binding buffer and  **150 µL** 100% ethanol.

- 19 

1m

Add the two diluted reactions to one column and spin at  **16000 x g, 00:01:00**.

- 20 

Wash each column with wash buffer and spin twice (see below).

- 20.1  

1m

Wash with  **500 µL** wash buffer (spin at  **16000 x g, 00:01:00**). (1/2)




- 20.2  

1m

Wash with  **500 µL** wash buffer (spin at  **16000 x g, 00:01:00**). (2/2)

- 21 

1m

Elute in  **25 µL** of RNase-free H<sub>2</sub>O per column, incubate for  **00:01:00** at  **Room temperature**.

- 22 

1m

Spin at  **16000 x g, 00:01:00** into fresh  **1.5 mL** tubes and pool elution together into one tube.

23 Measure the RNA concentration using a Nanodrop1000c spectrophotometer.

The yield from the total of 8 reactions is usually ~ **200 µg**.

24 Store the mRNA at **-80 °C** and avoid repeated freezing and thawing.