



Aug 31, 2020

CDx autoXpress detection of COVID 19

Mickey Shah¹, Gary D Niehaus¹¹Crystal Diagnostics

In Development

dx.doi.org/10.17504/protocols.io.bki4kugw

Mickey Shah

ABSTRACT

The Crystal Diagnostics AutoXpress provides rapid detection of COVID 19 at 10 - 100 CFU per ul of clinical sample (nasal swab, sputum, or saliva). The patented Liquid Crystal technology is combined with antibody-coated paramagnetic microspheres to selectively capture and detect COVID 19.

DOI

dx.doi.org/10.17504/protocols.io.bki4kugw

PROTOCOL CITATION

Mickey Shah, Gary D Niehaus 2020. CDx autoXpress detection of COVID 19. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bki4kugw>

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 31, 2020

LAST MODIFIED

Aug 31, 2020



PROTOCOL INTEGER ID

41276

Automated

- 1 Transfer **600 µl** of test sample to 96 well plate. Each sample will be transferred to two wells, one for POS and one for NEG.
- 2 Mix the magnetic 3µm microspheres for POS and NEG using a pipette.
- 3 Dispense **100 µl** of the magnetic 3µm microspheres to each well.
- 4 Shake the plate at **1200 rpm, Room temperature , 00:00:10** , followed by **500 rpm, Room temperature , 00:15:00** .

- 5 Transfer the plate to the magnet block and wait for ⌚ 00:05:00
- 6 Aspirate all the liquid from each well without disturbing the pellet.
- 7 Add 🧴 900 µl of PBS to each well and wait for ⌚ 00:03:00
- 8 Aspirate all the liquid from each well without disturbing the pellet.
- 9 Transfer plate to shaker from the magnet block and shake the plate at 🌀 1000 rpm, Room temperature , 00:00:10
- 10 Mix the 10um poly microspheres using a pipette.
- 11 Dispense 🧴 100 µl of the 10um poly microspheres to each well.
- 12 Shake the plate at 🌀 900 rpm, Room temperature , 00:00:30 , followed by 🌀 500 rpm, Room temperature , 00:15:00
- 13 Transfer the plate to magnet block and wait for ⌚ 00:05:00
- 14 Aspirate all the liquid from each well without disturbing the pellet.
- 15 Dispense 🧴 300 µl of PBS to each well and wait for ⌚ 00:03:00
- 16 Aspirate all the liquid from each well without disturbing the pellet.
- 17 Dispense 🧴 300 µl of PBS to each well and wait for ⌚ 00:03:00

- 18 Aspirate  **280 µl** of liquid from each well without disturbing the pellet.
- 19 Transfer plate to shaker from the magnet block and shake the plate at  **1200 rpm, Room temperature , 00:01:00**
- 20 Transfer the plate from shaker to Reader side
- 21 Mix the sample using a pipette. Transfer  **10 µl** of the sample to  **60 µl** of liquid crystal.
- 22 Gently mix the sample with liquid crystal and transfer the mix to BioSlide.
- 23 To disperse the liquid crystal, apply positive pressure to the BioSlide.
- 24 After  **00:15:00** , BioSlide is read using the camera reader. The image and event count is obtained for each sample.