



Jul 11, 2020

# S30-S30A-S30-Buffers- Haseloff Lab

Fernando FGC Guzman Chavez<sup>1</sup>, Jim Haseloff<sup>1</sup><sup>1</sup>University of Cambridge**1** Works for me [dx.doi.org/10.17504/protocols.io.bigckbsw](https://dx.doi.org/10.17504/protocols.io.bigckbsw)Fernando Guzman Chavez  
University of Cambridge

## ABSTRACT

Following this recipe, you will obtain 1L of **S30**, **S30A** and **S30B Buffers**

**S30 Buffer** is a standard solution to wash cell pellets previous to lysis. Likewise, it is used during the cell lysis.

## NOTE :

The protocol described here is an adaptation from these papers:

- **Adam D. Silverman**, Nancy Kelley-Loughnane, Julius B. Lucks, and Michael C. Jewett (2019). *Deconstructing Cell-Free Extract Preparation for in Vitro Activation of Transcriptional Genetic Circuitry*. ACS Synthetic Biology, 403-414. DOI: 10.1021/acssynbio.8b00430.
- **Yang WC**, Patel KG, Wong HE, Swartz JR. (2012). *Simplifying and streamlining Escherichia coli-based cell-free protein synthesis*. Biotechnol Prog. 28(2):413-420. DOI:10.1002/btpr.1509.
- <https://www.protocols.io/researchers/anibal-arce-medina>

## DOI

[dx.doi.org/10.17504/protocols.io.bigckbsw](https://dx.doi.org/10.17504/protocols.io.bigckbsw)

## PROTOCOL CITATION

Fernando FGC Guzman Chavez, Jim Haseloff 2020. S30-S30A-S30-Buffers- Haseloff Lab. **protocols.io**  
[dx.doi.org/10.17504/protocols.io.bigckbsw](https://dx.doi.org/10.17504/protocols.io.bigckbsw)



## KEYWORDS

S30 buffer, cell-free

## LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Jul 10, 2020

## LAST MODIFIED

Jul 11, 2020

## PROTOCOL INTEGER ID

39140

## 1. S30 Buffer

## 1 To prepare S30 buffer, following compounds have to be ready-to-use:

- Tris-Acetate
- Magnesium acetate tetrahydrate
- Potassium acetate
- 1M DTT
- 5M KOH

### 1.1 1.1 Prepare for autoclaving

- Volumes indicated are sufficient for **1 L of S30 buffer**

Compound, Sum Formula	MW [g/mol]	Required amount [g or mL]	Concentration [mM]
Tris-Acetate	181.19	<b>1.81g</b>	10
Magnesium acetate	214.45	<b>3.0g</b>	14
Potassium acetate	98.14	<b>5.89g</b>	60

Preparation: Dissolve in demi-H<sub>2</sub>O, fill up to 1000 ml, autoclave 15 min at 121°C  
**\*adjust to pH=8.2 with 5M KOH (approx. 900µL)**

- **Before use** add 2 mL of **1 M DTT** (2mM, final concentration)
- **For dialysis**, S30 buffer is prepared at **5mM** Tris-Acetate and 1mM DTT. Rest of the components remain at same concentrations.
- To prepare **TBST Buffer** (pH=7.6), dissolve 1 tablet of Tris buffered saline with Tween 20 in 500mL of water and autoclave 15 min at 121°C. Catalogue Number: 91414-10TAB SIGMA.

## 2. S30A and S30B Buffer

## 2 To prepare S30A and S30B buffers, following compounds have to be ready-to-use:

- 100 mM Tris (solution)
- 1:1 Acetic Acid
- 1M DTT

### 2.1 2.1 Prepare for autoclaving

- Volumes indicated are sufficient for **1 L of S30 Buffer A and S30 Buffer B**

#### S30 Buffer A

Compound, Sum Formula	MW [g/mol]	Required amount [g or mL]	Concentration [g/L]
100 mM Tris Base	121.14	<b>500 mL</b>	6.057 [50 mM]
Mg-Glutamate	388.61	<b>5.44g</b>	5.44 [14mM]
K-Glutamate	203.23	<b>12.194g</b>	12.194 [60mM]

Dissolve in demi-H<sub>2</sub>O, fill up to 1000 ml, autoclave 15 min at 121°C  
**\*adjust to pH=7.7 with acetic acid**

- Before use add 2 mL of 1 M DTT

### S30 Buffer B

Compound, Sum Formula	MW [g/mol]	Required amount [g or mL]	Concentration [g/L]
100 mM Tris Base	121.14	50 mL	0.6057 [5 mM]
Mg-Glutamate	388.61	5.44g	5.44[14mM]
K-Glutamate	203.23	12.194g	12.194[60mM]

Dissolve in demi-H<sub>2</sub>O, fill up to 1000 ml, autoclave 15 min at 121°C  
**\*adjust to pH=8.2 with acetic acid**

- Before use add 2 mL of 1 M DTT

### 3. Literature/ References

3



Silverman AD, Kelley-Loughnane N, Lucks JB, Jewett MC (2019). Deconstructing Cell-Free Extract Preparation for in Vitro Activation of Transcriptional Genetic Circuitry.. ACS synthetic biology.  
<https://doi.org/10.1021/acssynbio.8b00430>



Yang WC, Patel KG, Wong HE, Swartz JR (2012). Simplifying and streamlining Escherichia coli-based cell-free protein synthesis.. Biotechnology progress.  
<https://doi.org/10.1002/btpr.1509>

- <https://www.protocols.io/researchers/anibal-arce-medina>

### 4. Change history

- 4 V1→V1.1 TBST recipe added.  
V1.1→V1.2 Recipes to prepare S30A and S30B buffers were added.