



Version 1

Sep 21, 2021

Construction of individual ddRAD libraries V.1

Claire Daguin Thiebaut¹, Stephanie Ruault¹, Charlotte Roby¹, Thomas Broquet¹,
Frédérique Viard¹, Alan Brelsford²

¹Station Biologique de Roscoff, France; ²Dept of Biology, U. of California, Riverside, USA

Claire Daguin Thiebaut: CNRS, Sorbonne Université, UMR7144 Adaptation et Diversité en Milieu Marin;

Stephanie Ruault: CNRS, Sorbonne Université, UMR7144 Adaptation et Diversité en Milieu Marin

Charlotte Roby: present address: Vegenov, Saint-Pol de Léon, France

Thomas Broquet: CNRS, Sorbonne Université, UMR7144 Adaptation et Diversité en Milieu Marin

Frédérique Viard: present address: ISEM, University of Montpellier, France

1 Works for me

Share

dx.doi.org/10.17504/protocols.io.bv4tn8wn

DYDIV



Claire Daguin Thiebaut
Station Biologique de Roscoff, France

ABSTRACT

This protocol describes a double digested restriction-site associated DNA (ddRADseq) procedure, that is a variation on the original RAD sequencing method ([Davey & Blaxter 2011](#)), which is used for *de novo* SNP discovery and genotyping.

This protocol differs from the original ddRADseq protocol ([Peterson et al 2012](#)), in which the samples are pooled just after the ligation to adaptors (i.e. before size selection and PCR).


The present ddRAD protocol has been slightly adapted from Alan Brelsford's protocol published in the supplementary material of this paper:

Brelsford, A., Dufresnes, C. & Perrin, N. 2016. High-density sex-specific linkage maps of a European tree frog (*Hyla arborea*) identify the sex chromosome without information on offspring sex. *Heredity* **116**, 177–181 (2016).
<https://doi.org/10.1038/hdy.2015.83>


In the present protocol, all samples are treated separately, in a microplate, until final PCR amplification performed before pooling. Despite being slightly more costly and time-consuming in the lab, it allows for fine adjustment of each sample representation in the final library pool, ensuring similar number of sequencing reads per sample in the final dataset.

Briefly, genomic DNA from the samples are individually digested with 2 restriction enzymes (one rare-cutter and one more frequent cutter) then ligated to a barcoded adaptor (among 24 available) at one side, and a single adaptor at the other side, purified with magnetic beads, and PCR-amplified allowing the addition of a Illumina index (among 12 available) for multiplexing a maximum of 288 sample per library. Samples are then pooled in equimolar conditions after visualisation on an agarose gel. Purification and size selection is then performed before final quality control of the library and sequencing.

Double-digested restriction site associated DNA sequencing (ddRADseq) at a glance



individual ddRAD at a glance



Sequences of adaptors used in the present protocol (Alicia Mastretta-Yanes, pers. comm.).

DOI

dx.doi.org/10.17504/protocols.io.bv4tn8wn

PROTOCOL CITATION

Claire Daguin Thiebaut, Stephanie Ruault, Charlotte Roby, Thomas Broquet, Frédérique Viard, Alan Brelsford 2021. Construction of individual ddRAD libraries. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bv4tn8wn>
Version created by **Claire Daguin Thiebaut**

KEYWORDS

ddRAD, SNP, population genomics, DNA, RADseq, genomic

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited




CREATED

Jun 25, 2021

LAST MODIFIED

Sep 21, 2021

OWNERSHIP HISTORY

Jun 25, 2021		Stephanie Ruault	Station Biologique de Roscoff, France
Sep 07, 2021		Stephanie Ruault	Station Biologique de Roscoff, France
Sep 07, 2021		Claire Daguin Thiebaut	Station Biologique de Roscoff, France

PROTOCOL INTEGER ID

51059

GUIDELINES

This protocol has been used successfully for several marine macro-algae ([Guzinski et al 2018](#), [Le Cam et al 2019](#)) and marine vertebrate (ivory gull, Charbonnel et al, *in revision Heredity*), invertebrate species (isopods, tunicates (Le Moan et al, *in revision, Mol Ecol*), mollusks (Tran Lu Y et al, *in revision, Mol Ecol*). Efficiency may rely on specific genome composition (GC rate, genome size, genome complexity). The main advantage of this protocol is to ensure homogeneous representation of the different samples in the final library pool.

The choice of restriction enzymes depends on the genome size and content. In the case of very large genomes, SbfI may replace PstI, using the same P1 adaptors set. For genomes with low GC content, use MspI instead of MseI.

The number of individual per library may be adjusted according to the expected yield of the sequencing run and the desired depth of sequencing per individual.

Starting DNA amount may also be adjusted to the number of individual multiplexed in a pool. This protocol has been successfully used with as little as 50 ng of genomic DNA for each individual. Strict normalization of all DNA extracts' concentrations is not required because of the normalizing PCR included at the end of the protocol.

Several individuals should be included in triplicates, from the digestion (or better, DNA extraction) step, allowing better adjustment of reads assembly parameters (for more explanations see [Mastretta-Yanes et al 2015](#) and [Paris et al 2017](#)).

MATERIALS TEXT

Reagents :

- 100 % ethanol (molecular biology grade)
- Agarose
- TBE 1X buffer

- Ethidium Bromide Solution [10mg/ml]
- 100bp ladder
- Gel loading dye
- DNase-Rnase free ultrapure water
- PstI-HF® ref NEB : R3140S
- MseI ref Neb : R0525S
- T4 DNA Ligase ref NEB : M0202L
- ATP 10mM ref NEB : P0756
- Annealing buffer stock (10X): 100 mM Tris HCl, pH 8, 500 mM NaCl, 10 mM EDTA
- dNTP set 4x0,25mlx100mM ref Neb : N0446S
- Q5® Hot Start High-Fidelity DNA Polymerase ref NEB : M0493L
- Ampure XP beads Beckman coulter (A63881) or NGS clean up and size selection (Macherey Nagel REF 744970.50)
- Elution solution/buffer for magnetic beads (Tris 10mM)
- Freshly diluted 70% ethanol
- Pippin Prep reagents and cassettes (1,5 % DF Marker K, ref SageScience CDF1510)
- NEBNext Library Quant kit for illumina ref Neb : E7630
- Qubit 1X dsDNA BR (Thermo, Q33266)

Primer and adaptor sequences :

NB : All primers and adaptors are NGS grade oligonucleotides, HPLC purified, delivered at 100µM.

▪ Adaptors

Adapter P1: annealed oligos 1.1 and 1.2 (1 pair per barcode, 24 6-base barcodes) :

PstI adapter p1.1

ACACTCTTTCCCTACACGACGCTCTTCCGATCTnnnnnnTGCA

PstI adapter p1.2

nnnnnnAGATCGGAAGAGCGTCGTGTAGGGAAAGAGTGT

 [P1-PstI_sequences.pdf](#)

Adapter P2: annealed oligos 2.1 and 2.2 :

MseI adaptor 2.1: GTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT

MseI adaptor 2.2: /5Phos/**TA**AGATCGGAAGAGCGAGAACAA

or for AT rich genomes:

MspI adaptor 2.1 : GTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT

MspI adaptor 2.2: /5Phos/**CG**AGATCGGAAGAGCGAGAACAA

▪ PCR primer oligos:

“ILLPCR1” : forward primer common to all reactions; “ILLPCR2” reverse primers contain unique 6-bases indices that allow for multiplexing in a single pool in an illumina lane. Each individual is tagged with one inline barcode, and one Illumina index. Therefore, with 24 barcodes and 12 indices, 288 individuals can be multiplexed in a single pool. One index read is necessary during the illumina sequencing run (either single read or paired-end).

ILLPCR1 primer:

aaTGATACGGCGACCAACGAGATCTACACTCTTCCCTACACGACG

Truseq Indexed illPCR2 Primers: each unique pcr2 primer lets you reuse the 24 barcodes in the same lane. If using two indexed primers, Illumina recommends to use 6 and 12. If three primers, use 4, 6, 12. If four primers, use 1, 4, 7, 9 or 6, 8, 11, 12. If six primers, 2, 4, 5, 6, 7, 12.

ILLPCR2_ind01

caAGCAGAAGACGGCATAACGAGAT**CGTGAT**GTGACTGGAGTTCAGACGTGTgc

ILLPCR2_ind02

caAGCAGAAGACGGCATAACGAGAT**ACATCG**GTGACTGGAGTTCAGACGTGTgc

ILLPCR2_ind03

caAGCAGAAGACGGCATAACGAGAT**GCCTAA**GTGACTGGAGTTCAGACGTGTgc

ILLPCR2_ind04

caAGCAGAAGACGGCATAACGAGAT**TGGTCA**GTGACTGGAGTTCAGACGTGTgc

ILLPCR2_ind05

caAGCAGAAGACGGCATAACGAGAT**CACTGT**GTGACTGGAGTTCAGACGTGTgc

ILLPCR2_ind06
 caAGCAGAAGACGGCATAACGAGAT**ATTGGCGT**GACTGGAGTTCAGACGTGTgc
 ILLPCR2_ind07
 caAGCAGAAGACGGCATAACGAGAT**GATCTG**GACTGGAGTTCAGACGTGTgc
 ILLPCR2_ind08
 caAGCAGAAGACGGCATAACGAGAT**TCAAGT**GACTGGAGTTCAGACGTGTgc
 ILLPCR2_ind09
 caAGCAGAAGACGGCATAACGAGAT**CTGATC**GACTGGAGTTCAGACGTGTgc
 ILLPCR2_ind10
 caAGCAGAAGACGGCATAACGAGAT**AAGCTA**GACTGGAGTTCAGACGTGTgc
 ILLPCR2_ind11
 caAGCAGAAGACGGCATAACGAGAT**GTAGCC**GACTGGAGTTCAGACGTGTgc
 ILLPCR2_ind12
 caAGCAGAAGACGGCATAACGAGAT**TACAAG**GACTGGAGTTCAGACGTGTgc

small letters indicate PTO modifications during synthesis of primers.

Consumables :

PCR plates
 8-tubes and 12-tubes *PCR* strips; Volume: 0.2mL
 thermosealable films
 2 mL low binding microcentrifuge tubes
 Filter tips (better in low retention quality)
 qPCR white plate 96
 Agilent DNA High sensitivity DNA kit (5067-4626)
 Qubit microtubes

Equipment :

Manual pipettors, monochannel and multichannel
 Centrifuge for microcentrifuge tubes and plates
 Vortex mixer
 Personal protection equipment (lab coat, gloves, goggles)
 Plate sealer
 Microwave
 Thermocycler
 Gel Electrophoresis Systems
 magnetic stand 96 (ref Thermo AM10027)
 Pippin prep (Sage Science)
 Qubit fluorimeter
 qPCR equipment (e.g. LightCycler 480 Roche)
 Bioanalyzer (Agilent)

SAFETY WARNINGS





Fully respect the strict separation between pre- and post-PCR labs. Use low retention filter tips for all pipeting steps.

Be very careful to avoid cross contamination between samples. Include non-template controls in each steps.

This protocol uses Ethidium Bromide (EtBr), which is a known mutagen. Wear a lab coat, eye protection and gloves when working with this chemical.

Preparation of double-stranded barcoded P1 adaptors 1µM 5m

5m

- 1 In a PCR plate wells, combine each oligo 1.1 with its complementary oligo 1.2 :
 -  **5 µl** adaptor P1-1 10µM
 -  **5 µl** adaptor P1-2 10µM
 -  **5 µl** annealing buffer 10X (100 mM Tris-HCl, pH 8, 500 mM NaCl, 10 mM EDTA Na₂)
 -  **35 µl** nuclease free water

2

		25			9				17		
		2			10				18		
		3			11				19		
		4			26				20		
		5			13				21		
		6			14				22		
		7			15				23		
		8			16				24		

Preparation of double-stranded P2 Adaptor P2_Mse 10µM	5m
---	----

Double digestion	10s
------------------	-----

- Remove the NucleoMag® NGS Bead Suspension (or Beckman Ampure XP) from the fridge. Let for approximately 30 min to bring the bead suspension to room temperature.
- Then, vortex this Bead Suspension stock solution carefully until homogeneization and fill all wells of a 8-tube PCR strip with this solution, or use a reagent reservoir.
- Add **60 µl** of well-homogenized bead suspension into the first 4 columns of a new PCR microplate with a multichannel pipette.
- Spin down the plate containing the digested-ligated products. Fill each well of a 8-PCR tube strip with Tris-HCl pH8 10mM.
- In the DNA digestion plate, add **5 µl** of Tris 10mM with a multichannel pipette, to be sure that you will be able to take 40µl of digested-ligated DNA.

11 Binding

This step binds DNA fragments 200bp and larger to the magnetic beads.

- For the first 4 columns of the digestion-ligation plate, Transfer **40 µl** of the DNA samples in the first 4 columns of the purification plate containing beads and carefully mix by pipetting up and down 10 times.
- Incubate the purification plate at room temperature for 5 min.

12 Separation

5m

- Place the purification plate onto the 96-well magnetic separator.
- Wait at least **00:05:00** until all the beads have been attracted by the magnets or until the liquid appears clear.
- The supernatant contains unwanted low molecular weight contaminants and unwanted smaller DNA fragments.
- Remove and discard the supernatant (ca 100 µl) by pipetting.
- Do not disturb the attracted beads while aspirating the supernatant. Remove the supernatant with the multichannel from the opposite side of the well.

13 1st wash with 80 % ethanol

- Place 80% (or 70%) ethanol in a 25ml reservoir.
- With a multichannel pipette, dispense **200 µl** 80% (or 70%) ethanol into the purification plate without disturbing the bead pellet.
- Incubate the purification plate at room temperature for at least 30 s.
- Carefully and completely remove and discard ethanol by pipetting.

14 2nd wash with 80 % ethanol

- With a multichannel pipette, dispense **200 µl** 80% (or 70%) ethanol into the purification plate without disturbing the bead pellet.
- Incubate the purification plate at room temperature for at least 30 s.
- Carefully remove and discard ethanol by pipetting.

15 Dry the beads

Let the purification plate on the magnetic separator and incubate at room temperature **for maximum 5 min** in order to allow the remaining traces of ethanol to evaporate.

NOTE: take care not to over dry the bead pellet (bead pellet appears cracked in this case) as this will significantly decrease elution efficiency.

16 Elute DNA fragments library

10m

- Take the purification plate from the magnetic stand, and add **40 µl** of Tris 10mM pH 8 with a multichannel pipette to resuspend the bead pellet by pipetting up and down 10 times.
- Incubate the purification plate at room temperature for **00:05:00**.
- Separate the magnetic beads against the side of the wells by placing the 96-well plate on the magnetic separator.

- 17 In this PCR, Illumina indexed primers are incorporated in order to produce fragments compatible with Illumina sequencing, and to insert an index allowing multiplexing of barcoded samples (see ligation step for details). This index will be read during the sequencing run.
- This PCR is expected to have a homogenizing effect. Primers are thus included in limiting quantity, in order to produce equalized amounts of PCR fragments among samples. The number of cycles is limited to a maximum of 15 (optimal with 12). After those cycles, a final PCR cycle is then performed after addition of primers in large excess.

The Reaction mixture for a total volume of 40μL is:

A	B	C	D
	conc initiale	conc final	vol 1X μ l
Q5 buffer NEB	5X	1x	8
dNTP mix	25mM each	0.2 mM	0.3
primer mix (illPCR1 and illPCR2index)	5 μ M each	0.17 μ M	1.4
Q5 hotstart hifi DNA polymerase NEB	2u/ μ l		0.4
High GC enhancer NEB	5X	1X	8
nuclease-free water			11.9
Adaptor-ligated purified template DNA			10
total mix			30
TOTAL reaction			40

10s

Spin down all reagents for approximately 00:00:05 and place on ice.

A	B	C
	Vol µl 1X µl	Vol µl for 1 plate (104 reactions)
Q5 buffer 5X	8	832
dNTP mix 25mM	0.3	31.2
Q5 hot start DNA polymerase 2u/µl	0.4	41.6
High GC enhancer 5X	8	832
nuclease free water	3.3	343.2
total	20	2080

Vortex mix all reagents in the mix and spin down.

Aliquot **130 µl** of the 1st mix into each well of a 8-tube *PCR* strip.

In a new *PCR* plate, dispense **20 µl** of 1st mix with a multichannel in each column.

20 2nd mix preparation (primer *PCR* mix) :

prepare 1 mix for each index in a 1.5 ml microtube, (4 mix for a full a plate).
for 1 index

A	B	C
	vol 1x µl	Vol µl for 3 columns
illPCR Primer mix 5µM	1.4	39.2
nuclease-free water	8.6	240.8
total	10	280

Second mix composition (primer *PCR* mix)

Vortex mix all reagents and spin down.

In a 12-tubes *PCR* strip, aliquot **93 µl** of the primer *PCR* mix into 3 consecutive wells according to the following scheme :



ind 1 Ind 4 ind 7 ind 9

In the *PCR* plate, dispense **10 µl** of 2nd mix in each line of the plate with a 12 multichannel pipette, onto the 20µl of the first mix.

21 DNA and mix combination

Spin down the ligation purified DNA plate.

With a multichannel pipette, transfer **10 µl** of purified ligated DNA into the *PCR* plate and mix by pipetting.

Finally, aliquot the **40 µl** by dispensing **10 µl** into 3 additional new empty *PCR* plates.

Seal the 4 *PCR* plates and spin down.

The 4 *PCR* will be performed in parallel In 4 different thermal cyclers, in order to reduce the *PCR* bias.

A	B	C	D
hot start initial denaturation	98°C	30s	
amplification cycles			
denaturation	98°C	20s	
annealing	60°C	60s	15 cycles
extension	72°C	40s	
final extension	72°C	10min	
hold	12°C		

PCR program for the Illumina indexing PCR

Seal the 4 PCR plates and spin down.

The 4 PCR will be performed in parallel in 4 different thermal cyclers, in order to reduce the PCR bias.

After PCR, pool back the 4 PCR plates into a single plate with a multichannel.

22 Final cycle (with primers in large excess)

Prepare 1 mix for each index in a 1.5 ml microtube

A	B	C	D	E
	conc stock	final concentration	vol µl 1X	vol µl for 3 columns
Q5 buffer	5X	1X	0.2	24
dNTP mix	25mM each	0.2mM	0.08	9.6
illPCR primer mix	5µM each	0.2µM	0.4	48
nuclease-free water			0.32	38.4
total			4	120

Vortex mix all reagents and spin down. In a 12-tube PCR strip, aliquot **40 µl** of the mix into 3 consecutive wells according to the following scheme:



ind 1 Ind 4 ind 7 ind 9

Dispense **4 µl** of final cycle mix in each line of the plate with a 12 multichannel pipette.

Seal the PCR plate and spin down.

In a thermocycler, run the final cycle as follows:

Warning :

If you perform the size selection on the pool without prior bead clean-up, you may observe a shift between observed and expected size ranges of the smear (approximately 250-750 obtained instead of 300-800).


If you see adaptor dimers on the gel after PCR, it is also preferable to perform a bead clean-up before size selection in the Pippin Prep.

In the case you do not have access to a Pippin prep, you can alternatively perform a double size-selection with beads (see the bead manufacturer's instructions for details), or by smear excision and purification from an agarose gel. From our experience, size selection with the Pippin prep is the most accurate and repeatable method, and has to be preferred if you have several pools of libraries (for example 3 pools of 192 samples) in order to increase the number of common loci between samples in the ddRAD sub sampling of the genome.

Quality control of the libraries

- 27 Control the quality of the library** with a Bioanalyzer (Agilent) in a High Sensitivity DNA chip. Dilute your pool 1:2 or more and load 1 µl of the pool before and after size selection, according to the manufacturer's instructions :

☐ [HighSensitivity_DNA_KG.pdf](#) ☒



29.1 *Suggestions to prepare library dilutions :*

in a 8-tube PCR strip:

Prepare intermediate dilutions (1:10 and 1:100) of the library with the dilution buffer supplied in the qPCR kit.

Then prepare the 4 library dilutions to be used in triplicate for qPCR analysis :

1:1000 : 10µl of 1/100 + 90µl buffer 1X

1:2000 : 50µl of 1/1000 + 50µl buffer 1X

1:4000 : 50µl of 1/2000 + 50µl buffer 1X

1:8000 : 50µl of 1/4000 + 50µl buffer 1X

You should get more than 10 nM, that is the library concentration usually required by the sequencing platform facilities.

The library is now ready for sequencing in single read or paired-end 150 bases in an Illumina sequencer, with one index read. Use the average size of the library size range as estimated from the Bioanalyzer profile to convert DNA concentration from nM to ng/µl.