

May 02, 2022

## Immunofluorescent staining of pancreatic sections

Caroline CT Tremblay<sup>1</sup>, Valentine Moull  <sup>1</sup>, Bader Zarrouki<sup>1</sup>, Julien Ghislain<sup>1</sup>, Vincent Poitout<sup>2</sup>

<sup>1</sup>Montreal Diabetes Research Center, CRCHUM, Montr  al, QC, Canada.;

<sup>2</sup>Montreal Diabetes Research Center, CRCHUM, and Department of Medicine, Universit   de Montr  al, Montr  al, QC, Canada.

1



[dx.doi.org/10.17504/protocols.io.bg6ujzew](https://dx.doi.org/10.17504/protocols.io.bg6ujzew)

### Laboratory of Vincent Poitout

Tech. support email: [julien.ghislain.chum@ssss.gouv.qc.ca](mailto:julien.ghislain.chum@ssss.gouv.qc.ca)

Julien Ghislain

This protocol describes the steps for performing fluorescence immunohistochemistry on fixed-frozen pancreatic tissue sections. It is suitable for pancreatic tissue isolated from rats and mice at postnatal to adult stages. We routinely apply this protocol to quantify total and proliferating or apoptotic beta-, alpha- and delta-cells. Briefly, pancreata are fixed in 4% paraformaldehyde solution and cryoprotected overnight in 30% sucrose. Tissue are then embedded, frozen, sectioned and mounted on slides. Antigen retrieval is performed using sodium citrate buffer prior to immunostaining.

DOI

[dx.doi.org/10.17504/protocols.io.bg6ujzew](https://dx.doi.org/10.17504/protocols.io.bg6ujzew)

Caroline CT Tremblay, Valentine Moull  , Bader Zarrouki, Julien Ghislain, Vincent Poitout 2022.

Immunofluorescent staining of pancreatic sections. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.bg6ujzew>

### CIHR

Grant ID: 0406014428

### NIH

Grant ID: 5R01DK058096-15

protocol

Glucose and fatty acids synergistically and reversibly promote beta cell proliferation in rats. Moull   VS, Vivot K, Tremblay C, Zarrouki B, Ghislain J, Poitout V. Diabetologia. 2017 May;60(5):879-888. doi: 10.1007/s00125-016-4197-8. Epub 2017 Jan 11. PMID: 28078385.

Epidermal growth factor receptor signaling promotes pancreatic  $\beta$ -cell proliferation in response to nutrient excess in rats through mTOR and FOXM1. Zarrouki B, Benterki I, Font  s G, Peyot ML, Seda O, Prentki M, Poitout V. Diabetes. 2014 Mar;63(3):982-93. doi: 10.2337/db13-0425. Epub 2013 Nov 5. PMID: 24194502.

immunohistochemistry, pancreas

protocol ,

Caroline Tremblay

Jun 04, 2020

May 02, 2022

37812

The slides should not be allowed to dry during the staining steps. Drying out will cause nonspecific antibody binding and therefore high background staining.

#### MATERIALS

[Hoechst 33342, Trihydrochloride, Trihydrate - 10 mg/mL Solution in Water](#) **Invitrogen - Thermo**

**Fisher Catalog #H3570**

[OCT \(Optimal Cutting Temperature compound\)](#) **Sakura**

**Finetek Catalog #4583**

[VECTASHIELD® Hardset™ Antifade Mounting Medium](#) **Vector**

**Laboratories Catalog #H-1400**

[Phosphate Buffered Saline 10x \(solution\)](#) **Bio Basic**

**Inc. Catalog #PD8117**

[Paraformaldehyde 10% buffered](#) **Newcomer**

**Supply Catalog #13301A**

[Donkey serum](#) **Sigma**

**Aldrich Catalog #S30-100ML**

[Triton X-100](#) **Sigma**

**Aldrich Catalog #T9284**

[Dako Pen Delimiting pen](#) **Agilent**

**Technologies Catalog #S2002**

[Bovine serum albumin](#) **Sigma**

**Aldrich Catalog #A7888**

[Sucrose Ultra](#)

**Pure Bioshop Catalog #SUC507.5**

[Sodium citrate crystals reagent](#) **ACP**

**Chemicals Catalog #S2990**

[Superfrost Plus Microscope Slides](#) **Fischer**

**Scientific Catalog #12-550-15**

[MilliQ water](#) **Contributed by users**

Surgipath® Clear Disposable Base Molds

Leica

75809-376

Well, 250ml. w/ lid, green, xylene resistant

TBS

SS-WLG

When working with PFA, always work in the chemical hood.

Make sure you have the right combination of secondary and primary antibodies with regards to the host species and the fluorophore.

## Preparation of cryosections

1d

### 1 Tissue fixation

19h



When working with PFA, always work in a chemical hood.

Harvest the pancreas and place it in a 50 ml Falcon tube containing **30 mL** of cold 4 % PFA ( **12 mL** of 10 % PFA + **18 mL** of PBS).

Fix for **04:00:00** at **4 °C** in the dark.

Then, working in a chemical hood, delicately remove the pancreas with forceps and place it on brown paper to absorb the fixative.

Place the organ in a new 50 ml Falcon tube containing **30 mL** of a 30 % sucrose solution ( **9 g** of sucrose and **30 mL** of PBS).

Store **15:00:00** at **4 °C** in the dark.

The next day, delicately remove the pancreas with forceps and place it on brown paper to absorb the sucrose solution. Place it in a mold, cover with OCT and freeze at **-80 °C** until ready for sectioning.

### 2 Preparing cryosections

2h

Set the cryostat and the pedestal temperature at **-20 °C**.

Gather all of the needed material (OCT, slides, pencil, blades, paintbrushes, aluminium foil, slide box, tissues).

Cut cryosections at 0.8 µm thickness and collect on Superfrost Plus microscopic slides.

Store the sections at **-80 °C** until ready for staining.

## Immunofluorescent staining

1d

### 3 Antigen retrieval

1h

Bring the slides to **Room temperature** by putting them in the slide holder in PBS until ready to perform antigen retrieval step.

Put the holder in a 250 ml well filled with milliQ water for **00:05:00** at **Room temperature**.

Transfer the slides to the sodium citrate solution ( **10 millimolar (mM)** **pH 6** ).

Heat the slides in the microwave with the following sequence:

Heat **00:03:00** (if the solution starts to boil, stop heating and add the remaining time to the pause section; the antigen retrieval step must last a total of **00:20:00** ).

Remove from the microwave and let stand **00:17:00**.

Heat **00:01:00** and let cool down for **00:30:00** at **Room temperature**.

Rinse the slides in milliQ water for **00:05:00**.

Transfer slides to PBS for **00:05:00**.

### 4 Blocking

1h

Dry the slides and outline the tissue with a hydrophobic barrier pen.

Place the slides with the tissue facing up in an humidified chamber and add a sufficient volume of blocking solution ( **0.1 % (v/v)** Triton, **1 % (v/v)** BSA plus **5 % (v/v)** normal donkey serum in PBS). Use between **150 µL** - **300 µL** depending on the size of the tissue.

Block for **01:00:00** at **Room temperature**.

## 5 Primary antibody staining

15h


At no time should the slides be allowed to dry. Drying out will cause nonspecific antibody binding and therefore high background staining.



Remove the blocking solution and replace with  150 µL -  300 µL antibody mix diluted in blocking solution.

Incubate in the humidified chamber  15:00:00 at  4 °C .

## 6 Secondary antibody staining

1h 30m





The following morning decant the primary antibody mixture and wash the slides three times in PBS for  00:05:00 using a slide holder and well.

Return slides to the humidified chamber and incubate with the secondary antibody diluted in blocking solution for  01:00:00 at  Room temperature .

Remove the antibody mixture and wash three times  00:05:00 in PBS.

## 7 Hoechst staining and mounting

30m

Label nuclear DNA by incubating the slides in the Hoechst solution (  10 µL of Hoechst in  250 mL of PBS) for  00:10:00 at  Room temperature using slide holder and well.

Wash the slides three times  00:05:00 in PBS.

Mount slides using VECTASHIELD HardSet Mounting Medium without DAPI.

Place the mounted slides in a cardboard slide tray holder in the dark until ready to acquire the images with a fluorescent microscope.