

NOV 12, 2023

## Preparation and Administration of Systemic 4-AP

Leonardo A Parra-

Rivas<sup>1</sup>

<sup>1</sup>University of California, San Diego

Subhojit Roy Lab



Leonardo A Parra-Rivas University of California, San Diego

**ABSTRACT** 

Preparation and Administration of Systemic 4-AP





DOI:

dx.doi.org/10.17504/protocol s.io.5qpvo38dxv4o/v1

**Protocol Citation:** Leonardo A Parra-Rivas 2023. Preparation and Administration of Systemic 4-AP. protocols.io https://dx.doi.org/10.17504/p rotocols.io.5qpvo38dxv4o/v1

License: This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status: Working** We use this protocol and it's working

Created: Nov 12, 2023

Last Modified: Nov 12,

2023

**PROTOCOL** integer ID:

90816

1

4-Aminopyridine was obtained from Sigma-Aldrich (Cat #A78403), and a stock solution was prepared in saline. Mice (8 weeks old) were given intraperitoneal injections of soluble 4-AP at a low dose of 2 mg/Kg or equivalent saline.

Animals were monitored every 20 min and presented no sign of convulsive seizure behavior. After 3 hrs, animals were anesthetized with isoflurane (VetOne Cat#502017), followed by cervical dislocation. Brains were rapidly removed and drop-fixed in 4% PFA solution for immunohistochemistry or placed on ice for biochemical assays.