




JAN 24, 2023

FBS heat inactivation

Andreas Sagen¹

¹University of Oslo



Andreas Sagen

University of Oslo, The National Institute of Occupational H...

ABSTRACT

The objective of heat inactivation is to destroy complement activity in the serum without affecting the growth-promoting characteristics of the product. Removal of complement activity from serum, such as fetal bovine serum, is not required for most cell cultures but may be necessary for cultures that are sensitive to the complement activity. Since heat inactivation of the serum may, to some extent, decrease the growth performance properties of the serum, this procedure should only be performed if required for optimal cell growth. Researchers should evaluate the applicability of heat inactivation regarding their own application.

GUIDELINES

Do not use heat inactivated FBS unless necessary.

MATERIALS

Circulating hot water bath

Timer

Thermometer

Ice bath

SAFETY WARNINGS



58 degrees Celsius is hot. Take care when working with these temperatures.

BEFORE START INSTRUCTIONS

Clean and sanitize everything

OPEN ACCESS

Protocol Citation: Andreas Sagen 2023. FBS heat inactivation. [protocols.io](https://protocols.io/view/fbs-heat-inactivation-cm85u9y6) <https://protocols.io/view/fbs-heat-inactivation-cm85u9y6>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Jan 24, 2023

Last Modified: Jan 24, 2023

PROTOCOL integer ID:
75773



Keywords: FBS, Fetal bovine serum, Inactivation



Inactivation

2d 3h 24m

1 Preparation

1d 0h 42m

- 1.1 Slow thaw FBS at  4 °C overnight 1d
- 1.2 Aliquot FBS into  50 mL or  100 mL units in a sterile environment 10m
- 1.3 Prepare an volumetrically equal blank container with distilled water 2m
- 1.4 Equalibriate FBS and distilled water for  00:30:00 at  Room temperature 30m
- 2 Prepare a water bath 30m
- 2.1 Clean and add distilled water to an appropriate hight 20m
- 2.2 Set temperature to  58 °C 10m
- 3 FBS inactivation
- 3.1 Set FBS aliquots and distilled water blank into preheated water bath at  58 °C for  00:30:00 Use a thermometer to track temperature in distilled water blank, and start timer 35m

when temperature reach  56 °C . Agitate aliquots every  00:05:00 by mixing

3.2 Remove aliquots from heat bath, and place FBS  On ice or  4 °C for  00:30:00 

3.3 Store aliquots at  -20 °C