



Mar 14, 2022

🌐 SP3 protocol for proteomic analysis of tendon cryosections

Stefania Wunderli¹, Jess Snedeker², The Tendon Seed Network²

¹ETHZ - ETH Zurich; ²University of Zurich

1



protocol .

Tendon Seed Network

Sarah Snelling

This protocol details the preparation of cryosection samples and laser capture microdissection (LCM) samples collected with a MMI Laser-Microdissection device for MS/MS analysis. The sample preparation includes tissue lysis, Sp3 technology (single-pot, solid-phase-enhanced sample preparation), trypsin digestion, C18-clean-up and peptide resuspension.

[ft3rbkj7p.docx](#)

Stefania Wunderli, Jess Snedeker, The Tendon Seed Network 2022. SP3 protocol for proteomic analysis of tendon cryosections. **protocols.io**
<https://protocols.io/view/sp3-protocol-for-proteomic-analysis-of-tendon-cryo-bzndp5a6>



Chan Zuckerberg Initiative
Grant ID: CZI2019-002426

Proteomics, Cryosection, Tissue Lysis

protocol ,

Nov 01, 2021

Mar 14, 2022

Nov 01, 2021







Maria Gladius

Nov 02, 2021

Sarah Snelling




54693

Reagents:

- Needle (1.2x40mm)
- RIPA buffer (Sigma)
- [M]100 millimolar (mM) TCEP
- [M]500 millimolar (mM) IAA
- MS water
- 80 % EtOH
- 100 % EtOH
- Carboxylated Magnetic Beads (hydrophilic and hydrophobic)
 - SpeedBead Magnetic Carboxylated, modified particles, Ref: 65152105050250 ( 15 mL , azide 0.05%, Sigma (??))
 - SpeedBead Magnetic Carboxylated, modified particles, Ref: 45152105050250 ( 15 mL , azide 0.05%, Sigma (??))
- Ammonium bicarbonate (ABC) [M]500 millimolar (mM)
 -  100 µL of [M]500 millimolar (mM) ABC +  900 µL MilliQ
- Trypsin ([M]0.1 µg/µL in [M]50 millimolar (mM) ABC) -> dissolved in  200 µL
- 5% TFA
- MS buffer (3% ACN, 0.1% formic acid)
- iRT peptides (Biognosys) ( 4 °C)
- Waters glass vial

Preparing cryosection samples

1h 11m

- 1 Fill a  200 µL PCR tube with  100 µL 100% EtOH.
- 2 Wetten the dry tissue section on the glass slide with ca.  50 µL 100% EtOH.
- 3 Use a 1.2x40mm needle (18Gx1.5) to gently scrape off the tissue section
- 4 Transfer the scraped off tissue to EtOH in PCR tube.

5



1m

Spin down ⌚ 00:01:00 at 🌀 12000 x g .

6



Remove 100% EtOH supernatant from tubes by using a 📏 200 µL pipette.

7

Heat the tubes on thermocycler at 🔥 50 °C for ca. ⌚ 00:10:00 .

10m

8

Continue with A) or B).

Step 8 includes a Step case.

A) Tissue lysis LCM samples

B) Tissue lysis cryosection samples

step case

A) Tissue lysis LCM samples

9



Add 📏 20 µL RIPA lysis buffer to each sample: Add to tube, close tubes and shake down the liquid onto cap.

9.1

10

Preparing cryosection samples

1h 21m 30s

11



1h

Incubate samples at 🔥 95 °C on thermocycler for ⌚ 01:00:00 .

11.1 LCM samples: upside-down, using magnetic rack for SP3 as a holder and stick tubes to rack with a tape.

12



10m

Centrifuge LCM samples for ⌚00:10:00 at 🌀2000 x g .

13



Take 📏5 µL of RIPA in each sample tube, pipet onto cap and scrape off remaining tissue from cap. Pipet RIPA back into tube. (LCM samples only).

14



30s

Spin down for ⌚00:00:30 at 🌀2000 x g .

15



Check under microscope, whether tissue pieces are still sticking to the cap.

16

Pool all the blank samples into one tube (ca. 📏40 µL in total).

Reduction and alkylation

30m

17

Final volume/concentrations:

- [M]**2 millimolar (mM)** TCEP: Add 📏**0.4 µL** of [M]**100 millimolar (mM)** TCEP solution to
📏**20 µL** of sample.

- [M]**15.5 millimolar (mM)** IAA: Add 📏**0.622 µL** of [M]**500 millimolar (mM)** IAA

solution to $20\ \mu\text{L}$ of sample.

Prepare a reduction/alkylation stock solution for TCEP and IAA with the final concentrations from above.

$8.8\ \mu\text{L}$ TCEP + $13.684\ \mu\text{L}$ IAA

18




Add $1.022\ \mu\text{L}$ of reduction/alkylation stock solution to each sample (if volume is $20\ \mu\text{L}$),
add $2.044\ \mu\text{L}$ to blank.

19



Vortex

20



30m

Incubate samples at $60\ ^\circ\text{C}$ for 00:30:00 at 600 rpm .

21

Spin down samples quickly.

Sp3 assisted protein capture and clean-up

45m

22

Hydrophilic and hydrophobic beads come at a stock concentration of $50\ \mu\text{g}/\mu\text{L}$. Use $1\ \mu\text{L}$ of stock solution ($10\ \mu\text{g}/\mu\text{L}$) for each μg of protein.

LCM samples: assuming a protein amount of $10\ \mu\text{g}$ per sample we need $10\ \mu\text{L}$ stock solution per sample: $15(+1)$ samples x $10\ \mu\text{L}$ = $160\ \mu\text{L}$ stock.

Cryosection samples: ca. $60\ \mu\text{L}$ total protein per sample, so we need $60\ \mu\text{L}$ stock solution: 4 samples (+1) x $40\ \mu\text{L}$ = $200\ \mu\text{L}$.

Prepare: $360\ \mu\text{L}$ at $10\ \mu\text{g}/\mu\text{L}$ = $3600\ \mu\text{g}$ beads in total ($1800\ \mu\text{g}$

hydrophilic beads + $1800\ \mu\text{g}$ hydrophobic beads).

Prepare the Carboxylated Magnetic Beads stock solution ($360\ \mu\text{L}$) with a concentration of $10\ \mu\text{g}/\mu\text{L}$.

22.1 Gently shake flasks with beads to resuspend bottom layer.

22.2 

Mix hydrophilic and hydrophobic beads at a 1:1 ratio: $36\ \mu\text{L}$ hydrophilic and $36\ \mu\text{L}$ hydrophobic beads (conc: $50\ \mu\text{g}/\mu\text{L}$).

22.3 Put beads on rack and remove liquid.

22.4 

Wash the beads 3 times with water at a concentration of $50\ \mu\text{g}/\mu\text{L}$:

$3600\ \mu\text{g}$ beads

at $5\ \mu\text{g}/\mu\text{L}$ = $720\ \mu\text{L}$ total volume.

22.5 Resuspend the beads in water at a working concentration of $10\ \mu\text{g}/\mu\text{L}$:

$3600\ \mu\text{g}$ beads

at $10\ \mu\text{g}/\mu\text{L}$ = $360\ \mu\text{L}$ water total.

23 Protein binding to beads.

23.1 

Add 100% EtOH to the sample to reach 80% EtOH (v/v).

- LCM samples: add $80\ \mu\text{L}$
- Cryosection samples: add $400\ \mu\text{L}$

23.2

Add magnetic beads stock to samples (stock concentration: **10 µg/µL** ,
100 µg or **600 µg** beads)

- LCM samples: add **10 µL**
- Cryosection samples: **60 µL**

23.3

45m

Incubate for **00:45:00** at **Room temperature** and **800 rpm** .

23.4 Spin down quickly.

24

Wash beads.

24.1 Insert the tubes into magnetic stand to collect the magnetic beads.

24.2 Remove supernatant and discard.

24.3

Add 80% EtOH (2 times the sample volume):

- LCM samples: use **200 µL**
- Cryosection samples: use **960 µL**

24.4

3m

Shake **00:03:00** at **Room temperature** , **800 rpm** .

24.5 

Collect beads in-between washed on magnetic stand, discard wash solution.

24.6 Repeat in total 3x.

Trypsin Digestion

25 

1. Make a trypsin stock solution in **50 millimolar (mM)** ABC for LCM and cryosection samples separately:

- **LCM samples:** final concentration = **0.005 µg/µL** , final volume = $15(+1) * 20 \text{ µL} = 320 \text{ µL}$

Mix **304 µL** ABC (**50 millimolar (mM)**) with **16 µL** trypsin (at **0.1 µg/µL** in **500 millimolar (mM)** ABC).

- **Cryosection samples:** final concentration = **0.025 µg/µL** , final volume = $4 * 60 \text{ µL} (+ 10 \text{ µL}) = 250 \text{ µL}$

Mix **187.5 µL** ABC (**50 millimolar (mM)**) with **62.5 µL** trypsin (at **0.1 µg/µL** in **500 millimolar (mM)** ABC).

26 

Add trypsin stock solution to each sample:

- **LCM samples:** add **20 µL** to LCM samples (contains ca. **0.1 µg** trypsin).
- **Cryosection samples:** add **60 µL** to reference samples (contains **1.2 µg** trypsin).


27 

Digest overnight at **37 °C** .

Peptide Extraction



15m


28 Spin the tubes down quickly.

29 Place tubes onto magnetic stage and transfer supernatant into a new  **1.5 mL** Eppendorf tube.

30 

Add ddH₂O to magnetic beads

- LCM samples: add  **30 µL**
- Cryosection samples: add  **70 µL**



31 Sonicate for  **00:15:00** .

15m

32 Spin down quickly.

33 Combine supernatant in the same Eppendorf tube as in step 1 above.

34 Quench digestion by adding 5% TFA

- LCM samples: add  **10 µL**
- Cryosection samples: add  **26 µL**

35 Check pH.



Drying and resuspension

10m

36 Snap freeze samples and place open collection tube in a vacuum evaporator (until completely dry).

37 Freeze samples or resuspend.

38 For resuspension: Prepare iRTs in MS buffer at a ratio of 1:20.

39 Resuspend samples in  **10 µL** (LCM) or  **20 µL** (cryo) MS buffer + iRT.

40 Sonicate for  **00:10:00** at  **Room temperature** .

10m

41 Spin down quickly.

42 Transfer the whole sample into a Waters glass vial.