



Sep 09, 2020

# Easy Covid LAMP

ayten<sup>1</sup><sup>1</sup>Hayat

1

Works for me

[dx.doi.org/10.17504/protocols.io.bk6gkzbw](https://dx.doi.org/10.17504/protocols.io.bk6gkzbw)

ahmet

DOI

[dx.doi.org/10.17504/protocols.io.bk6gkzbw](https://dx.doi.org/10.17504/protocols.io.bk6gkzbw)

PROTOCOL CITATION

ayten 2020. Easy Covid LAMP. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bk6gkzbw>



LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Sep 09, 2020

LAST MODIFIED

Sep 09, 2020

PROTOCOL INTEGER ID

41896

- 1 Collect sample into 1.5ml eppendorf tube.
- 2 Insert sample transfer plastic into the sample tube, leave for 5 seconds, then insert into pre-filled 0.2ml master mix PCR tube, shake for 5 seconds, discard.
- 3 Place PCR tube into pre-heated heat block device and run . 🔥 **65 °C (Heat block is pre-set)**  
🕒 **00:29:00**
- 4 Check the color of the PCR tube: positive if yellow, negative if pink or orange.