



Upload image

Aug 11, 2020

# Competitive ELISA to study the inhibition of the SpA-binding to Ab-3 by Ab-2.

Angel A Justiz-Vaillant<sup>1</sup><sup>1</sup>University of the West Indies St. Augustine**1** Works for me [dx.doi.org/10.17504/protocols.io.bjkbkksn](https://dx.doi.org/10.17504/protocols.io.bjkbkksn)[University of the West Indies](#) [angel.vaillant@sta.uwi.edu](mailto:angel.vaillant@sta.uwi.edu)Angel Justiz-Vaillant  
University of the West Indies St. Augustine

## ABSTRACT

This immunoassay tested the ability of anti-anti-SpA (Ab-2) to inhibit or interfere with the binding of anti-anti-SpA (Ab-3) to the original antigen- staphylococcal protein A (SpA). It is a confirmation test of the functional capacity of Ab-2. In this ELISA, both Ab-2 and the antigen SpA compete for binding to Ab-3. When Ab-2 is present, it blocks the binding of HRP-labeled SpA to Ab-3, resulting in inhibition of reaction color development.

## DOI

[dx.doi.org/10.17504/protocols.io.bjkbkksn](https://dx.doi.org/10.17504/protocols.io.bjkbkksn)

## PROTOCOL CITATION

Angel A Justiz-Vaillant 2020. Competitive ELISA to study the inhibition of the SpA- binding to Ab-3 by Ab-2..  
**protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bjkbkksn>

## LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Aug 11, 2020

## LAST MODIFIED

Aug 11, 2020

## PROTOCOL INTEGER ID

40291

## GUIDELINES

Use appropriate positive and negative controls, as well as blanks.

## MATERIALS

NAME	CATALOG #	VENDOR
Substrate (TMB)	EZMADP-60K(kit)	Millipore Sigma
Stop Solution (0.3M HCl)	EZMADP-60K(kit)	Millipore Sigma
96-Well Microtiter <sup>®</sup> ; Microplates, Polystyrene, 280 <sup>μ</sup> m; L, Nonsterile, V-bottom	2605	Thermo Fisher
Pierce <sup>®</sup> ; Recombinant Protein A, Peroxidase Conjugated	32400	Thermo Fisher

## SAFETY WARNINGS

Do not use any reagent where damages to the packaging has occurred.

- 1 Prepare materials and ELISA buffer solutions and reagents.
- 2 Pipette 53 µl of purified Ab-3 (100 µg/ml) mixed with 2.8 ml of coating buffer into each well.
- 3 Incubate the microplate at 37°C for 4 hrs.
- 4 Aspirate the contents of the wells.
- 5 Fill each well with an appropriately diluted washing solution and aspirate.
- 6 Wash the microplate 3 times.
- 7 Pipette 53 µl of serial dilutions of pooled Ab-2 (100 µg/ml) in triplicates.
- 8 Incubate the microplate at RT for 1 hr.
- 9 Rewash the microplate filling each well with 100 µl of washing buffer.
- 10 Pipette 53 µl of commercially available peroxidase-labeled SpA conjugate (Sigma-Aldrich) diluted 1:5000 to each well.
- 11 Reincubate the microplate at RT for 1 hr.
- 12 Repeat step 5.
- 13 Pipette 53 µl of TMB (Sigma-Aldrich) to each well.

14 Incubate in the dark at RT for 14 min.

15 Measure absorbance at 450 nm in a microplate reader and analyze the results