



APR 11, 2023

🌐 Untargeted IMS Tentative Identification Lipidomics

Lukasz Migas¹, Madeline E. Colley², Katerina V Djambazova², Martin Dufresne², Angela R.S. Kruse², David Anderson², Olof Isberg², Jamie Allen², ali.zahraei², Melissa Farrow², Jeff Spraggins², Raf Van De Plas²

¹Delft University of Technology; ²Vanderbilt University

VU Biomolecular Multimodal Imaging Center / Spraggins Research Group

Human BioMolecular Atlas Program (HuBMAP) Method Development



Angela R.S. Kruse
Vanderbilt University

OPEN ACCESS

DOI:
dx.doi.org/10.17504/protocols.io.4r3l27j7qg1y/v1

Protocol Citation: Lukasz Migas, Madeline E. Colley, Katerina V Djambazova, Martin Dufresne, Angela R.S. Kruse, David Anderson, Olof Isberg, Jamie Allen, ali.zahraei, Melissa Farrow, Jeff Spraggins, Raf Van De Plas 2023. Untargeted IMS Tentative Identification Lipidomics. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.4r3l27j7qg1y/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Apr 07, 2023

Last Modified: Apr 11, 2023

PROTOCOL integer ID:
80191

ABSTRACT

The purpose of this protocol is to generate tentative annotations for lipids detected using IMS.

- 1 Following pre-processing, tentative identification is performed using an in-house developed annotation software - annotine.
- 2 Generate an average mass spectrum of the dataset (in profile mode).
- 3 Scale the mass spectrum between 0 and 1, and peak-pick its profile to retrieve a list of m/z features (commonly 100s to 1000s).
- 4 Filter the peak list to retain only peaks that have a relative intensity above 0.001 (sensitive mode) or 0.01 (standard mode). Peaks whose intensity value falls below this threshold are removed.
- 5 (*optional*) De-isotope the peak list to remove M+1, M+2, ... and other potential isotopes from further consideration.

- 6 Generate an internal database on the basis of a user-supplied list of molecular species databases and a set of user-supplied expected adduct types:

Databases:

- coreMetabolome, LMSD, SHexCer, HMDB5
- (optional) a local LC-MS database

Adducts:

- Positive mode: $[M+H]^+$, $[M+Na]^+$, $[M+K]^+$
- Negative mode: $[M-H]^-$, $[M-CH_3]^-$

Note

- coreMetabolome5 is retrieved from metaspace (<https://metaspace2020.eu/help>);
- LMSD database is retrieved from LipidMaps (<https://www.lipidmaps.org/databases>);
- HMDB5 database is retrieved from HMDB (<https://hmdb.ca/downloads>);
- SHexCer is a manually curated list, and
- an LC-MS database can be (optionally) included if prior LC-MS experiments have been conducted.

- 7 Perform tentative identification by comparing the peak list with the built database of species and adduct combinations. If a peak is within a ± 5 ppm window of an annotation in the database, that annotation is assigned to that peak. This process is repeated until each peak has been compared to the database.
- 8 Evaluate each tentative annotation using metrics.
- 9 (optional) To reduce the number of unlikely annotations, calculate a false discovery rate (FDR) for every tentative identification. Annotations can be filtered based on the FDR (or any other) score.
- 10 (optional) If LC-MS results are available (see LC-MS/MS lipidomics protocol below), associate these directly with the annotine results. This allows immediate highlighting of which tentative identifications have also been observed by LC-MS/MS, increasing the confidence of the identification.

Protocol



NAME

Bulk Untargeted LC-MS/MS Lipidomics

CREATED BY

Madeline E. Colley

PREVIEW