

MAY 22, 2023

OPEN ACCESS

DOI:
dx.doi.org/10.17504/protocols.io.n92ldpx89l5b/v1

Protocol Citation: Masayasu Misu, Tomoki Yoshikawa, Satoko Sugimoto, Yuki Takamatsu, Takeshi Kurosu, Yukiteru Ouji, Masahide Yoshikawa, Masayuki Shimojima, Hideki Ebihara, Masayuki Saijo 2023. RCA-NGS for RNA viruses with ONT V14 chemistry.
protocols.io
<https://dx.doi.org/10.17504/protocols.io.n92ldpx89l5b/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
 We use this protocol and it's working

Created: Apr 10, 2023

Last Modified: May 22, 2023


PROTOCOL integer ID:
 80235

RCA-NGS for RNA viruses with ONT V14 chemistry

Forked from [RCA-NGS for RNA viruses](#)

Tomoki
 Masayasu Misu^{1,2}, Yoshikawa¹,
 Takeshi
 Yuki Takamatsu¹, Kurosu¹, Yukiteru Ouji²,
 Masahide Masayuki Hideki
 Yoshikawa², Shimojima¹, Ebihara¹,
 Masayuki
 Saijo¹

¹Department of Virology I, National Institute of Infectious Diseases;
²Department of Pathogen, Infection and Immunity, Nara Medical University



Tomoki Yoshikawa

ABSTRACT

Note that this version of the protocol was adopted to V14 chemistry of ONT.


This RCA-NGS were optimized for an NGS machine, MinION. These methods do not require nucleic acid amplification with virus-specific PCR primers, physical viral particle enrichment, and RACE.

These methods enable whole RNA viral genome sequencing by combining the following techniques:

















- 1) removal of unwanted DNA and RNA other than the RNA viral genome by nuclease treatment
- 2) the terminal of viral genome sequence determination by barcoded linkers ligation
- 3) Amplification of the viral genomic cDNA using an isothermal DNA amplification technique, such as rolling circle amplification (RCA).


This method can be exploited to determine any whole RNA viral genomes (i.e., single-stranded, double-stranded, positive-stranded, negative-stranded, non-segmented or multi-segmented genomes).

MATERIALS

-  Micrococcal Nuclease - 320,000 gel units New England Biolabs Catalog #M0247S

Keywords: Oxford Nanopore Technology, RNA virus, Sequence method, MinION, Nanopore sequencing, RCA-NGS

-  High Pure Viral RNA Kit Roche Catalog #11858882001
-  Turbo DNA-free Kit Invitrogen - Thermo Fisher Catalog #AM1907
- **NucleoSpin RNA Clean-up XS - Takara, Catalog #740903.10**
-  T4 RNA Ligase 2, truncated KQ - 2,000 units New England Biolabs Catalog #M0373S
- **The barcode-polyA linker DNA (e.g., The cSP6-polyA linker DNA)**
-  Superscript IV Thermo Fisher Scientific Catalog #18090050
- **SP6 primer (e.g., 5' phosphorylated SP6 primer)**
-  Deoxynucleotide (dNTP) Solution Mix New England Biolabs Catalog #N0447S
-  Superase-In RNase Inhibitor Thermofisher Catalog #AM2694
- **Dr.GentLE Precipitation Carrier - Takara Catalog #9094**
-  RNase H - 250 units New England Biolabs Catalog #M0297S
-  Agencourt AMPure XP Beckman Coulter Catalog #A63880
- **CirLigase II ssDNA Ligase - Biosearch Technologies Catalog #CL9021K**
- **GenomiPhi V3 Ready-To-Go DNA Amplification Kit - Cytiva Catalog #25-6601-24**
-  T7 Endonuclease I - 250 units New England Biolabs Catalog #M0302S
-  NEBNext FFPE DNA Repair Mix - 24 rxns New England Biolabs Catalog #M6630S
-  NEBNext Ultra II End Repair/dA-Tailing Module - 24 rxns New England Biolabs Catalog #E7546S
-  Blunt/TA Ligase Master Mix - 50 rxns New England Biolabs Catalog #M0367S
-  NEBNext Quick Ligation Module - 20 rxns New England Biolabs Catalog #E6056S
- **Native Barcoding Kit 24 V14 - Oxford Nanopore Technologies Catalog #SQK-NBD114.24**
-  Qubit 4 Fluorometer Thermo Fisher Scientific Catalog #Q33238
-  Qubit 1X dsDNA HS Assay Kit Thermo Fisher Scientific Catalog #Q33230
-  DNA LoBind Tube 1.5ml Eppendorf Catalog #022431021

-  0.2 ml PCR Tube strips Eppendorf Catalog #0030124359
- 100 % ethanol
- 70 % ethanol
- TE(pH8.0)
- nuclease-free H₂O

SAFETY WARNINGS




- ! Follow your facility's regulations and biosafety practices.




BEFORE START INSTRUCTIONS

This method was only confirmed to work with the working stocks that contain isolated RNA viruses at least 3.0×10^5 TCID₅₀ per ml.

It is recommended to check no bacterial contamination(e.g., *Mycoplasma* spp.).


Preparation for virus supernatant

- 1 Centrifuge the working stock virus to remove debris. 10m
 6000 x g, Room temperature, 00:10:00
 - 2 Transfer  180 µL virus supernatant to a 1.5ml screw cap tube.
 - 3 Unwanted DNA and RNA mainly originating from the virus-infected cells are digested using
 Micrococcal Nuclease - 320,000 gel units New England Biolabs Catalog #M0247S
- 3.1 Total 201 µl reaction 1h

-  180 μL virus supernatant
-  20 μL 10X Micrococcal Nuclease Reaction Buffer
-  1 μL Micrococcal nuclease

Mix by pipetting and spin down.


 37 °C water bath

 01:00:00

The viral genomic RNA extraction

4 The viral genomic RNA extraction is performed using


 High Pure Viral RNA Kit Roche Catalog #11858882001 .

4.1 Add  400 μL of binding buffer (with  4 μL PolyA carrier RNA).

10m


Mix gently by ~5 times pipetting and flicking thoroughly the tube, and spin down.

 Room temperature


 00:10:00

4.2 Transfer the sample to a High Pure Filter Tube.


1m

 8000 x g, Room temperature,
00:01:00

Discard the flow-through liquid and Collection Tube, and insert the Filter Tube into a new Collection Tube.

4.3 Add  500 μL of inhibitor removal bo transfer the sample to a High Pure Filter Tube.


1m

 8000 x g, Room temperature,
00:01:00

Discard the flow-through liquid and Collection Tube, and insert the Filter Tube into a new Collection Tube.

4.4 Add  450 μL of wash buffer.

1m

 8000 x g, Room temperature,
00:01:00


Discard the flow-through liquid and Collection Tube, and insert the Filter Tube into a new Collection Tube.

4.5 Add  450 µL of wash buffer.

1m


 13000 x g, Room temperature,
00:01:00 and discard the flow-through liquid.

Discard the Collection Tube and insert the Filter Tube into a 1.5 ml tube -

 DNA LoBind Tube 1.5ml Eppendorf Catalog
#022431021

4.6 Add  50 µL Elution Buffer.

1m


 13000 x g, Room temperature,
00:01:00

Note

The eluted RNA can be stored at -80°C.




Remove unwanted DNA

5 Unwanted DNA mainly from the virus-infected cells in the RNA sample is digested using a

 Turbo DNA-free Kit Invitrogen - Thermo Fisher Catalog
#AM1907

5.1 Total 56 µl reaction


30m

-  50 µL the eluted RNA
-  5 µL 10X reaction buffer
-  1 µL DNase I

Mix gently by pipetting and spin down.


 37 °C  00:30:00


6 The viral RNA is purified using **NucleoSpin RNA Clean-up XS - Takara, Catalog #740903.10.**

6.1 Add an equal volume  56 μL of Buffer RCU and mix gently.


6.2 Transfer the sample to a NucleoSpin RNA XS Column.

1m


 11000 x g, Room temperature,
00:01:00

6.3 Wash the column by  400 μL Buffer RA3.


1m

 11000 x g, Room temperature,
00:01:00

Discard the flow-through liquid and Collection Tube, and insert the NucleoSpin RNA XS Column into a new Collection Tube.

6.4 Wash the column by  200 μL Buffer RA3.


2m

 11000 x g, Room temperature,
00:02:00


Discard the flow-through liquid and Collection Tube, and insert the NucleoSpin RNA XS Column into a Nuclease-free Collection Tube(1.5 ml).

6.5 Add  10 μL RNase-free H_2O .

1m

 11000 x g, Room temperature,
00:01:00


Transfer the sample to a 0.2 ml PCR tube -

 0.2 ml PCR Tube strips Eppendorf Catalog
#0030124359

cSP6-polyA Linker DNA ligation

7

The viral RNA is ligated with cSP6-polyA Linker DNA using

 T4 RNA Ligase 2, truncated KQ - 2,000 units New England Biolabs Catalog
#M0373S






The RNA is ligated to the 3' end with the barcoded(complementary sequence of SP6 (cSP6)) polyA linker DNA. It is able to identify the 3' terminal viral genome sequence. The PolyA sequence is required for reverse transcription for ONT kit (SQK-PBK004/ PCS109).

Note

The cSP6-polyA linker DNA (5'-5rApp-CTATAGTGTCACCTAAATCAAAAAAAAAAAAAAAAAAAAA-3ddC-3'), which is pre-adenylated at the 5' terminal (5rApp), and consists of the complementary sequence of SP6 (CTATAGTGTCACCTAAATC), oligo (dA) 20, and dideoxycytidine (3ddC) at the 3' terminal, was synthesised for 3' linker ligation by Integrated DNA Technologies (Coralville, IA).

7.1 Total 20 µl reaction

15m

-  10 µL Purified RNA
-  1 µL 10 µM the cSP6-polyA linker DNA
-  2 µL 10X T4 RNA Ligase Reaction Buffer
-  6 µL 50% PEG8000 solution
-  1 µL T4 RNA Ligase 2, truncated KQ


Mix gently by pipetting and spin down.

Incubate  25 °C  00:15:00

8 The linker-ligated viral RNA is purified using **NucleoSpin RNA Clean-up XS - Takara, Catalog #740903.10**




Fill the sample to 100 µl with 80 µl TE (pH 8.0) and add 100 µl (equal volume) of Buffer RCU.

Eluted the RNA in  10 µL RNase-free H₂O and transfer the sample to a 0.2 ml PCR tube.

Reverse transcription

9 The viral RNA is reverse transcribed using

 Superscript IV Thermo Fisher Scientific Catalog #18090050

5' phosphorylated SP6 primer is used for reverse transcription.

Note

SP6 primer (5' phosphorylated SP6 primer); 5' [Phos]GATTTAGGTGACACTATAG 3'
5' phosphorylation is due to circularization.

9.1 Set up pre-mixture

6m

- 10 μ L RNA (~ 50ng)
- 1 μ L 50 μ M SP6 primer
- 1 μ L nuclease-free H₂O
- 1 μ L 10mM dNTP -



Deoxynucleotide (dNTP) Solution Mix New England Biolabs Catalog #N0447S

Mix gently by flicking the tube, and spin down.

65 °C

00:05:00

and

4 °C

00:01:00

9.2 Total 20 μ L reaction

20m

- 13 μ L pre-mixture sample
- 4 μ L 5X SSIV Buffer
- 1 μ L 100mM DTT
- 1 μ L RNase OUT -



Suprase-In RNase Inhibitor Thermofisher Catalog #AM2694

- 1 μ L SuperScript IV Reverse Transcriptase

Mix gently by flicking the tube, and spin down.

55 °C

00:10:00

80 °C

00:10:00


RNase H treatment

20m

10

Add  1 μL  RNase H - 250 units New England Biolabs Catalog
#M0297S


20m

 37 °C 00:20:00

cDNA purification using AMPure XP

11

cDNA is purified using

 Agencourt AMPure XP Beckman Coulter Catalog
#A63880

Prepare AMPure XP beads for use; resuspend by vortexing.


Transfer amplified DNA sample to 1.5ml low binding tube.

11.1

Add  36 μL (X1.8
volume)


AMPure XP reagent and mix by pipetting.

Incubate on rotor mixer.


 00:05:00 Room temperature

11.2

Spin down and pellet on a magnet.

Wait for  00:01:00 and pipette off the supernatant.

11.3

Wash twice by  100 μL 70 % ethanol and remove the ethanol using a pipette and discard.


11.4

Spin down and pipette off any residual ethanol.

11.5

Resuspend pellet in  12 μL TE(pH 8.0).

Incubate on a rotor mixer.

 00:05:00 Room temperature

5m


11.6

Spin down and pellet the beads on the magnet until the elute is clear and colourless.

11.7 Remove and retain  12 μL elute into a new tube.

(Optional step) Short cDNA fragment removal

12 Short cDNA fragment is removed from the viral RNA sample using

 Agencourt AMPure XP Beckman Coulter Catalog
#A63880


Prepare AMPure XP beads for use; resuspend by vortexing.
Transfer amplified DNA sample to 1.5ml low binding tube.

Note

If a significant proportion of the reads obtained from an NGS run fail to match with the NCBI-nr database (i.e., no hits), it could indicate a large number of short cDNA fragments in the sample. In such instances, re-performing the optional step following step 11 by adding x1.8 volume of AMPure XP could significantly enhance the outcomes.

12.1 Add  8 μL of TE to the  12 μL of elute to adjust  20 μL .


13 Size selection of the cDNA sample is performed using

 Agencourt AMPure XP Beckman Coulter Catalog
#A63880

X0.8 volume of AMPure beads recovers more than 200 bp of nucleic acids.


13.1 Add  16 μL (X0.8 volume) AMPure beads and mix by pipetting.


Incubate on rotor mixer.

 00:05:00

 Room temperature

13.2 Spin down and pellet on a magnet.


Wait for  00:01:00 and pipette off the supernatant.


13.3 Wash twice by  100 µL 70 % ethanol and remove the ethanol using a pipette and discard.

13.4 Spin down and pipette off any residual ethanol.

13.5 Resuspend pellet in  12 µL nuclease-free water.

Incubate on a rotor mixer.

 00:05:00

 Room temperature

5m

13.6 Spin down and pellet the beads on the magnet until the elute is clear and colourless.

13.7 Remove retain  12 µL elute into a new tube.






Circularization of cDNA

1h 10m

14 The cDNA is circularized using **CircLigase II ssDNA Ligase - Biosearch Technologies Catalog #CL9021K**.

14.1 Total 20 µl reaction

1h 10m

-  12 µL cDNA
-  2 µL 10X reaction buffer
-  1 µL 50 mM MnCl₂
-  4 µL 5M Betaine
-  1 µL CircLigase II

Mix by pipetting and spin down.

60 °C

01:00:00

80 °C

00:10:00

15

cDNA is purified using



Agencourt AMPure XP Beckman Coulter Catalog #A63880

go to step #11

(Add



36 µL (X1.8 volume)

AMPure beads)

Elute the pellet in



10 µL

nuclease-free H₂O.

Amplification of cDNA by rolling circle amplification (RCA)

16

cDNA is amplified by Rolling circle amplification (RCA) using **GenomiPhi V3 Ready-To-Go DNA Amplification Kit - Cytiva Catalog #25-6601-24.**

16.1

Total 20 µL reaction

3m

- 10 µL cDNA
- 10 µL 2X denaturation buffer

Mix by pipetting and spin down.

95 °C

00:03:00

4 °C on ice

16.2

Add 20 µL denatured sample to Ready to go GenomiPhi cake.

4h 10m

30 °C

04:00:00

65 °C

00:10:00

17

The cDNA is purified using



Agencourt AMPure XP Beckman Coulter Catalog #A63880

go to step #11

(Add



36 µL (X1.8 volume)

AMPure beads)

Elute the pellet in



40 µL



nuclease-free H₂O.

18

DNA concentration is measured using a Qubit 4 Fluorometer with



Qubit 1X dsDNA HS Assay Kit Thermo Fisher Scientific Catalog
#Q33230

-  199 μL 1X working solution
-  1 μL DNA

Mix by vortexing.

Incubate  00:02:00  Room temperature and measure.

Note

Confirm the total amplified cDNA to be over 1500 ng, as confirmed using, for instance, a Qubit 4 Fluorometer and Qubit 1X dsDNA HS Assay Kit.

T7 endonuclease treatment

19 The amplified cDNA by RCA is digested using



T7 Endonuclease I - 250 units New England Biolabs Catalog
#M0302S





to remove

branching.

The following protocol is modified based on the Native barcoding amplicons (with EXP-NBD104, EXPNBD114, and SQK-LSK109) protocol (NBA_9093_v109_revA_12Nov2019) provided by Oxford Nanopore Technologies website.


19.1 Total 30 μL reaction

30m

-  x μL (1.0 μg) DNA
-  3 μL NEBuffer 2
-  1.5 μL T7 endonuclease I
-  25-x μL nuclease-free H_2O

Mix by pipetting and spin down.

 37 $^{\circ}\text{C}$

 00:30:00

20

The cDNA is purified using



Agencourt AMPure XP Beckman Coulter Catalog
#A63880

➡ go to step #11

(Add

30 µL (X1.0
volume)

AMPure beads)

Resuspend pellet in 13 µL nuclease-free H₂O.

DNA repair and end-prep

21 The purified cDNA is end-prepped using



NEBNext FFPE DNA Repair Mix - 24 rxns New England Biolabs Catalog
#M6630S

and



NEBNext Ultra II End Repair/dA-Tailing Module - 24 rxns New England Biolabs Catalog
#E7546S

21.1 Total 15 µl reaction

35m

- 12 µL DNA
- 0.875 µL NEB Next FFPE DNA repair buffer
- 0.5 µL NEB Next FFPE DNA repair Mix
- 0.875 µL Ultra II end-prep reaction buffer
- 0.75 µL Ultra II end-prep reaction Mix

Mix by pipetting and spin down.

20 °C 00:30:00

65 °C 00:05:00

22 The cDNA is purified using



Agencourt AMPure XP Beckman Coulter Catalog
#A63880

➡ go to step #11

(Add

15 µL (X1.0
volume)

AMPure beads)

Resuspend pellet in 10 µL nuclease-free H₂O.

23 DNA concentration is measured using a Qubit 4 Fluorometer with

2m



Qubit 1X dsDNA HS Assay Kit Thermo Fisher Scientific Catalog
#Q33230

- 199 µL 1X working solution
- 1 µL DNA

Mix by vortexing.

Incubate  00:02:00  Room temperature and measure.

Note


Confirm the purified cDNA to be approximately 700 ng or more using, for instance, Qubit 4 Fluorometer with a Qubit 1X dsDNA HS Assay Kit.

Note

The cDNA can be stored at 4°C overnight.





Native barcode ligation

- 24** The end-prepped cDNA is ligated with native barcode using **Native Barcoding Kit V14 - Oxford Nanopore Technologies Catalog #SQK-NBD114.24** and

 Blunt/TA Ligase Master Mix - 50 rxns New England Biolabs Catalog #M0367S

- 24.1** Total 20 µl reaction

20m

-  x µL DNA (about 400ng)
-  1.5 µL native barcode
-  10 µL Blunt/TA ligase master mix
-  8.5-x µL nuclease-free H₂O


Mix by pipetting and spin down.

 25 °C  00:20:00

- 25** Add  20 µL TE (pH 8.0).


26

The cDNA is purified using

 Agencourt AMPure XP Beckman Coulter Catalog #A63880

 go to step #11

(Add

 20 μL (X1.0 volume)

AMPure beads)



Resuspend pellet in  12 μL nuclease-free H_2O .

27

DNA concentration is measured using a Qubit 4 Fluorometer with



Qubit 1X dsDNA HS Assay Kit Thermo Fisher Scientific Catalog #Q33230

-  199 μL 1X working solution
-  1 μL DNA

Mix by vortexing.

Incubate  00:02:00  Room temperature and measure.

Convert nanogram(ng) into femtomole(fmol) by a calculator.

Note

The molar concentration of the cDNA sample can be converted based on the length of the major band confirmed by electrophoresis after T7 endonuclease treatment. Typically, the fragment lengths are around 2000 bases pairs.

Adaptor ligation

20m

28

Pool each barcoded sample into a 0.2ml PCR tube (Total 100–200 fmol).

Note

Even when the V14 kit was used, we have good results with using total 100 to 200 fmol of samples.






29

Adaptor Ligation with pooled samples is performed using
Ligation Sequencing Kit - Oxford Nanopore Technologies Catalog #SQK-NBD114.24 and




29.1 Total 20 µl reaction

20m

-  x µL DNA (100-200 fmol)
-  2 µL Native Adapter (NA)
-  4 µL NEB Next Quick Ligation Reaction Buffer(5X)
-  2 µL Quick T4 DNA ligase
-  12-x µL nuclease-free H₂O

mix gently and incubate.

 25 °C

 00:20:00

30 The adaptor-ligated cDNA is purified using



Agencourt AMPure XP Beckman Coulter Catalog
#A63880


Prepare AMPure XP beads for use; resuspend by vortexing.

Transfer amplified DNA sample to 1.5ml low binding tube.

30.1

Add  10 µL (X0.5 volume) AMPure XP reagent and mix by pipetting.

Incubate on a rotor mixer.


 00:05:00

 Room temperature

30.2




Spin down and pellet on a magnet. Wait for  00:01:00 and pipette off the supernatant.

30.3

- Wash twice by  100 µL Short Fragment Buffer(SFB) and remove the SFB using a pipette and discard.


30.4



Spin down and pipette off any residual SFB.

30.5 ■ Resuspend pellet in  13 µL Elution Buffer (EB)
 00:05:00  Room temperature and tapping occasionally.

30.6 Spin down and pellet the beads on the magnet until the elute is clear and colourless.

30.7 Remove retain  13 µL elute into a new tube.

31 DNA concentration is measured using a Qubit 4 Fluorometer with
 Qubit 1X dsDNA HS Assay Kit Thermo Fisher Scientific Catalog
#Q33230

-  199 µL 1X working solution
-  1 µL DNA

Mix by vortexing.

Incubate  00:02:00  Room temperature and measure.

Sequencing by MinION

32 Sequencing according to the manufacturer's instructions.