

Aug 22, 2024 Version 2

# Immunofluorescence protocol for FFPE human post-mortem brain sections to detect alpha-synuclein and tau pathology V.2

DOI

[dx.doi.org/10.17504/protocols.io.n92ld8nxxv5b/v2](https://dx.doi.org/10.17504/protocols.io.n92ld8nxxv5b/v2)

Felicia Suteja<sup>1</sup>, Hongyun Li<sup>1</sup>, YuHong Fu<sup>1</sup>

<sup>1</sup>University of Sydney

ASAP Collaborative Rese...

Team Kirik



courtney.wright Wright

University of Sydney

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.n92ld8nxxv5b/v2](https://dx.doi.org/10.17504/protocols.io.n92ld8nxxv5b/v2)

**Protocol Citation:** Felicia Suteja, Hongyun Li, YuHong Fu 2024. Immunofluorescence protocol for FFPE human post-mortem brain sections to detect alpha-synuclein and tau pathology. protocols.io <https://dx.doi.org/10.17504/protocols.io.n92ld8nxxv5b/v2> Version created by [courtney.wright Wright](#)

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** August 18, 2024

**Last Modified:** August 22, 2024

**Protocol Integer ID:** 106180

**Keywords:** ASAPCRN, human, post-mortem, formalin-fixed, paraffin embedded, alpha-synuclein, tau, immunofluorescence

**Funders Acknowledgement:**

**Michael J Fox Foundation**

**Grant ID:** ASAP-000497



## Disclaimer

The **protocols.io** team notes that research involving animals and humans must be conducted according to internationally-accepted standards and should always have prior approval from an Institutional Ethics Committee or Board.

## Abstract

This protocol details the method for preparing and staining human formalin-fixed, paraffin embedded post mortem brain tissue to detect alpha-synuclein and tau pathology using immunofluorescence.

## Materials

### Materials

- HistoChoice® Clearing agent - *Sigma* #H2779
- Series 2 adhesive microscope slides - *Trajan*, #472042491
- Sodium borohydride - *Sigma* #71320, Lot #STBJ4218
- Human FC blocker - *BD Pharmingen* Human BD Fc Block, clone Fc1, #564220
- TrueBlack® Background Suppressor- *Biotium* #23012A
- TrueBlack® Lipofuscin Autofluorescence Quencher - *Cell Signaling Technology* #23007)
- Biotium's CoverGrip™ Coverslip Sealant - *Biotium* #23005)
- EverBrite™ Mounting Medium - *Biotium* #23001
- 100%, 95%, 70% ethanol
- 70% formic acid
- 0.1 M sodium citrate buffer pH 6 (house made)
- 0.1 M PBS (house made)

### Antibodies

A	B	C
SNCA (Gt)	R&D, AF1338	1:200
Ubq (Rb)	Abcam, ab7780	1:100
P62 (Ms IgG1)	BD #610833	1:100

#### Midbrain and cortical primary antibody details

A	B	C
Dn anti Gt-AF800	Thermo, #A32930	1:250
Dn anti Rb-AF647	Thermo, #A31573	1:200
Dn anti Ms-AF488	Thermo, #A21202	1:200
Hoescht-405	Sigma, #B2261	1:500

#### Midbrain and cortical secondary antibody details

A	B	C	D
Midbrain conjugated antibody			
AT8(Ms, IgG1)	Thermo, #M N1020	Zenon™ Ms IgG1 Labeling Kits (AF594)	Thermo, #Z2 5007
TH (Rb)	Merck, AB152	Zenon™ Rb Labeling Kits (AF532)	Thermo, #Z2 5303

Zenon™ midbrain conjugated antibody details

A	B	C	D
Cortical conjugated antibody			
AT8(Ms, IgG1)	Thermo, #M N1020	Zenon™ Ms IgG1 Labeling Kits (AF594)	Thermo, #Z2 5007
HuD (Rb)	Abcam ab302514	Zenon™ Rb Labeling Kits (AF532)	Thermo, #Z2 5303

Zenon™ cortical conjugated antibody details

## Hardware

- Leica HM325 rotary microtome
- Aptum Bio Retriever 2100, Aptum Biologics Ltd, UK

## Safety warnings

- ⚠ For hazard information and safety warnings, please refer to the SDS (Safety Data Sheet) for each of the raw materials used.



## Before start

Formalin-Fixed Paraffin-Embedded (FFPE) sections of the human midbrain and cortical regions were cut at 6µm with a Leica HM325 rotary microtome and mounted on Series 2 adhesive microscope slides.



## Day I - Tissue Prep

- 1 1. Place slides with FFPE human brain tissues on a slide rack and bake in a 60 °C oven for 01:00:00

1h

### De-paraffinize

- 2 2. Continue using the slide rack, submerge slides in the following solution to de-paraffinize:

19m

a. HistoChoice: 2 x 00:07:00

b. 100% ethanol: 2 x 00:03:00

c. 95% ethanol: 00:03:00

d. 70 % ethanol: 00:03:00

3. Submerge in MiliQ/ultrapure water for 00:03:00 .

4. Tap gently sideways on flat surface covered with paper towels to remove excess water.

### Antigen retrieval

- 3 5. Submerge slides in 70% Formic Acid for 00:20:00 .

2h 31m

6. Wash in MiliQ/ultrapure water for 2x 00:05:00 .

7. Remove slides from the rack and submerge slides directly into 0.01 Molarity (M) sodium citrate buffer ( 6 ) and incubate sections in a programmable antigen retrieval cooker

8. Let the pressure cooker reach its peak of 121 °C before gradually cooling for a total of 02:00:00 .

9. Wash with MiliQ/ultrapure water for 00:01:00 , followed by washing with 1 X PBS for 00:05:00 .

### Quenching aldehyde group



- 4 10. Prepare 0.1% Sodium borohydride ( $\text{NaBH}_4$ ) in 1xPBS for ⌚ 00:30:00 of quenching. The solution must always remain chilled in ice.

35m



11. Wash in 1x PBS for 2x ⌚ 00:05:00

## Human Fc blocker treatment

- 5 12. Add Human FC blocker in 1X PBS with a ratio of 1:50, incubate slides with Human FC blocker at 🌡 Room temperature for ⌚ 00:05:00

5m

## Background Suppressor system treatment

- 6 Component A (Background Suppressor) may become turbid or form a gel at 🌡 4 °C this does not affect performance. Warm the buffer to room temperature or 🌡 37 °C until clear (light blue) and completely liquid before use.

1h 30m



13. Add enough TrueBlack® Background Suppressor to completely cover sample.

🌡 Room temperature ⌚ 00:30:00

14. Remove the Background Suppressor and add IF Blocking Buffer (home-made).

🌡 Room temperature ⌚ 01:00:00

## Primary antibody incubation

- 7 5. Prepare 150 ul per sample of primary antibody solution consisting of selected primary antibody diluted in home-made blocking buffer.


2d

A	B
SNCA (Gt)	1:200
Ubq (Rb)	1:100
P62 (Ms IgG1)	1:100




## Primary antibody table and dilutions

16. Prepare humidified chamber and remove Blocking buffer by tapping on paper towel.

17. Add primary antibody diluted in Blocking buffer and incubate for  48:00:00 in

 4 °C

**Day 4 - Secondary antibodies**


8 1. Wash slides in PBST 3x  00:10:00 . Slides must be protected from light from this step.

2h 20m

2. Prepare 150 ul per sample of secondary antibody solution consisting of selected secondary antibody diluted appropriately in home-made blocking buffer.


A	B
Dn anti Gt-AF800	1:250
Dn anti Rb-AF647	1:200
Dn anti Ms-AF488	1:200
Hoescht-405	1:500

## Midbrain and cortical secondary antibodies and dilutions

3. Incubate sample in secondary antibody in  Room temperature for  02:00:00

4. Wash slides in PBST 3x  00:05:00 , followed by 1x PBS 1x  00:05:00

**Conjugated antibodies**

9 5. Prepare conjugated antibodies using Zenon™ labeling kits within  00:30:00 before applying the antibodies.

2h 35m

A	B
AT8(Ms, IgG1)	Zenon <sup>®</sup> Mouse IgG1 Labeling Kits (AF594)



A	B
TH (Rb)	Zenon <sup>®</sup> Rb Labeling Kits (AF532)

Midbrain conjugated antibodies

A	B
AT8(Ms, IgG1)	Zenon <sup>®</sup> Ms IgG1 Labeling Kits (AF594)
HuD (Rb)	Zenon <sup>®</sup> Rb Labeling Kits (AF532)

Cortical conjugated antibodies

6. Incubate sample in conjugated antibodies in 🌡 Room temperature for 🕒 02:00:00

7. Wash slides in 1x PBS 3x 🕒 00:05:00

## Post-treatment with TrueBlack<sup>®</sup> Lipofuscin Autofluorescence Quencher

10 8. Prepare 150 ul per sample of 1xTrueBlack<sup>®</sup> Lipofuscin Autofluorescence Quencher in 70% ethanol. 5m

9. Tap slides to paper towels to remove excess washing solution and then place in humidified chamber.


10. Add diluted TrueBlack<sup>®</sup> solution and incubate for 30-60 seconds.

11. Rinse slides in 1xPBS 3x 🕒 00:05:00


## Mounting and cover slipping

11 12. Remove as much moisture without drying tissue using Kimwipes



13. Retrieve Mounting with EverBrite™ Mounting Medium to warm at  Room temperature before dispensing approximately 10 - 20 ml on slides.

14. Place coverslip gently on the slides and wait for the Mounting Medium to dry before applying the perimeter of the coverslip with Biotium's CoverGrip™ Coverslip Sealant.

15. Store in a protected slide box away from light at  4 °C