



DEC 27, 2023

🌐 Fixative preparation for Baermann funnel extraction

Casey A Schlenker¹¹The Ohio State University

OSU Soybean Pathology and Nematology



Casey A Schlenker

ABSTRACT

Fixative preparation for Baermann funnel extraction

MATERIALS

Materials for fixative preparation:

- Para-formaldehyde powder (PFA)
- Phosphate-buffered saline (PBS) tablets
- 1M NaOH
- Glycerol

OPEN  ACCESS

Protocol Citation: Casey A Schlenker 2023. Fixative preparation for Baermann funnel extraction.

protocols.io

<https://protocols.io/view/fixative-preparation-for-baermann-funnel-extractio-c6t4zeqw>











License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited


Protocol status: Working
We use this protocol and it's working

Created: Dec 27, 2023**Last Modified:** Dec 27, 2023**PROTOCOL integer ID:**
92764


Fixative Preparation


50m

- 1 Weigh  1.6 g PFA powder
- 2 Dissolve in  20 mL of DI water
- 3 Incubate in water bath at  60 °C to help dissolve PFA powder for  00:40:00 40m
- 4 Add  20 µL of 1 M NaOH to further help dissolve PFA powder
- 5 Place back in  60 °C water bath for  00:10:00 10m
- 5.1 If still not dissolved, add  10 µL more of 1M NaOH
- 6 Pull out of water bath and cool to  Room temperature
- 7 Add two tablets of PBS to a  100 mL volumetric flask

8 Mix  20 mL of each solution to get 4% formaldehyde

9 Add  400 μ L of glycerol to fixative

10 Ensure  7.5

11 Place fixative in a hot water bath till  60 °C