



Aug 30, 2024

JAX-Sen: Mouse kidney dissociation for single-cell RNA sequencing

DOI

dx.doi.org/10.17504/protocols.io.j8nlk833wl5r/v1

Ramalakshmi Ramasamy¹, Juliana Alcoforado Diniz¹, Jessica Garofalo¹, Patrick Fleming¹, Paul Robson^{1,2}

¹The Jackson Laboratory for Genomic Medicine, Farmington, CT, USA;

²Department of Genetics and Genome Sciences, University of Connecticut School of Medicine, Farmington, CT, USA

Cellular Senescence Net...



Harshpreet Chandok

The Jackson Laboratory

OPEN  ACCESS



DOI: **dx.doi.org/10.17504/protocols.io.j8nlk833wl5r/v1**

Protocol Citation: Ramalakshmi Ramasamy, Juliana Alcoforado Diniz, Jessica Garofalo, Patrick Fleming, Paul Robson 2024. JAX-Sen: Mouse kidney dissociation for single-cell RNA sequencing. **protocols.io** **<https://dx.doi.org/10.17504/protocols.io.j8nlk833wl5r/v1>**

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: July 02, 2024

Last Modified: August 30, 2024

Protocol Integer ID: 102773

Keywords: JAX-Sen



Funders Acknowledgement:

**National Institute on Aging
(NIA) JAX-Sen Senescence
Tissue Mapping Center
Grant ID: ID: U54 AG079753**

Abstract

These samples are part of the JAX-Sen project in the SENNET Consortium. We aim to study and characterize senescence in the C57Bl/6 mouse kidney. We compare young (6 months old) and old (24 months old) mouse kidneys using scRNAsequencing. This protocol describes the single-cell dissociation of the kidneys before library preparation and sequencing.



Reagents and Materials

- 1
 - 5 ml Protein lobind tubes
 - 1.5 or 2 ml Protein lobind tubes
 - Micro scissors/ scalpel
 - Ice-cold PBS
 - Wide-bore pipette tips
 - GentleMACS dissociator
 - GentleMACS C tubes
 - DMEM + 10% FBS
 - ACK lysis buffer – 2 ml/sample
 - SS_buffer– 5 ml/sample
 - 40um cell strainer
 - 70uM strainer

SS_Buffer:

A	B	C	D	E	F	G	H	I
	Stock co nc.	Working conc.	50 ml	100 ml	250 ml	500 ml	20 ml	
EDTA	0.5 M	2 mM	0.2	0.4	1	2	0.08	ml
BSA	10%	2%	10	20	50	100	4	ml
FBS	100%	15%	7.5	15	37.5	75	3	ml
DPBS	1X	make up	32.2	64.4	161	322	12.88	ml

Procedure

- 2 Collect shipment from the warehouse.
- 2.1 Fill up ice buckets and prep them with reagents:
Bucket 1 (small): Enzyme mix, PBS, SS_buffer
Bucket 2 (large): For mincing tissues and sample tubes
- 2.2 Confirm if the sample is cold. Take pictures. Transfer sample to 5 ml tubes.
- 2.3 Prepare enzyme mix by adding 2.35 mL of serum-free RPMI 1640 or DMEM, 100 µL of Enzyme D, 50 µL of Enzyme R, and 12.5 µL of Enzyme A of the Multi Tissue Dissociation Kit 1 into a



gentleMACS C Tube for up to 0.5 g of tissue.

- 2.4 Mince mouse kidneys (2-4 mm³) with sterile micro scissors.
- 2.5 After cutting, transfer the sample into labeled 5 ml tubes.
- 2.6 Tightly close C Tube and attach it upside down onto the sleeve of the gentleMACS Dissociator.
- 2.7 Run the gentleMACS Program Multi_B. If using the heating function of the gentleMACS Octo Dissociator with Heaters run program 37C_Multi_E.
- 2.8 After termination of the program, detach C Tube from the gentleMACS Dissociator.
- 2.9 Resuspend sample and apply the cell suspension to a MACS® SmartStrainer (70 µm) placed on a 5 mL tube.
- 2.10 Wash MACS SmartStrainer (70 µm) with **4 mL of ice-cold PBS**.
- 2.11 Centrifuge at **300 x g, 4°C, 7 min** to get rid of the supernatant. Aspirate supernatant.
- 2.12 Resuspend the cell pellet in 2 mL **ACK Lysing buffer**, and incubate on **ice for 3 min** to exclude RBCs.
- 2.13 Add 5 mL **SS_buffer** to stop ACK Lysing buffer reaction, then collect cells by centrifugation **500 x g, 4°C, 5 min**.
- 2.14 Resuspend the pellet in **1ml of ice cold PBS**.
- 2.15 Filter the cell suspension through **Flowmi 40um cell strainers** tube to remove debris and non-dissociated tissue fragments.
- 2.16 Collect cells by centrifugation **500 x g, 4°C, 5 min** and proceed to Flex protocol.

