



Jul 03, 2021

THP-1 differentiation

Juliane C R Fernandes¹¹Faculdade de Medicina da Universidade de São Paulo

1 Works for me



Share

This protocol is published without a DOI.



Juliane C Fernandes

Faculdade de Medicina da Universidade de São Paulo, Institut...

DISCLAIMER

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

ABSTRACT

Differentiation of THP-1 monocytic cell line in macrophage-like cells

PROTOCOL CITATION

Juliane C R Fernandes 2021. THP-1 differentiation. **protocols.io**
<https://protocols.io/view/thp-1-differentiation-bvx7n7rn>



FUNDERS ACKNOWLEDGEMENT

FAPESP 2017/21906-9

KEYWORDS

THP-1, monocyte, macrophage, differentiation, PMA, cell line

LICENSE

— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jun 19, 2021

LAST MODIFIED

Jul 03, 2021

PROTOCOL INTEGER ID

50911

PARENT PROTOCOLS

In steps of

[Leishmania infectivity assessment using imaging flow cytometry](#)

GUIDELINES

Differentiation of THP-1 monocytic cell line in macrophage-like cells.

MATERIALS TEXT

THP-1 monocytic cell

line ATCC Catalog #TIB202

RPMI 1640 LGC

biosearch Catalog #BR30011-05

FBS Gibco,

ThermoFisher Catalog #12657-029

phorbol-12-myristate-13-acetate (PMA) Sigma-

aldrich Catalog #1002714643

Trypan Blue Solution, 0.4% Thermo

Fisher Catalog #15250061

Falcon Tube (50 mL) Fischer Scientific Step 1 PenStrep Invitrogen - Thermo Fisher

Neubauer Improved Haemocytometer
Counting Chamber
Counting Chamber

Hawksley AC1000 [Link](#)

Void depth: 0.1 mm
Counting Area: 1mm²



SAFETY WARNINGS

Recommended taking biosafety precautions.

DISCLAIMER:

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Plating and starting differentiation process

3d

10m

1

THP-1

monocytes ATCC Catalog #TIB-202

cultivated in RPMI 10% FBS Pen/Strep in a humidified

incubator at 34 °C 5% CO₂. Carefully check for any contamination before initiating the protocol.

THP-1

Centrifuge [monocytes ATCC Catalog #TIB-202](#)

at [200 x g, Room temperature](#), 00:08:00

using [Falcon Tube \(50 mL\) Fischer Scientific](#)

Centrifuge
Benchtop Centrifuge

Eppendorf 5405000441 [↗](#)

Any benchtop centrifuge will suffice



1.1 Estimate the number of bottles necessary to achieve enough cells for the experiment. After the first round of centrifuging, discard culture media and centrifuge the content of new bottles in the same tube. [↗ go to step #1](#)

1.2 Resuspend pelleted cells in a volume sufficient for counting (usually 5 mL for each centrifuged bottle)^{2m}

1.3 Mix [10 µl](#) of the resuspended cells with [10 µl](#) of Trypan Blue. 2m

2 Count and evaluate cell viability placing [10 µl](#) of the 1:1 cell:Trypan mixture in a Neubauer Chamber 10m
Calculate the concentration of cell per mL (Medium of cells/chamber × 2 × 10⁴)

2.1 Dilute cells/mL according to the preferred plate for the follow-up experiment

24 well plate - 1 × 10⁶ cells/well - [1 mL](#) /well (flow cytometry)

6 well plate - 5 × 10⁶ cells/well - [4 mL](#) /well (RNA extraction)

3 Add PMA.
Add PMA from the diluted stock (1 µg/mL). First, dilute 100x to obtain an aliquot of 10 ng/µL. Final concentration for THP-1 differentiation: 30 ng/mL diluted in RPMI 10% FBS.

3.1 Store for [72:00:00](#) at the incubator [34 °C 5% CO2](#) 3d

Resting 3d

4 Aspirate cell culture media with a vacuum pump. Wash cells 3 times with PBS 1X to remove residual PMA.

- 5 Incubate cells with 1 or 4 mL - respectively for 24 well and 6 well plates - of fresh RPMI 10% FBS medium for further ^{3d}
🕒 72:00:00 🌡️ 34 °C 5% CO2

THP-1 macrophages

- 6 Proceed to the assay of interest with fully differentiated THP-1 macrophages.

This protocol was tested to the following assays: miRNA inhibition, *Leishmania* infection, flow cytometry, and RNA extraction.

After differentiation cells stop duplicating and adhere to the bottom of plate wells (THP-1 monocytic cell line grow in suspension).