



JAN 30, 2024

Rotarod test

In 1 collection

mariangela.massarocenere^{1,2,3}

¹Department of Experimental Neuroscience, Santa Lucia Foundation IRCCS, Rome, Italy;

²Department of Systems Medicine, University of Rome Tor Vergata, Rome, Italy;

³Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, United States



n.berretta

ABSTRACT

Protocol to assess motor coordination and balance by using an accelerating rod

MATERIALS

Rotarod apparatus custom for rats

OPEN ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.kqdg3xo81g25/v1

Protocol Citation: mariangela.massarocenere 2024. Rotarod test . protocols.io
<https://dx.doi.org/10.17504/protocols.io.kqdg3xo81g25/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: Jan 25, 2024

Last Modified: Jan 30, 2024

- 1 Habituate the animal for 1 h to the testing room before the test
- 2 Place the rat on the elevated rod (47 cm from the floor) rotating at low speed (4 rpm) for at least 30 seconds before the first session of the test
- 3 Place the animal on an accelerating rod from 4 to 40 rpm in 300 s
- 4 Repeat step 3 for 3 times with at least 5 min of rest between trials
 - 4.1 Record the time the rat fall from the rod (maximum 300 seconds)
- 5 Clean the apparatus every time between animals and trials
- 6 Repeat steps 4-5 for 4 times in two consecutive days (with an inter-session interval of 3.5 h)