

MAR 18, 2024

🌐 Immunohistochemistry protocol optimized at IMM-JLA

Ana M Biscaia Santos¹, Ana M Cristóvão Pinto¹, Ana R Pires¹,
Joana G Antunes¹

¹Comparative Pathology Unit, Instituto de Medicina Molecular - João Lobo Antunes, Lisboa, Portugal

Ana M Biscaia Santos: Biomedical Scientist Specialist

Ana M Cristóvão Pinto: Biomedical Scientist

Ana R Pires: Biomedical Scientist Specialist

Joana G Antunes: Biomedical Scientist - corresponding author: joana.antunes@medicina.ulisboa.pt

OPEN  ACCESS



Joana G Antunes

Instituto de Medicina Molecular - João Lobo Antunes

DISCLAIMER

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](https://dx.doi.org/10.17504/protocols.io.j8nlkon2wv5r/v1) may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](https://dx.doi.org/10.17504/protocols.io.j8nlkon2wv5r/v1), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

DOI:

dx.doi.org/10.17504/protocols.io.j8nlkon2wv5r/v1

Protocol Citation: Ana M Biscaia Santos, Ana M Cristóvão Pinto, Ana R Pires, Joana G Antunes 2024.

Immunohistochemistry protocol optimized at IMM-JLA.

protocols.io

<https://dx.doi.org/10.17504/protocols.io.j8nlkon2wv5r/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: Oct 31, 2023

Last Modified: Mar 18, 2024

PROTOCOL integer ID: 90197

Keywords:

immunohistochemistry, histology,
immunostaining, FFPE, frozen,
Envision HRP

ABSTRACT

The immunohistochemistry (IHC) protocol is used to showcase the location of specific proteins in relation to the architecture of the tissue.

In this method, we use a primary antibody which binds to a specific antigen in the protein of interest. This antibody has a specific host and can be ready to use or undiluted. The secondary antibody is an anti-host IgG which will bind to the primary antibody ¹. It's also conjugated with an Horseradish Peroxidase (HRP).

The HRP will react with the Hydrogen Peroxide in the substrate of the DAB solution, oxidizing the DAB, creating a brown deposit on the location of the protein of interest.

Depending on location of the protein, we can have staining in the nuclei, the cytoplasm or the cellular membrane.

Possible variables are:

the concentration of the antibody;

type of antigen retrieval;

the pH of the antigen retrieval solution (for Heat Induced Epitope Retrieval (HIER));

the time of incubation;

the tissue used ².

You can use this protocol for Formalin-Fixed Paraffin-Embedded samples (FFPE) and Fixed Frozen samples (Gelatin and OCT).

MATERIALS

Coplin Jars

Incubation chamber

Slide Dipper

Slide dishes with Lids

Wash bottles with nozzle

Magnets

Magnetic mixer

Tweezers

Micropipettes p20, p200 and p1000 and respective tips

Hydrophobic PAP pen

Eppendorfs 1.5mL

Plastic pipettes

Equipment Pretreatment ModuleTM Deparaffinization and Heat-Induced Epitope Retrieval by Thermo Scientific

Reagents:

Ethanol 70%

Ethanol 96%

Ethanol 100%

Xylene

Distilled Water

PROTOCOL MATERIALS

- ⊗ Protein Block Serum free **Dako Catalog #Ref. X0909** Step 17
- ⊗ Hydrogen Peroxide Solution **Merck MilliporeSigma (Sigma-Aldrich) Catalog #7722-84-1**
- Step 13
- ⊗ Entellan ® **Merck MilliporeSigma (Sigma-Aldrich) Catalog #1.00869** Step 29
- ⊗ Envision™ Dako DAB **Dako Catalog #Ref. K3468** Step 22
- ⊗ Harris hematoxylin for histology **Bio-optica Catalog #05-06004E** Step 24
- ⊗ Antibody Diluent **Dako Catalog #Ref. 8006** Step 18
- ⊗ EnVision FLEX Target Retrieval Solution Low pH **Dako Catalog #K8005** Step 10
- ⊗ EnVision FLEX Target Retrieval Solution High pH **Dako Catalog #K8004** Step 10
- ⊗ Envision™ Flex Wash Buffer **Dako Catalog #Ref. K8007** In 5 steps
- ⊗ EnVision™ Flex conjugated w/ HRP **Dako Catalog #Ref. K4061** Step 20
- ⊗ EMSURE Methanol **Supelco Catalog #1.06009.2500** In 2 steps

SAFETY WARNINGS

! This protocol uses the following hazardous chemicals:

Xylene
Methanol
Hydrogen peroxide
DAB
Resinous Mounting Medium




Please comply with all safety measures and risk assessment protocols.

Latex gloves can be used for the entire protocol, except when manipulating the reagents previously stated, during which you should use nitrile gloves. Type IV biological disposable bin should be readily available.






BEFORE START INSTRUCTIONS

If FFPE - begin protocol after adhering sections to the slides for 1h at 60°C

Preparation of Wash Buffer

- 1 Measure  20 mL of 20X  Envision™ Flex Wash Buffer **Bio-optica Catalog #Ref. K8007**
- 2 Make up with distilled water to  1 L of total volume

Sample Preparation after sectioning - FFPE

- 3 Deparaffinize slides with Xylene for  00:10:00 10m
- 4 Change Xylene and deparaffinize for another  00:10:00 10m
- 5 Place slides in Ethanol 100% for  00:05:00 5m
- 6 Place slides in Ethanol 96% for  00:05:00 5m
- 7 Place slides in Ethanol 70% for  00:05:00 5m

8 To complete the hydration, place slides in distilled Water for 00:05:00

5m

Sample Preparation after sectioning - Frozen

9 **If frozen samples:**

10m

Remove from storage and let dry.

If embedded in OCT, put in distilled water for 00:05:00 at Room temperature .

If embedded in gelatin, put in PBS for 00:05:00 at 37 °C .

Antigen Retrieval

10 Perform antigen retrieval using PT Module by Thermo Scientific with a cycle of pre-heat to 65°C and heat to 95°C for 20 minutes.

Note

This equipment, when optimized, allows you to skip the step of deparaffinization in this section, by putting the slides directly through the antigen retrieval step.

This equipment allows two buffers with different pH, in which we use

EnVision FLEX Target Retrieval Solution High pH **Dako Catalog #K8004** and

EnVision FLEX Target Retrieval Solution Low pH **Dako Catalog #K8005** , depending on the optimization of the antibody.

Blocking Endogenous Peroxidase and Total Proteins

35m





11 Wash in 1X EnvisionTM Flex Wash Buffer **Bio-optica Catalog #Ref. K8007** 3 times for 00:05:00 each

5m

12 Put slides in 200 mL of 100% EMSURE Methanol **Bio-optica Catalog #1.06009.2500**

Note

Extemporaneous solution. 100% methanol can be used up to 3 times.



- 13 Put  6 mL of  Hydrogen Peroxide Solution **Bio-optica Catalog #7722-84-1** at  4 °C in the 100%  EMSURE Methanol **Bio-optica Catalog #1.06009.2500** with the slides

- 14 Cover the slide dish and incubate slides in the solution for  00:30:00

30m

Note

Hydrogen Peroxide must be added each time when immersing the slides, to obtain a concentration of 3%.




- 15 Wash in 1X  Envision™ Flex Wash Buffer **Bio-optica Catalog #Ref. K8007** 3 times for  00:05:00 each

- 16 Draw a limit on the edge of the tissue with the hydrophobic PAP Pen

- 17 Put  Protein Block Serum free **Bio-optica Catalog #Ref. X0909** on the tissue

Incubation and Washing



1h 54m 10s




- 18 Incubate with Primary antibody in the optimized concentration (if needed, dilute in  Antibody Diluent **Bio-optica Catalog #Ref. 8006**) for  01:00:00 at  Room temperature

1h

Note



Normally, around 100µL of diluted antibody per section, depending on the size of the tissue.




19 Wash in 1X  Envision™ Flex Wash Buffer **Bio-optica Catalog #Ref. K8007** 3 times for  00:05:00 each

20 Incubate with secondary antibody 30m
 EnVision™ Flex conjugated w/ HRP **Bio-optica Catalog #Ref. K4061** for  00:30:00
 Room temperature

Note


Be sure to match the secondary antibody to the host of the primary antibody (ex: if using mouse anti-CD3 primary, use anti-mouse secondary)

21 Wash in 1X  Envision™ Flex Wash Buffer **Bio-optica Catalog #Ref. K8007** 3 times for  00:05:00 each

22 Incubate with  Envision™ Dako DAB+ **Bio-optica Catalog #Ref. K3468** for  00:02:00 to 7m
 00:05:00

Safety information

After incubation, pour DAB excess on a paper towel and dispose on type IV waste, with tips and eppendorfs used, avoiding spills at any cost.
 Clean the incubation chamber and surfaces with Bleach.


23 Wash in distilled water for  00:05:00


5m








Counterstain and mount

17m 10s

- 24 Counterstain with  Harris hematoxylin for histology **Bio-optica Catalog #05-06004E** for 10s

 00:00:10

Note

filter Hematoxylin before use
- 25 Put slides in lukewarm running water for  00:05:00 5m
- 26 Dehydrate with 96% ethanol for  00:02:00 2m
- 27 Place in 100% ethanol for  00:02:00 2m
- 28 Clear in Xylene for  00:10:00 10m
- 29 Mount with appropriate mounting medium (i.e.  Entellan® **Bio-optica Catalog #1.00869**)

Note

Let the slides dry in a fume hood until complete polymerization of the mounting medium is achieved