

FEB 28, 2023

OPEN BACCESS

יוסם

dx.doi.org/10.17504/protocol s.io.36wgqj9zyvk5/v1

Protocol Citation: Jason D Limberis 2023. uDumBell - Circularization of rv0678 for genotypic bedaquiline resistance testing of Mycobacterium tuberculosis. protocols.io

https://dx.doi.org/10.17504/protocols.io.36wgqj9zyvk5/v1

License: This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working We use this protocol and it's working

Created: Sep 16, 2022

Last Modified: Feb 28, 2023

PROTOCOL integer ID:

70166

UDumBell - Circularization of rv0678 for genotypic bedaquiline resistance testing of Mycobacterium tuberculosis

Jason D Limberis¹

¹University of California, San Francisco



Jason D Limberis University of California, San Francisco

ABSTRACT

The ligation of dumbbell (hairpin) oligos to linear dsDNA produces pseudo-circular DNA. Including deoxyUridine in the PCR primer sequences causes Q5 and other high-fidelity polymerases to arrest elongation. This results in overhangs that were successfully ligated to a complementary hairpin structure. The deoxyUridine reduced the PCR product by approximately two-thirds, but this was ameliorated by increasing the Q5 DNA polymerase concentration three-fold.

ATTACHMENTS

222.png

MATERIALS

Rv0678 amplification primers (you can use any primer set here, and multiplex them, these are not yet optimized, the random sequence is for blunt-end cloning to detect *chimeras*)

A	В
Forward primer	/5Phos/GUCTATTTTCTGTTGGTG CTGATATTGC
Reverse primer	/5Phos/GUCTATACTTGCCTGTCG CTCTATCTTC

А	В
uDu mBel I	/5Phos/ATAGACCGAGACAGTAGAAGACCATGAACAAGCAGCACACGATA AACTAGACACCCTACTGTCTCG

Preferably PAGE purified

- Q5 Hot Start High-Fidelity DNA Polymerase 500 units **New England Biolabs Catalog #M0493L**
- X T4 DNA Ligase 20,000 units **New England Biolabs Catalog #M0202S**

Optinal

- 🔀 Exonuclease I (E.coli) 15,000 units New England Biolabs Catalog #M0293L
- 🔀 Exonuclease VIII truncated New England Biolabs Catalog #M0545S

Amplicon PCR

1

A	В
Component	Volume (ul)
5X Reaction Buffer	10
5X Q5 High GC Enhancer	10
10 mM dNTPs	1
Forward primer	2.5

A	В
Reverse primer	2.5
DNA (5ng)	2
Q5 High-Fidelity DNA Polymerase	1.5
Nuclease-Free Water	20.5

PCR using primer set

A	В	С	D
Step	Temp (C)	Time (s)	Cycles
Denaturation	98	30	1
Denaturation	98	10	
Annealing	62	10	34
Extension	72	20	
Extension	72	2	1

Cycle parameters

Adapter ligation

Prepare the dumbell (hairpin) by incubating at \$80 °C followed by cooling to room temperature over 00:30:00 (this only needs to be done once)

30m

3

A	В
Component	Volume (ul)
T4 DNA Ligase Buffer (10X)	2
PCR product (upto 1ug), as low as 50ng, probably much lower possible)	10
dumbell adapter	3
Ligase (add last, don't vortex)	1

A	В
H20	4

Incubate as below, with the lid temperature set to $$40 ^{\circ}C

A	В
Temp	Minutes
22	30
15	120
4	120
65	5

4 Incubate at 8 Room temperature for \$\infty\$ 00:05:00

Place on amagnetic rack

Aspirate supernatant

Add \perp 200 μ L [M] 70 % (V/V) ethanol

Wait for (5) 00:00:30

Aspirate and discard the supernatant

Add \perp 200 μ L [M] 70 % (v/v) ethanol

Wait for (5) 00:00:30

Aspirate and discard the supernatant

Resuspend beads in Z 20 µL of H20

Incubate for 👏 00:02:00

Transfer to a clean PCR tube

Exonuclease treatment - optional

1h

1h

8m

A B
Component Volume (ul)

NEBuffer 4 (10x) 1

Exonuclease VIII (truncated) 1

DNA 18

Incubate at 4 37 °C for 6 00:30:00

Stop reaction by adding EDTA to at least 11 mM.

protocols.io |

5

Heat Inactivation 4 70 °C C for 00:30:00