

VERSION 1

DEC 11, 2023

OPEN ACCESS



Protocol Citation: Satheeswaran Balasubramanian, Ekambaram Perumal 2023. Cartilage staining. **protocols.io** <https://protocols.io/view/cartilage-staining-c5hxy37n>Version created by Satheeswaran Balasubramanian

MANUSCRIPT CITATION: Balasubramanian S, Rangasamy S, Vivekanandam R, Perumal E. Acute exposure to tenorite nanoparticles induces phenotypic and behavior alterations in zebrafish larvae. Chemosphere. 2023 ;339:139681. <https://doi.org/10.1016/j.chemosphere.2023.139681>

🌐 Cartilage staining V.1

Satheeswaran Balasubramanian¹, Ekambaram Perumal¹

¹Bharathiar University

Fish behavior and physiology



Satheeswaran Balasubramanian
Bharathiar University

ABSTRACT

The Alcian blue staining technique is widely used among developmental biologists to observe the development of cartilage and bone in zebrafish embryos and larvae. Alcian blue is a positively charged dye that stains the cartilage through an electrostatic interaction with negatively charged acidic mucopolysaccharides in the presence of magnesium ions (Scott et al. 1996).

GUIDELINES

This protocol is optimized for 5 days post-fertilization (dpf) of zebrafish larvae. It can be adapted for higher dpf.

MATERIALS

A	B	C
Abbreviation	Chemical Name	Molecular Weight (g/mol)
Na ₂ HPO ₄	Disodium hydrogen phosphate	141.96
MgCl ₂	Magnesium chloride	95.21
KCl	Potassium chloride	74.55
KH ₂ PO ₄	Potassium dihydrogen phosphate	136.09
KOH	Potassium hydroxide	56.11
NaCl	Sodium chloride	58.44
dH ₂ O	Distilled water	-
EtOH	Ethanol	46.07
PFA	Paraformaldehyde	30.03
PBS	Phosphate-buffered saline	-

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Nov 25, 2023

Last Modified: Dec 11, 2023

PROTOCOL integer ID:
91415

Keywords: Cartilage, zebrafish, alcian blue

Funders

Acknowledgement:

ICMR-SRF

Grant ID: No.45/51/2020-Nan/BMS

A	B	C
PBST	PBS plus 0.1% Tween 20	-
MeOH	Methanol	32.04

1. 10X PBS stock; pH 7.4: To prepare 1 L stock, add 1.37 M (80 g) NaCl, 27 mM (0.2 g) KCl, 100 mM (14.4 g) Na₂HPO₄, and 18 mM (0.24 g) KH₂PO₄ in 800 ml of dH₂O. Make up the volume to 1000 ml and autoclave. Store at room temperature. To prepare 1X PBS, add 100 ml of 10X PBS to 900 ml dH₂O.
2. 1X PBST: Add 0.1 ml of Tween 20 to 100 ml of 1X PBS and mix.
3. 4 % PFA; pH 7.4: Dissolve 4 g of PFA in 80 ml of PBS at 60 °C, until the solution clears. Cool the solution and make it to 100 ml. To avoid repeated thawing and freezing, PFA can be aliquoted in 5 ml tubes and stored at -20 °C.
4. Bleach solution: For 1 ml, add 100 µl of 30% solution of H₂O₂, 250 µl of 2% KOH and 650 µl of dH₂O.
5. Staining solution: For 1 ml, add 100 µl of 10X alcian blue (0.01%), 200 µl of 1 M MgCl₂ (60 mM), and 700 µl of 100% EtOH (70%). a) 10X of alcian blue solution (8 GX; 05500, Sigma Aldrich): For 10 ml solution, add 10 mg alcian blue and 10 ml of methanol. b) 300 mM of MgCl₂: Add 0.29 g of MgCl₂ to 10 ml of dH₂O. c) 100% EtOH.
6. Storage solution: Keep 100% glycerol and 1% KOH (100 mg in 10 ml) separately.
 - a. 25% solution: For 1 ml, mix 250 µL glycerol, 250 µl of KOH, and 500 µl of dH₂O.
 - b. 50% solution: For 1 ml, mix 500 µL glycerol, 250 µl of KOH, and 250 µl of dH₂O.

SAFETY WARNINGS



Make sure to read all Safety Data Sheets for the reagents. Hydrogen peroxide and paraformaldehyde solution must be prepared under a fume hood at all times.

BEFORE START INSTRUCTIONS

1. Every protocol needs to be validated in the laboratory when first introduced. The present protocol describes validation steps that were taken in the Molecular Toxicology lab, at Bharathiar University.
2. Always wear personal protective equipment.
3. All the steps are optimized in a 24-well plate for multiple groups such as control and treated.

Larval fixation

3h

- 1 At the desired stage, take the zebrafish larvae and wash in PBS for 5 minutes.

5m

- 2 Euthanize the larvae using the cold shock method (Keep at 4°C for 5-10 minutes). 10m
- 3 Transfer the euthanized larvae into 4% PFA and keep it in the rocker for 2 hours at room temperature (2h required, it can be kept overnight at 4°C).
- 4 Wash the fixed larvae with PBST for 5 minutes, two times. 10m
- 5 Dehydrate the larvae using 50% ethanol and keep it in the rocker for 10 minutes at room temperature 10m

Staining

8h

- 6 Add 1 ml of alcian blue staining solution to the dehydrated larvae and keep it in a rocker overnight at room temperature. 8h

Bleaching

40m

- 7 Rinse the stained larvae once using PBST. 1m
- 8 Transfer the larvae to a bleach solution and incubate in dark at room temperature. Assess the bleach process (using a microscope) and stop the reaction when the larvae become transparent (approximately 20-30 min for 4 dpf and above). 30m

Clearing

9h

9 Remove the bleach solution and rinse once using PBST.

1m

10 Clear the beached larvae with 25% glycerol and 0.25% KOH. Keep them in a rocker at room temperature for 30 minutes to overnight.

30m

11 Replace the solution with 50% glycerol and 0.25% KOH. Keep them in a rocker at room temperature for hours to overnight.

8h

Storage

12 Store larvae in a solution of 50% glycerol and 0.1% KOH at 4°C. Avoid long-term storage as it will diminish the stain. Capture images of the cartilage using any bright field microscope (preferably a stereo microscope).

Results

13

Expected result




Figure 1. Representative staining images of 5 dpf larvae. a) Whole larvae; b) head region (taken in molecular toxicology lab, Bharathiar University); c) Schematic representation of cartilaginous structures in the head region of zebrafish (adapted from Raterman et al. 2020).