



VERSION 2
AUG 09, 2023

ITS1 Amplicon Prep V.2

aglazer¹

¹TGen North



aglazer

ABSTRACT

This protocol was provided by IDT in the the xGen ITS1 Amplicon Panel manual.

ATTACHMENTS

[xgen-16s-amplicon-panel-v2-and-xgen-its1-amplicon-panel-protocol.pdf](#)

OPEN  ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.n92ldmzd8l5b/v2

Protocol Citation: aglazer
2023. ITS1 Amplicon Prep.
protocols.io
<https://dx.doi.org/10.17504/protocols.io.n92ldmzd8l5b/v2>
Version created by [aglazer](#)

MANUSCRIPT CITATION:

xGen 16S Amplicon Panel v2
xGen ITS1 Amplicon Panel.
(n.d.). Retrieved July 30,
2023, from
https://sfvideo.blob.core.windows.net/sitefinity/docs/default-source/protocol/xgen-16s-amplicon-panel-v2-and-xgen-its1-amplicon-panel-protocol.pdf?sfvrsn=acace007_16

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development
We are still developing and optimizing this protocol

Created: Aug 08, 2023


Last Modified: Aug 09, 2023

PROTOCOL integer ID:
86191

Prepare panel-specific Multiplex PCR Reaction Mix

- 1 Before mixing reagents, calculate the total volume of the Master Mix based on the number of reactions required, with 15% overage for pipetting.
- 2 Make the Multiplex PCR Reaction Mix according to the table below and **pipette-mix** the Master Mix thoroughly.

Note

Vortex components G1 (ITS1 Amplicon Panel Primer) and G2 and pulse-spin tubes to collect contents before use. Gently shake Enzyme G3 at room temperature for  00:05:00 minutes or until all solutes appear to be in solution. Place back on ice for the remainder of the protocol.

Note

Keep prepared Master Mix and all components **ON ICE** until ready to use.

Reagents	Volume per sample (uL)
Reagent G1 - ITS1 Amplicon Panel Primer	2
Reagent G2	3
Enzyme G3	15
Total Volume	20

Multiplex PCR Master Mix Recipe

- 3 Use a multichannel to pipette out  20 µL the PCR Mastermix into each well of a mid-skirt plate as needed and keep ON ICE until template addition.

Perform Multiplex PCR

- 4 Pre-program the thermal cycler for multiplex PCR according to the table below and **allow the block to reach 98°C before loading samples**



Note

Confirm that lid heating is turned **ON** and is set to 105°C

A	B	C
Preheating Step (SKIP TO RUN)	98C	Infinity
	98C	30 sec
4 cycles	98C	10 sec
	63C	5 min

A	B	C
	65C	1 min
14 cycles	98C	10 sec
	64C	1 min
	65C	1 min
	4C	Hold

Multiplex PCR Thermal Cycler Program

- 5 Use a multichannel to add  10 µL of template to the  20 µL Multiplex PCR Reaction Mix and pipette-mix thoroughly.
- 6 Place the mid-skirt plate containing template and Master Mix in the thermal cycler and run the Multiplex PCR Program.

Note

Near the completion of the thermal cycler run, prepare the Indexing Reaction Mix in the post-PCR area. **Assemble this reaction mix on ice and keep cold until adding it to samples in the indexing step.** All components except indexes may be added to the Master Mix when running multiple samples in parallel.

Prepare Indexing PCR Master Mix (while the samples are on..

- 7 Before mixing reagents, calculate the total volume of master mix based on the number of reactions required, with 15% overage for pipetting. It is recommended to prepare at least 10 reactions at any one time to maintain a volume of Enzyme I3 that can be accurately pipetted.
- 8 Prepare the Indexing Master Mix **ON ICE** according to the following table

Note

All components **except indexes** may be added to the Master Mix when running multiple samples in parallel.


Component	Volume per Sample (uL)
Reagent I1	3.3
Enzyme I2	0.5
Enzyme I3	0.1
Enzyme I4	25
Total Volume	28.9

Indexing Master Mix Recipe

- 9 Set aside the Indexing MMX **on ice** until template addition.



Perform Post-Multiplex PCR Cleanup

10m 30s


- 10 Make sure that beads and samples are at room temperature. Spin down the sample plate and briefly vortex beads to homogenize before use.
- 11 Add  30 μ L of magnetic beads (ratio: 1.0x) to each 30 uL sample after thermal cycling and mix by **vortexing**. Pulse-spin the samples in a microcentrifuge tube to collect contents.



Note

Make sure no bead-sample suspension droplets are left on the sides of the tube. If droplets are visible, pulse-spin the samples in a microcentrifuge to collect contents.

12 Incubate the samples for  00:05:00 minutes at room temperature, off the magnet. After incubation, place the sample tubes on a magnetic rack until the solution clears and a pellet is formed (about  00:05:00 minutes).

10m

13 While the sample is on the magnet, remove and discard the supernatant without disturbing the pellet (about  5 μL may be left behind) using a clean pipette tip. Leave the plate on the magnet.



14 Add  180 μL of **freshly prepared** 80% ethanol solution to the pellet while it is still on the magnet. Use care not to disturb the pellet. Incubate for  00:00:30 seconds, then carefully remove the ethanol solution.

30s


15 Repeat the previous step, for a second wash with the ethanol solution.

16 Pulse-spin the samples in a microcentrifuge, place back onto the magnet, and remove any residual ethanol solution from the bottom of the tube with a small volume tip.

17 Resuspend each bead pellet in  19.4 μL post-PCR TE Buffer off-magnet and pipette mix.

18 Incubate the samples for  00:05:00 minutes at room temperature, off the magnet. After incubation, place the sample tubes on a magnetic rack until the solution clears and a pellet is formed (about  00:05:00 minutes).

10m

19 While the sample is on the magnet, remove  17.4 μL of the supernatant without disturbing the pellet using clean pipette tips and transfer to a new plate.

Perform Indexing PCR

- 20 Program the thermal cycler with the Indexing PCR program described in the table below and **allow the block to reach 37°C before loading samples.**

Note





Confirm lid heating is turned **ON** and set to reach 105°C.

Note

The PCR cycle number can be increased for samples that may give low yields.

A	B	C
Preheating Step (SKIP TO RUN)	37C	Infinity
	37C	20 min
	98C	30 sec
7 cycles	98C	10 sec
	60C	30 sec
	66C	1 min
	4C	Hold

Indexing PCR Program


- 21 Add  3.7 µL of the pre-mixed xGen Amplicon UDI primer pair to each sample ( 17.4 µL).
- 22 Add  28.9 µL of the cold indexing PCR reaction mix, stored on ice, to each sample. Mix thoroughly by **pipetting** (total volume will be  50 µL).

- 23 Spin down the sample plate, place in the thermal cycler, and run the Indexing PCR program.

Perform Post-Indexing PCR Cleanup



12m 30s

- 24 Make sure the beads are at room temperature. Spin down the sample plate and briefly vortex the beads to homogenize before use.


- 25 Add  42.5 µL of fresh magnetic beads (ratio: 0.85x) to each 50 µL sample after thermal cycling and mix by **vortexing**.



Note

Make sure no bead-sample suspension droplets are left on the sides of the tube. If droplets are visible, pulse-spin the samples in a microcentrifuge to collect contents.

- 26 Incubate the samples for  00:05:00 minutes at room temperature, off the magnet. Pulse-spin the samples in a microcentrifuge. Afterwards, place the sample tubes on a magnetic rack until the solution clears and a pellet is formed (about  00:05:00 minutes).




10m

- 27 While the sample is on the magnet, remove and discard the supernatant without disturbing the pellet (about  5 µL may be left behind) using a clean pipette tip. Leave the plate on the magnet.

- 28 Add  180 µL of **freshly prepared** 80% ethanol solution to the pellet while it is still on the magnet. Use care not to disturb the pellet. Incubate for  00:00:30 seconds, then carefully remove the ethanol solution.

30s

- 29 Repeat the previous step, for a second wash with the ethanol solution.

- 30 Pulse-spin the samples in a microcentrifuge, place back onto the magnet, and remove any residual ethanol solution from the bottom of the tube with a small volume tip.
- 31 Immediately add  22 µL of post-PCR TE buffer and resuspend the pellet by pipetting up and down until homogenous. Incubate for  00:02:00 minutes off the magnet. 2m
- 32 Place the samples back on the magnet and transfer the clean  20 µL library eluates to a fresh plate. Make sure that the eluate does not contain magnetic beads (indicated by brown color in eluate). If magnetic beads are present, place the plate on the magnet, wait for pellets to re-form, then transfer eluate again.
- 33 **SAFE STOP:** Store freshly prepared libraries at -20°C.