



FEB 16, 2024

OPEN ACCESS



DOI:

[dx.doi.org/10.17504/protocols.io.j8nlkoqdxv5r/v1](https://dx.doi.org/10.17504/protocols.io.j8nlkoqdxv5r/v1)

**Protocol Citation:** Humaira Noor, Richard Wade-Martins 2024. Quantification of glutamate released from human induced pluripotent stem cells (iPSC) derived cortical neurons (CNs).

protocols.io

<https://dx.doi.org/10.17504/protocols.io.j8nlkoqdxv5r/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working

# Quantification of glutamate released from human induced pluripotent stem cells (iPSC) derived cortical neurons (CNs)

In 1 collection

Humaira Noor<sup>1,2,3,4</sup>, Richard Wade-Martins<sup>1,2,4</sup>

<sup>1</sup>Oxford Parkinson's Disease Centre and Department of Physiology, Anatomy and Genetics, University of Oxford, South Park Road, Oxford OX1 3QU, United Kingdom;

<sup>2</sup>Kavli Institute for Neuroscience Discovery, University of Oxford, Dorothy Crowfoot Hodgkin Building, South Park Road, Oxford OX1 3QU, United Kingdom;

<sup>3</sup>Nuffield Department of Medicine (NDM), University of Oxford, Henry Wellcome Building for Molecular Physiology, Old Road, Oxford OX3 7BN, United Kingdom;

<sup>4</sup>Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815, USA

Team Cragg



Cláudia C. Mendes

## ABSTRACT

This protocol describes the procedure to quantify glutamate released from induced pluripotent stem cells derived cortical neurons (iPSC-CN) cultured in adherent monolayer.

## MATERIALS

### Reagents:

- [Glutamate Assay Kit](#) (Abcam, CAT# ab83389) - [Instructions for Use](#)
- [Hanks' Balanced Salt Solution \(HBSS++\)](#) (ThermoFisher, CAT# 14025092)
- [Potassium chloride \(KCl\)](#) (Sigma Aldrich, CAS# 7447-40-7)
- [PierceTM BCA Protein Assay Kit](#) (ThermoFisher, CAT# 23225)

### Equipment:

- 96-well microreader plate
- [PHERAstar® FSX microplate reader](#) (BMG LabTech)

**Created:** Feb 05, 2024

## BEFORE START INSTRUCTIONS

**Last Modified:** Feb 16, 2024

**PROTOCOL integer ID:** 94701

**Keywords:** Glutamate release, human neurons, tonic and evoked release

**Funders Acknowledgement:**

Aligning Science Across  
Parkinson's (ASAP)  
Grant ID: ASAP-020370

Conditioned media used in this assay can come from any human cell type capable of producing and releasing glutamate. However, the protocol has been optimised for iPSC-derived cortical neurons differentiated according to the following **Protocol: Differentiation of human cortical neurons (CNs) from induced pluripotent stem cells (iPSCs)**, and plated at 10,000 cells per well in half area 96-well plates and cultured in 100 µL of media.

Prior to use, media must be warmed preferentially to 37°C as the cells are very temperature sensitive.

Plate cells in a layout so that you have at least triplicate readouts for each cell line and/or condition.

## Collection of Condition Media

### 1 Collection of conditioned media for tonic release:

- 1.1 Remove the cell culture media and add 60 µl of pre-warmed HBSS++ (2.4 mM KCl) media to the cells.
- 1.2 Incubate the cells for 10 min before collecting the media.
- 1.3 Collect the media carefully and slowly making sure not to disturb the cells while adding and removing media.

### 2 Collection of conditioned media for evoked release:

- 2.1 Prepare 40 mM KCl HBSS++ solution by adding 29.82 mg KCl to 10ml of HBSS++. Filter sterilise prior to use.
- 2.2 Add 60  $\mu$ L of pre-warmed HBSS++ (40 mM KCl) media to the cells and incubate for 5 min. Collect the conditioned media.
- 2.3 Supernatants collected in both instances can be stored in -20°C following snap freezing until further use.

## Measuring Glutamate in Conditioned Media

### 3 Preparation of standards:

- 3.1 Always prepare fresh set of standards on the day of use. The standards are unstable and needs to be used within 4 hours following preparation.
- 3.2 Prepare 1 mM of glutamate standard by diluting 5  $\mu$ L of 0.1 M glutamate with 495  $\mu$ L of the supplied Assay Buffer. Prepare a range of standards starting from 0 nmol to 10 nmol following the instruction booklet.
- 3.3 Prepare a range of lower concentrations of the standards in addition the concentrations mentioned in the manufacturer's instructions (**Step 3.2**) as often the glutamate released by the human neurons are in the lower concentration range (depending on the cell type and the seeded density).

## 4 Assay Procedure:

- 4.1 Bring all the materials and prepared reagents to room temperature before use.
- 4.2 Set up the reaction by adding 50  $\mu$ L of the standards and samples to wells of a 96-well microreader plate according to your plate map.
- 4.3 For calculating background for the samples, use 50  $\mu$ L of fresh HBSS++.
- 4.4 Add 100  $\mu$ L of reaction mix to each well and mix by pipetting up and down.
- 4.5 Incubate the plate at 37°C for 30 minutes protected from the light.
- 4.6 Read OD at 450 nm using a microplate reader for e.g. PHERAstar® microplate reader (BMG Labtech).

## Analysis

- 5 Subtract blank values and average the data obtained in triplicates.

- 6** Plot a standard curve using the data obtained for the standards.
- 7** Calculate glutamate levels in nmol for each sample from the standard curve.
- 8** Normalise the glutamate release data for difference in plating density/cell number by dividing by the total protein values obtained by BCA.