



MAR 15, 2024

## Isolation of NeuN+ cells from brain tissue (for CUT and RUN)

anita.adami<sup>1</sup>

<sup>1</sup>Laboratory of Molecular Neurogenetics, Department of Experimental Medical Science, Wallenberg Neuroscience Center and Lund Stem Cell Center, BMC A11, Lund University, 221 84 Lund, Sweden.

ASAP Collaborative Research Network

Jakobsson



Raquel Garza  
Lund University

### ABSTRACT

This protocol describes the steps to isolate NeuN+ cells from brain tissue in preparation for CUT and RUN

OPEN ACCESS



#### DOI:

[dx.doi.org/10.17504/protocols.io.4r3l27pejg1y/v1](https://dx.doi.org/10.17504/protocols.io.4r3l27pejg1y/v1)

**Protocol Citation:** anita.adami 2024. Isolation of NeuN+ cells from brain tissue (for CUT and RUN). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.4r3l27pejg1y/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
We use this protocol and it's working

**Created:** Mar 03, 2023

Last Modified: Mar 15, 2024

PROTOCOL integer ID: 78067

Keywords: ASAPCRN

Funders Acknowledgement:

Aligning Science Across  
Parkinson's through the Michael  
J. Fox Foundation for Parkinson's  
Research  
Grant ID: ASAP-000520  
Swedish Research Council  
Grant ID: 2018-02694  
Swedish Brain Foundation  
Grant ID: FO2019-0098  
Cancerfonden  
Grant ID: 190326  
Barncancerfonden  
Grant ID: PR2017-0053  
NIHR Cambridge Biomedical  
Research Centre  
Grant ID: NIHR203312  
Swedish Society for Medical  
Research  
Grant ID: S19-0100  
National Institutes of Health  
Grant ID: HG002385  
Swedish Research Council  
Grant ID: 2021-03494  
Swedish Research Council  
Grant ID: 2020-01660

Isolation of NeuN+ cells

45m

- 1 Nuclei were isolated from frozen tissue as described [here](#).
- 2 Before FACSing, nuclei were incubated with Recombinant Alexa Fluor® 488 Anti-NeuN antibody [EPR1276 - Neuronal Marker (Abcam, Cat# ab190195, RRID:AB\_2716282) at a concentration of 1:500 for 

00:30:00 On ice

 as according to [Spalding, K.L. et al. \(2013\)](#). 

30m

3 The nuclei were run through the FACS at 4 °C with a low flow rate using a 100 mm nozzle and 300.000 nuclei Alexa Fluor – 488 positive nuclei were sorted.

4 The sorted nuclei were pelleted at 1,300 x g for 00:15:00 and resuspended in 1 mL of ice-cold nuclei wash buffer (20 mM HEPES, 150 mM NaCl, 0.5 mM spermidine, 1x cOmplete protease inhibitors, 0.1% BSA) and 10 µL per antibody treatment of ConA-coated magnetic beads (Epicyphe<sup>15m</sup>) added with gentle vortexing (Pipette tips for transferring nuclei were pre-coated with 1% BSA), to perform CUT&RUN.