

NOV 05, 2023

## Sleep-Wake Recording in Dan lab

Xiaolin Huang<sup>1</sup>

<sup>1</sup>University of California, Berkeley



Xiaolin Huang

**ABSTRACT** 

Surgery and Setups for Recording 24-Hour Natural Sleep/Wake Cycles in Mice.

# OPEN ACCESS



#### DOI:

dx.doi.org/10.17504/protocol s.io.x54v9pkd4g3e/v1

Protocol Citation: Xiaolin Huang 2023. Sleep-Wake Recording in Dan lab. protocols.io https://dx.doi.org/10.17504/p rotocols.io.x54v9pkd4g3e/v1

License: This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working We use this protocol and it's working

Created: Oct 27, 2023

Last Modified: Nov 05,

2023

#### **PROTOCOL** integer ID:

89994

#### Funders Acknowledgement:

**ASAP** 

Grant ID: 020551



- 1 Mice were anesthetized with 1.5% isoflurane and placed in a stereotaxic frame.
- 2 Body temperature was maintained using a heating pad.
- 3 The skin was incised to expose the skull after asepsis and connective tissue was removed.
- To implant electroencephalogram (EEG) and electromyogram (EMG) recording electrodes, two stainless steel screws were inserted into the skull 2.5 mm from the midline and 3 mm posterior to the bregma, and two EMG electrodes were also inserted into the neck muscles.
- **5** A reference screw for grounding was placed on top of the left cerebellum.
- 6 Insulated leads from the EEG and EMG electrodes were soldered to a pin header, which was secured to the skull using dental cement.

### Recording

7 After a minimum of 7 days of undisturbed post-surgery recovery, behavioral

experiments were carried out in home cages placed in sound-attenuating boxes.

- **8** EEG and EMG electrodes were connected to flexible recording cables via a mini-connector.
- **9** Recordings started after 1-2 days of habituation and continued for 3-4 days on a 12-h dark/12-h light cycle.
- Mouse behavior was recorded using a video camera at 10 frames per second (fps), along with EEG and EMG recordings.
- 11 The signals were recorded with a TDT RZ10x/PZ5-32 or LR10-SK1 system.