



Jun 06, 2022

EMP 16S rRNA PCR

Allyson Hirsch¹, George Testo¹¹The Pathogen & Microbiome Institute, Northern Arizona University

1



protocol .

**George Testo**
The Pathogen & Microbiome Institute

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

The 16S protocol detailed here is designed to amplify prokaryotes (bacteria and archaea) using paired-end 16S community sequencing on the Illumina platform. Primers 515F–806R target the V4 region of the 16S SSU rRNA.

Allyson Hirsch, George Testo 2022. EMP 16S rRNA PCR. **protocols.io**
<https://protocols.io/view/emp-16s-rna-pcr-cakescte>



_____ protocol ,

Jun 03, 2022

Jun 06, 2022

63846

Reagents

- PCR Water
- PCR Reagents (buffer, Forward Primer, dNTPs, BSA, Taq)

Supplies

- 50mL conical(s)
- Unskirted plate(s)
- Green temporary seal(s)
- PCR seal(s)
- Foil seal(s)
- Sharpie

Equipment

- Plate sealer
- Conical rack
- Plate holder(s)
- Boat & covers
- 200uL pipette & tips
- 1000uL pipette & tips
- 10uL 12-channel pipette & tips
- 200uL 12-channel pipette & tips

:

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

The primer sequences without linker, pad, barcode, or adapter are as follows:

Updated sequences: 515F (Parada)–806R (Aprill), forward-barcoded:
FWD:GTGYCAGCMGCCGCGGTAA; REV:GGACTACNVGGGTWTCTAAT

Original sequences: 515F (Caporaso)–806R (Caporaso), reverse-barcoded:
FWD:GTGCCAGCMGCCGCGGTAA; REV:GGACTACHVGGGTWTCTAAT

Note: Primers used for PCR reaction(s) are determined by project and or lab.

- 1 Prepare PCR cocktail calculations and a plate map for all PCR reactions using barcoded primers.

The following calculations must be made based on the number of samples to be amplified:

1	PCR Cocktail for 16S ribosomal amplification			
2	16S Illumina amplicons			
3	Final Volume	Number of samples (incl controls)		
4	25	300		
5	Component	Conc	Volume/ reaction	[final]
6	PCR buffer	10X	2.5	1X
7	Forward primer	10uM	1	.4uM
8	Reverse primer	10uM	1	.4uM
9	BSA	20mg/ml	1	0.56mg/ul
10	dNTP mix	2.5mM each	2	200 uM each
11	Hot Start ExTaq	5U/ul	0.125	0.625
12	Template	10-100 ng	2	n/a
13	Total Reagents		9.625	
14	H2O		15.375	4843.13
15	Total volume		25	7875
16	Volume per well/tube		23	
17				
18	Date			
19	Name			

Note: The numbers used above are only for demonstration purposes. If calculating the amount of PCR reagents necessary for a set of PCR reactions, calculate using the demonstrated concentrations and volumes per reaction. Typically, every 6 triplicates (or 6 samples + NTCs) requires an input of 25 samples. Therefore, 1 unskirted plate of PCR reactions equates to 100 samples (as the required input to the PCR cocktail calculator as shown above).

Preparing PCR Hood 30m


- 2 Decontaminate a PCR hood with DNA Away and Ethanol.
- 3 Place materials inside of the PCR Hood (refer to materials section).

- 4 Turn on the UV setting for at least 00:30:00 .

30m

Preparing Master Mix 30m

- 5 Label unskirted plate(s) with the following: Project, Extraction Name, Initials, Date, & Plate Number


- 6 Pour PCR water into a 50mL conical and make the master mix (specific to project and or number of samples).
- 7 Using a 200uL 12-channel pipette, pipette  **22 µL** of master mix into each well on each unskirted plate.


Note: Samples will be amplified in triplicate. This means that each sample will be amplified in 3 replicate 25uL PCR reactions. This ensures that enough amplicons are created for subsequent preparations. Triplicate PCR reactions for each sample will be pooled together during EGels.

Adding Reverse Primer 30m

8



Once finished, retrieve EMP  **10 micromolar (µM)** Reverse Primer Plate and DNA plate. Let thaw and spin down.

- 9 Add  **1 µL** of Reverse Primer according to the plate map specific for the project and sample(s).

Adding Sample(s) 30m

- 10 Add  **2 µL** of sample according to the plate map specific for the project and sample(s).

- 11 Seal everything with PCR and or freezer seals and clean the PCR hood.

11.1 Return sample(s), reagents, and EMP Reverse Primer plate to the appropriate freezer.

12



Spin down PCR plate(s).

Running Thermocycler

30m

13

Place plate(s) on the thermocycler and start an amplification protocol.

13.1

Use the following thermocycler conditions as a part of this earth microbiome project (EMP) protocol:

A	B	C
98C	3 mins	1x
98C	45 secs	30x
50C	60 secs	30x
72C	90 secs	30x
72C	10 mins	1x
4C	hold	

14

Once the PCR has finished on the thermocycler, place PCR plate(s) in a **-20 °C** freezer or run a gel.