

May 10, 2024

# Protein Extraction for Amyloid Beta Fractionation

DOI

**dx.doi.org/10.17504/protocols.io.j8nlk8ky1l5r/v1**

Lisa Blackmer-Raynolds<sup>1</sup>, Tim Sampson<sup>1</sup>

<sup>1</sup>Emory University



Lisa Blackmer-Raynolds

Emory University

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.j8nlk8ky1l5r/v1](https://dx.doi.org/10.17504/protocols.io.j8nlk8ky1l5r/v1)

**Protocol Citation:** Lisa Blackmer-Raynolds, Tim Sampson 2024. Protein Extraction for Amyloid Beta Fractionation. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.j8nlk8ky1l5r/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** May 08, 2024

**Last Modified:** May 10, 2024

**Protocol Integer ID:** 99473

**Keywords:** amyloid beta, fractionation, solubility, protein

**Funders Acknowledgement:**

NIH

Grant ID: F31AG076332-02

ASAP

Grant ID: ASAP-020527



## Abstract

This protocol can be used to isolate protein from murine brain tissue (volumes listed are for hippocampal samples) by solubility. It creates three protein fractions: tris soluble, triton soluble, and formic acid soluble. This protocol can be helpful when evaluating amyloid beta pathology as different amyloid beta structures can be found within each fraction. For example, insoluble amyloid beta associated with amyloid plaques will be found within the formic acid soluble fraction, whereas more soluble non-plaque associated amyloid beta will be found within the other two fractions.



## Preparation

- 1 **5x Homogenization Buffer (store 4°C):**  
125 ml of 1.0 M Tris  
30 ml of 0.5 M MgCl<sub>2</sub>  
25 ml of 0.1 M EDTA (pH to 7.2 and store at 4°C for up to 6 months)  
20ml of MiliQ Water
- 2 **1x Homogenization Buffer (dilute day of):**  
Dilute 5x Homogenization buffer to 1x in MiliQ Water and add 1 tablet protease inhibitor (cat#11697498001 Roche) per 10 mL
- 3 **1x Homogenization Buffer with 1% Triton X (dilute day of):**  
Add 1% Triton X-100 to 1x Homogenization buffer made above
- 4 **70% Formic Acid (store RT) \*\*Caution Acidic\*\***  
7mL formic acid  
3mL water
- 5 **1.5x Homogenization Buffer pH 14 (store 4°C) \*\*Caution Corrosive\*\***  
Dilute 5x Homogenization buffer to 1.5x in MiliQ Water and pH to 14

## Tissue Processing

- 6 Weigh frozen tissue sample and record for later normalization
- 7 Homogenize the tissue on ice using quick bursts of a sonicator in 400µl 1x Homogenization buffer
- 8 Centrifuge maximum speed at 4°C for 1hr
- 9 Collect supernatant as tris soluble fraction and store -80°C for later use
- 10 Homogenize the pellet on ice using quick bursts of a sonicator in 200µl 1x Homogenization Buffer with 1% Triton X-100



- 11 Centrifuge maximum speed at 4°C for 1hr
- 12 Collect supernatant as triton soluble fraction and store -80°C for later use
- 13 Homogenize the pellet on ice using quick bursts of a sonicator in 50µl 70% formic acid \*\*Be careful with acid\*\*
- 14 Centrifuge maximum speed at 4°C for 1hr
- 15 Collect supernatant as formic acid soluble fraction and neutralize with 1.5x Homogenization buffer pH14 on ice \*Be careful with acid and base\*\*. Store -80°C for later use