



# Extraction of Oxycodone from Plasma for Mass Spec Analysis

Sierra Simpson<sup>1</sup>, Olivier George<sup>1</sup>

<sup>1</sup>University of California, San Diego

1 Works for me This protocol is published without a DOI.

George Lab  
Tech. support email: [olgeorge@ucsd.edu](mailto:olgeorge@ucsd.edu)



Sierra Simpson  
University of California, San Diego

## ABSTRACT

From the Momper Lab.

## PROTOCOL CITATION

Sierra Simpson, Olivier George 2020. Extraction of Oxycodone from Plasma for Mass Spec Analysis.  
**protocols.io**  
<https://protocols.io/view/extraction-of-oxycodone-from-plasma-for-mass-spec-bamhic36>

## LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Dec 17, 2019

## LAST MODIFIED

Jul 23, 2020

## PROTOCOL INTEGER ID

31113

## MATERIALS TEXT

Quantitative determination of oxycodone in plasma and brain was accomplished by the use of high-performance liquid chromatography with tandem mass spectrometry detection (LC-MS/MS). Oxycodone was precipitated from 20 µL of plasma with 60 µL of 500ng/mL ISTD ([Oxycodone-D6](#)) in acetonitrile. Brain tissues were homogenized and filtered. 20µl supernatant was injected directly onto a C-18 reversed phase HPLC column (MacMod Ace-5, 2.1 x 150 mm). The LC mobile phase consisted of HPLC grade water with [0.1% formic acid](#) (elute A) and ACN with 0.1% formic acid (elute B), and eluted with a gradient program of 0-25 min / 95%B, flow rate of 0.4 mL/min; 1.00-3.00 min / 95 % B, flow rate of 0.6 mL/min; 3.10-6.00 min / 5 % B, flow rate of 0.3 mL/min. MS/MS detection is made in positive electrospray ionization mode, at mass transitions of 316→241 m/z (Oxycodone), 322→247 m/z (ISTD). The method has a dynamic range of 9.8–313 ng/mL. For analytes in plasma and brain, calibration standards are used to generate a curve using a linear regression algorithm to plot the peak area ratio versus concentration with 1/x weighting, over the full dynamic range of analyte concentrations.

### Harvesting Blood

- 1 Collect blood in a tube coated with EDTA. Heparin leads to hemolysis which can effect downstream processing.
- 2 Spin blood at 3000rpm for 10 minutes to separate erythrocytes from plasma

- 3 Aliquot plasma into 50-100ul aliquits to avoid freeze thaw. 20ul is processed per sample.

#### Plasma Preparation

- 4 Add 20ul plasma to 60ul Acetonitrile spiked with 500ng/ml Oxycodone-d-6 ( Heavy oxy for internal standard)

This will crash out all the extra protein within the sample.

- 5 Spin the sample at 13K for 10 minutes, then extract supernatant and put it into a mass spec vial for sampling.

#### Running the samples

- 6 Load the samples into the autosampler

20µl supernatant is injected directly onto a C-18 reversed phase HPLC column (MacMod Ace-5, 2.1 x 150 mm).

- 7 The LC mobile phase consisted of HPLC grade water with [0.1% formic acid](#) (elute A) and ACN with 0.1% formic acid (elute B), and eluted with a gradient program of 0-.25 min / 95%B, flow rate of 0.4 mL/min; 1.00-3.00 min / 95 % B, flow rate of 0.6 mL/min; 3.10-6.00 min / 5 % B, flow rate of 0.3 mL/min.

MS/MS detection is made in positive electrospray ionization mode, at mass transitions of 316→241 m/z (Oxycodone), 322→247 m/z (ISTD).