

Sep 16, 2020

Quantitative Real-Time Polymerase Chain Reaction (qRT-PCR)

In 1 collection

Yingchao Xue^{1,2}, Xiping Zhan³, Shisheng Sun⁴, Senthilkumar S. Karuppagounder^{5,6,7}, Shuli Xia^{2,5}, Valina L Dawson^{5,6,7,8,9}, Ted M Dawson^{5,6,7,8,10}, John Laterra^{2,5,8,11}, Jianmin Zhang¹, Mingyao Ying^{2,5}

¹Department of Immunology, Research Center on Pediatric Development and Diseases, Institute of Basic Medical Sciences, Chinese Academy of Medical Sciences and School of Basic Medicine, Peking Union Medical College, State Key Laboratory of Medical Molecular Biology;

²Hugo W. Moser Research Institute at Kennedy Krieger; ³Department of Physiology and Biophysics, Howard University;

⁴College of Life Sciences, Northwest University; ⁵Department of Neurology, Johns Hopkins University School of Medicine;

⁶Neuroregeneration and Stem Cell Programs, Institute for Cell Engineering, Johns Hopkins University School of Medicine;

⁷Adrienne Helis Malvin Medical Research Foundation; ⁸Department of Neuroscience, Johns Hopkins University School of Medicine;

⁹Department of Physiology, Johns Hopkins University School of Medicine;

¹⁰Department of Pharmacology and Molecular Sciences, Johns Hopkins University School of Medicine;

¹¹Department of Oncology, Johns Hopkins University School of Medicine

1 Works for me

This protocol is published without a DOI.

Neurodegeneration Method Development Community

Tech. support email: ndcn-help@chanzuckerberg.com



Anita Broellochs
protocols.io

ABSTRACT

This protocol explains Quantitative Real-Time Polymerase Chain Reaction (qRT-PCR) of lines ND1014, N1, and ND27760 from *Synthetic mRNAs Drive Highly Efficient iPS Cell Differentiation to Dopaminergic Neurons*.

EXTERNAL LINK

<https://www.ncbi.nlm.nih.gov/pmc/articles/PMC6344911/>

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Synthetic mRNAs Drive Highly Efficient iPS Cell Differentiation to Dopaminergic Neurons. Xue Y, Zhan X, Sun S, Karuppagounder SS, Xia S, Dawson VL, Dawson TM, Laterra J, Zhang J, Ying M. *Stem Cells Transl Med*. 2019 Feb;8(2):112-123. doi: 10.1002/sctm.18-0036. Epub 2018 Nov 1. PMID: 30387318

PROTOCOL CITATION

Yingchao Xue, Xiping Zhan, Shisheng Sun, Senthilkumar S. Karuppagounder, Shuli Xia, Valina L Dawson, Ted M Dawson, John Laterra, Jianmin Zhang, Mingyao Ying 2020. Quantitative Real-Time Polymerase Chain Reaction (qRT-PCR). **protocols.io**
<https://protocols.io/view/quantitative-real-time-polymerase-chain-reaction-q-9u8h6zw>

MANUSCRIPT CITATION please remember to cite the following publication along with this protocol

Synthetic mRNAs Drive Highly Efficient iPS Cell Differentiation to Dopaminergic Neurons. Xue Y, Zhan X, Sun S, Karuppagounder SS, Xia S, Dawson VL, Dawson TM, Laterra J, Zhang J, Ying M. *Stem Cells Transl Med*. 2019 Feb;8(2):112-123. doi: 10.1002/sctm.18-0036. Epub 2018 Nov 1. PMID: 30387318

EXTERNAL LINK

<https://www.ncbi.nlm.nih.gov/pmc/articles/PMC6344911/>

COLLECTIONS ⓘ




Protocols for Synthetic mRNAs Drive Highly Efficient iPS Cell Differentiation to Dopaminergic Neurons

KEYWORDS

ND1014, N1, ND27760, ipsc, SNCA, Atoh2, Ngn2, qRT-PCR

LICENSE

 This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Nov 27, 2019

LAST MODIFIED

Sep 16, 2020

OWNERSHIP HISTORY

Nov 27, 2019  Liz Brydon [Protocols.io](#)

Sep 16, 2020  Anita Broellochs [protocols.io](#)

PROTOCOL INTEGER ID

30336

PARENT PROTOCOLS

Part of collection

[Protocols for Synthetic mRNAs Drive Highly Efficient iPS Cell Differentiation to Dopaminergic Neurons](#)

GUIDELINES

Primers for qRT-PCR

Primer	Sequence (5'-3')
NEUROD1	AGACACTCGTCTGTCCAGCTT; GCCCCAGGGTTATGAGACTA
FOXA2	GGAACACCACTACGCCTTCAAC; AGTGCATCACCTGTTCTAGGC
NURR1	AAACTGCCCAGTGGACAAGCGT; GCTCTTCGGTTTCGAGGGCAAA
LMX1A	CATCGAGCAGAGTGTCTACAGC; TGTCGTCGCTATCCAGGTCATG
TH	GCTGGACAAGTGTCTACCTG; CCTGTAAGGGAAGGCGATCTCA
IFNA	ACCCACAGCCTGGATAACAG; ACTGGTTGCCATCAAATCC
IFNB	CATTACCTGAAGGCCAAGGA; CAGCATCTGCTGGTTGAAGA
IFIT1	AAAAGCCACATTTGAGGTG; GAAATTCCTGAAACCGACCA
OAS1	CGATCCCAGGAGGTATCAGA; TCCAGTCTCTTCTGCCTGT
PKR	TCGCTGGTATCACTCGTCTG; GATTCTGAAGACCGCCAGAG
RIGI	GTTGTCCCCTGCTGTTCTT; GCAAGTCTTACATGGCAGCA
18s	ACCCGTTGAACCCATTCTGTA; GCCTCACTAAACCATCCAATCGG

MATERIALS

NAME	CATALOG #	VENDOR
SYBR Green		Life Technologies
RNeasy Mini Kit	74104	Qiagen

SAFETY WARNINGS

Please refer to the Safety Data Sheets (SDS) for safety and environmental hazards.

BEFORE STARTING

Obtain approval to work with human stem cells from an appropriate Institutional Review Board.

RNA Extraction

- 1 Extract total RNA using a RNeasy Mini Kit.

qRT-PCR

- 2 Reverse transcribe using murine leukemia virus reverse transcriptase and oligo(dT) primers.
- 3 Set up qRT-PCR using SYBR green PCR Master Mix in the IQ5 RT-PCR detection system. Primer sequences are listed in "Guidelines".
- 4 Normalize relative expression to 18S rRNA.