

Jun 18, 2024

🌐 Post GEM–RT Cleanup and cDNA Amplification

📁 In 1 collection

DOI

dx.doi.org/10.17504/protocols.io.bp2l62d3dgqe/v1



Heidi Monroe¹, Nayra Cardenes², Melanie Königshoff², koenigshoffm², Robert Lafyatis¹

¹University of Pittsburgh;

²Division of Pulmonary, Allergy, Critical Care and Sleep Medicine, School of Medicine, University of Pittsburgh

TriState SenNet

Cellular Senescence Net...

1 more workspace



Nayra Cardenes

Division of Pulmonary, Allergy, Critical Care and Sleep Med...

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.bp2l62d3dgqe/v1

Protocol Citation: Heidi Monroe, Nayra Cardenes, Melanie Königshoff, koenigshoffm, Robert Lafyatis 2024. Post GEM–RT Cleanup and cDNA Amplification. [protocols.io https://dx.doi.org/10.17504/protocols.io.bp2l62d3dgqe/v1](https://dx.doi.org/10.17504/protocols.io.bp2l62d3dgqe/v1)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: June 12, 2024

Last Modified: June 18, 2024

Protocol Integer ID: 101648

Keywords: Post GEM–RT Cleanup, cDNA amplification, Dynabeads, cDNA quantification, SenNet, TriState, snRNAseq, Lung, PCLS, Frozen tissue

Funders Acknowledgement:
TriState SenNET (Lung and Heart) Tissue Map and Atlas consortium - NIA
Grant ID: U54AG075931

Abstract

The Chromium Single Cell Gene Expression Solution upgrades short read sequencers to deliver a scalable microfluidic platform for 3' digital gene expression by profiling 500-10,000 individual cells per sample.

A pool of ~3,500,000 10x Barcodes are sampled separately to index each cell's transcriptome. It is done by partitioning thousands of cells (or nuclei) into nanoliter-scale Gel Beads-in-emulsion (GEMs), where all generated cDNA share a common 10x Barcode. Dual Indexed libraries are generated and sequenced from the cDNA and 10x Barcodes are used to associate individual reads back to the individual partitions.

After the GEMs Generation and Barcoding, GEMs are broken and pooled fractions are recovered. Silane magnetic beads are used to purify the first-strand cDNA from the post GEM-RT reaction mixture, which includes leftover biochemical reagents and primers. Barcoded, full-length cDNA is amplified via PCR to generate sufficient mass for library construction.

This protocol details the post GEM-RT cleanup and cDNA amplification, cleanup - SPRIselect and quantification.

Attachments






[snRNAseq_ProtocolsIO..](#)

⋮
66KB

Image Attribution

Nayra Cardenes, PhD

Materials

-  Dynabeads MyOne Silane **10x Genomics Catalog #2000048**
-  Reducing Agent B **10x Genomics Catalog #2000087**
- cDNA Primers (2000089)
-  Dynabeads MyOne Silane **10x Genomics Catalog #2000048**
- Amp Mix (2000047/2000103)
- Cleanup Buffer (2000088)

Dynabeads Cleanup Mix:

A	B	C	D
Reagents	1X (μl)	4X+10% (μl)	8X+10% (μl)
Cleanup Buffer	182	801	1602
Dynabeads MyOne SILANE	8	35	70
Reducing Agent B	5	22	44
Nuclease-free Water	5	22	44
Total	200	880	1760

Elution Solution I:

A	B	C
Reagents	1X (μl)	10X (μl)
Buffer EB	98	980
10% Tween 20	1	10
Reducing Agent B	1	10
Total	100	1000

cDNA Amplification Reaction Mix:

A	B	C	D
Reagents	1X (μl)	4X+10% (μl)	8X+10% (μl)
Amp Mix	50	220	440
cDNA Primers	15	66	132
Total	65	286	572



Safety warnings





The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to **protocols.io** is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with **protocols.io**, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Before start

- **Equilibrate to Room temperature (RT)** – Reducing Agent B (2000087), cDNA Primers (2000089) and Dynabeads MyOne SILANE (2000048)
- **Place on ice** – Amp Mix (2000047/2000103)
- **Thaw at 65°C**- Cleanup Buffer (2000088)


Dynabeads

29m

1 Add  125 µL Recovery Agent to each sample at  Room temperature .




Note

DO NOT pipette mix or vortex the biphasic mixture and wait  00:02:00 .

The resulting biphasic mixture contains Recovery Agent/Partitioning Oil (pink) and aqueous phase (clear), with no persisting emulsion (opaque).

Note

If biphasic separation is incomplete: Firmly secure the cap on the tube strip, ensuring that no liquid is trapped between the cap and the tube rim. Mix by inverting the capped tube strip 5x, centrifuge briefly, and proceed to step 2.

2 Slowly remove and discard  125 µL Recovery Agent/Partitioning Oil (pink) from the bottom of the tube.


Note

DO NOT aspirate any aqueous sample.

3 Prepare Dynabeads Cleanup Mix:

A	B	C	D
Reagents	1X (µl)	4X+10% (µl)	8X+10% (µl)
Cleanup Buffer	182	801	1602
Dynabeads MyOne SILANE	8	35	70
Reducing Agent B	5	22	44
Nuclease-free Water	5	22	44
Total	200	880	1760

Calculations for Dynabeads Cleanup Mix preparation.

4 Vortex and add  200 μ L to each sample. Pipette mix 10x (pipette set to 200 μ L).




5 Incubate  00:10:00 at  Room temperature .

15m

Note



Keep caps open.

Pipette mix again at ~  00:05:00 after start of incubation to resuspend settled beads.


6 Prepare Elution Solution I. Vortex and centrifuge briefly.



Elution Solution I

A	B	C
Reagents	1X (μ l)	10X (μ l)
Buffer EB	98	980
10% Tween 20	1	10
Reducing Agent B	1	10
Total	100	1000

Calculations for Elution Solution I preparation.

7 At the end of  00:10:00 incubation, place on a 10x Magnetic Separator. High position (magnet.High) until the solution clears.

10m



8 Remove the supernatant (aqueous phase and Recovery Agent).

9 Add  300 μ L 80% ethanol to the pellet while on the magnet. Wait  00:00:30 .

30s



10 Remove the ethanol.

- 11

Add

🧪

200 μ L

 80% ethanol to pellet. Wait

⌚

00:00:30

 .

30s
- 12

Remove the ethanol.
- 13

Centrifuge briefly. Place on the magnet. Low.

⚙️
- 14

Remove remaining ethanol. Air dry for

⌚

00:01:00

 .

1m
- 15

Remove from the magnet. Immediately add

🧪

35.5 μ L

 Elution Solution I.
- 16

Pipette mix (pipette set to 30 μ l) without introducing bubbles.
- 17

Incubate

⌚

00:02:00

 at

🌡️

Room temperature

 .

2m
- 18

Place on the magnet. Low until the solution clears.
- 19

Transfer

🧪

35 μ L

 sample to a new tube strip.

cDNA amplification

3d

- 20

Prepare cDNA Amplification Mix

🌡️

On ice

 . Add reagents in the order listed. Vortex and centrifuge briefly.

⚙️

cDNA Amplification Reaction Mix

A	B	C	D
Reagents	1X (μ l)	4X+10% (μ l)	8X+10% (μ l)
Amp Mix	50	220	440
cDNA Primers	15	66	132

A	B	C	D
Total	65	286	572

Calculations for cDNA Amplification Reaction Mix preparation.

21 Add  65 µL cDNA Amplification Reaction Mix to  35 µL sample.



22 Pipette mix 15x (pipette set to 90 µl). Centrifuge briefly.






23 Incubate in a thermal cycler with the following protocol:



A	B	C
Lid Temperature	Reaction Volume	Run Time
105°C	100 µl	~30-45 min
Step	Temperature	Time
1	98°C	3 min
2	98°C	15 sec
3	63°C	20 sec
4	72°C	1 min
5	Go to Step 2 – 11 cycles	
6	72°C	1 min
7	4°C	Hold

Thermocycler protocol.


24 Store at  4 °C for up to  72:00:00 or  -20 °C for ≤1 week or proceed to the next step.

3d



cDNA Cleanup – SPRIselect:

3d 0h 9m 30s

25 Vortex to resuspend the SPRIselect reagent. Add  60 µL SPRIselect reagent (0.6X) to each sample and pipette mix 15x (pipette set to 150 µl).







26 Incubate  00:05:00 at  Room temperature .

5m



27 Place on the magnet. High until the solution clears.

28 Remove the supernatant.

29 Add  200 μ L 80% ethanol to the pellet. Wait  00:00:30 .

30s




30 Remove the ethanol.

31 Repeat steps 29 and 30 for a total of 2 washes.

32 Centrifuge briefly and place on the magnet. Low.



33 Remove any remaining ethanol. Air dry for  00:02:00 .



2m

Note

DO NOT exceed 2 min as this will decrease elution efficiency.

34 Remove from magnet. Add  40.5 μ L Buffer EB. Pipette mix 15x.



35 Incubate for  00:02:00 at  Room temperature .




2m



36 Place the tube strip on the magnet. High until the solution clears.




37 Transfer  40 μL sample to a new tube strip.

38 Store at  4 °C for up to  72:00:00 or at  -20 °C for up to 4 weeks, or proceed to the next step.

3d



cDNA Quantification

39 Run  1 μL sample (Dilution Factor 1:10) on an Agilent Bioanalyzer High Sensitivity chip.

Protocol references

https://cdn.10xgenomics.com/image/upload/v1660261285/support-documents/CG000505_Chromium_Nuclei_Isolation_Kit_UG_RevA.pdf