

Aug 28, 2024

Preparation of qPCR standards for AAV titering

DOI

dx.doi.org/10.17504/protocols.io.bp2l61jndvqe/v1

Roberta Marongiu^{1,2}, Santiago Unda^{1,2}, Michael G. Kaplitt^{1,2}

¹Department of Neurological Surgery, Weill Cornell Medical College, New York, NY 10065;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD 20815, USA

ASAP Collaborative Rese...



Eileen Ruth Torres

Weill Cornell Medicine

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.bp2l61jndvqe/v1

Document Citation: Roberta Marongiu, Santiago Unda, Michael G. Kaplitt 2024. Preparation of qPCR standards for AAV titering. protocols.io <https://dx.doi.org/10.17504/protocols.io.bp2l61jndvqe/v1>

License: This is an open access document distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Created: August 01, 2022

Last Modified: August 28, 2024

Document Integer ID: 68024

Keywords: ASAPCRN, AAV, qPCR

Funders Acknowledgement:

Aligning Science Across

Parkinson's

Grant ID: 020608

Abstract

qPCR standards for AAV titting

document annex to protocol "AAV Purification Protocol with Iodixanol gradient"



1 ug of 1 kb DNA = 9.1×10^{11} molecules

Use any plasmid with known size with WPRE (so you can use WPRE primers)

Digestion 100ul reaction – do not use any enzyme with star activity and make sure to use a single restriction site enzyme.

# of reactions	volume
1	
DNA (10ug)	x
10x CutSmart buffer	10
XbaI	5
H2O	Up tp 9 5 ul
V/reaction	100

Mix digestion and incubate (37°C) for minimum 2 hours up to O/N.

PCR purification and elute in 30 ul EB.

Spec very well and perform multiple readings. Spec values (ug/ml)

Calculations:

If using our pAAV.CBA.mcherry.WPRE = 6132 bases=6.132 kb

1 ug => for 1000 bp dsDNA = 9.1×10^{11} molecules

1 ug => for pAAV.CBA.mcherry.WPRE = $9.1 \times 10^{11} / 6.132 = 1.484 \times 10^{11}$ molecules

Standard $10^8 = 2 \times 10^7$ molecules

2×10^7 molecules weight 0.135 ng = 2×10^7 molecules / 1.484×10^{11} molecules

2×10^7 molecules / ul = 0.135 ng / ul = 135 ng/ml

When we use 5 ul in qPCR => 0.675 ng in 5 ul @ 2×10^7 / ul

In 1 ml TE – 135 ng

Calculate the DNA volume to dilute in 1 ml of TE based on spec'ed concentration

