



AUG 08, 2023

## Passing of Organoids in Matrigel

Annika  
Fendler<sup>1</sup>

<sup>1</sup>Charité

Annika Fendler: Review



Annika Fendler

### ABSTRACT

This protocol is for passaging prostate cancer organoids in matrigel. Protocol modified from Drost et al. (see citation below)

OPEN  ACCESS



**Protocol Citation:** Annika Fendler 2023. Passing of Organoids in Matrigel. **protocols.io**  
<https://protocols.io/view/passaging-of-organoids-in-matrigel-cycexste>

**MANUSCRIPT CITATION:**  
Drost J, Karthaus WR, Gao D, Driehuis E, Sawyers CL, Chen Y, Clevers H. Organoid culture systems for prostate epithelial and cancer tissue. Nat Protoc. 2016 Feb;11(2):347-58. doi: 10.1038/nprot.2016.006. Epub 2016 Jan 21. PMID: 26797458; PMCID: PMC4793718.  
<https://www.nature.com/articles/nprot.2016.006>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



**Protocol status:** In development  
We are still developing and optimizing this protocol

**Created:** Aug 08, 2023

**Last Modified:** Aug 08, 2023

**PROTOCOL integer ID:**  
86118

## Before start


- 1 Prepare sufficient medium for passaging and prewarm immediately before use
- 2 Thaw  Growth Factor Reduced (GFR) Matrigel® Corning Catalog #354230 on ice o/n
- 3 Warm  TrypLE™ Express Enzyme Thermo Fisher Scientific Catalog #12604013 to 37°C before use
- 4 Coat 15 ml Falcons with 1% BSA o/n
- 5 Coat tips before use in 1% BSA by pipetting up and down a few times
- 6 Take pictures of the organoids to count organoids per well and measure the size of the organoids










### Note


Add results from counting here:


## Dissociate organoids

45m

- 7 Break up matrigel lense by pipetting up and down with a P1000 for approximately 10 times and transfer all content of a well to a coated 15 ml Falcon.  
Repeat with all remaining wells
- 8 Once the content of all wells is transferred to the Falcon, pipette the whole suspension up and down a few more times to further break up the matrigel.
- 9 Add up to  10 mL of PBS

- 10  500 x g, Room temperature, 00:05:00 5m
- The higher speed helps to pellet the cells below the matrigel layer. Otherwise they get trapped in the matrigel and it is easy to suck them up with the supernatant
- 11 Mark the top of the matrigel pellet. Remove as much supernatant as possible without disturbing the matrigel.
- 12 Add  5 mL PBS
- 13  500 x g, Room temperature, 00:05:00 5m
- 14 Remove supernatant as described above
- 15 Add  1-2 mL prewarmed  TrypLE™ Express Enzyme Thermo Fisher Scientific Catalog #12604013 and pipette up and down to make sure that the cell pellet is well resuspended.
- 16 Place falcon in incubator for  00:10:00 to  00:15:00 25m
- 17 Pipette up and down for 10 times every 5 min and check the progress under the microscope
- 18 When only small clusters and single cells remain, add  5 mL PBS
- 19  400 x g, Room temperature, 00:05:00 5m

20 Remove supernatant and resuspend pellet in another  5 mL PBS

21  400 x g, Room temperature,  
00:05:00

5m

22 Remove as much supernatant as possible and keep the Falcon tube on ice until ready.

## Seed cells in organoids



15m

23 Depending in the passaging ratio, decide how many wells you want to seed the cells into and calculate the necessary amount of matrigel.

Plate size	Volume matrigel per well (ul)	Number of wells seeded	Split ratio	Total matrigel volume (ul)	Volume medium per well (ul)
96-well					
48-well	20	6	1:0.8	125	200
24-well					
12-well					
6-well					

24 Add matrigel to the cell pellet and carefully resuspend cells. Avoid formation of bubbles.

25 Carefully add a matrigel drop to the center of the well.

26 Flip the plate and let the matrigel solidify at  37 °C in the incubator for  00:15:00

15m

27 Carefully add prewarmed medium on top of matrigel lense.

**28** Incubate at  37 °C in the incubator to expand cells.