



SEP 23, 2023

OPEN ACCESS



DOI:
dx.doi.org/10.17504/protocols.io.5jyl8pmndg2w/v1

Protocol Citation: Elias Adriaenssens 2023. Immunoprecipitation of NAP1-GFP. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.5jyl8pmndg2w/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Jun 28, 2023

Last Modified: Sep 23, 2023

PROTOCOL integer ID:
84128

Keywords:
Immunoprecipitation

Immunoprecipitation of NAP1-GFP

Elias Adriaenssens¹

¹Sascha Martens lab, University of Vienna, Max Perutz Labs - Vienna

ASAP Collaborative Research Network

Sascha Martens lab, University of Vienna, Max Perutz Labs - Vienna



Elias Adriaenssens

Sascha Martens lab, University of Vienna, Max Perutz Labs - ...

ABSTRACT

This protocol describes immunoprecipitation of NAP1-GFP from HAP1 cells.

ATTACHMENTS

[757-1925.pdf](#)

MATERIALS

Materials

- NAP1-EGFP (Addgene xxx)
- Lipofectamine 2000 (Thermo Fisher)
- SDS-PAGE gels (NP0322BOX, Thermo Fisher)
- PageRuler Prestained protein marker (Thermo Fisher)

Lysis buffer

A	B
KCl	100 mM
MgCl ₂	2.5 mM
Tris-HCl pH 7.4	20 mM
NP-40	0.50%

Fixation solution

A
40% ethanol
10% acetic acid
50% dH2O



Pierce™ Detergent Compatible Bradford Assay Kit Thermo
Fisher Catalog #23246



ChromoTek GFP-Trap®
Agarose Proteintech Catalog # gta-20

Immunoprecipitation

1h

1




Transiently transfect FIP200 knockout HAP1 cells with pcDNA3.1 NAP1-EGFP (Addgene), or empty pcDNA3.1 vector as a negative

control, using Lipofectamine 2000 (Thermo Fisher).

Note

This cell line was selected as FIP200 deletion results in TBK1 hyperactivation and thus increased NAP1 phosphorylation.




- 2 After 48 h, collect the cells by trypsinization and wash the cell pellet with PBS once. Then, lyse the cells in lysis buffer.

- 3 Lyse the samples for  00:20:00  On ice . Then, clear the cell lysates by centrifugation at  20000 x g, 4°C, 00:10:00 .




30m

- 4 Determine the protein concentrations of the cleared protein lysates with the Pierce Detergent Compatible Bradford Assay Kit (23246, Thermo Fisher).

- 5 For both samples, negative control and NAP1-EGFP lysates, incubate  12 mg of cell lysate  Overnight with  20 µL or GFP-Trap agarose beads (GTA-20, Chromotek).



20m

- 6 In the morning, wash the samples three times in lysis buffer. Then, resuspend the beads in protein loading dye, supplemented with 100 mM DTT. Finally, boil the samples for  00:05:00 at  95 °C .



5m

- 7 Load the samples on 4-12% SDS-PAGE gels (NP0322BOX, Thermo Fisher) with PageRuler Prestained protein marker (Thermo Fisher).

- 8 After the run, stain the SDS-PAGE gel with Coomassie, fix for  00:10:00 in fixation solution, and destain it  Overnight in dH2O.



15m

- 9 Cut the band corresponding to NAP1-EGFP from the gel with a fresh scalpel and submit for mass spectrometry analysis.