



•





Feb 11, 2022

# Protein Expression Using BL21(DE3) (C2527) V.2

## New England Biolabs<sup>1</sup>

<sup>1</sup>New England Biolabs

1



dx.doi.org/10.17504/protocols.io.bdd4i28w

New England Biolabs (NEB)

Tech. support phone: +1(800)632-7799 email: info@neb.com



This protocol explains methods for Protein Expression Using BL21(DE3) Competent E. coli (C2527).

You can find the protocol for a large-scale Protein Expression here.

DOI

dx.doi.org/10.17504/protocols.io.bdd4i28w

https://www.neb.com/protocols/0001/01/01/protein-expression-using-bl21de3-c2527

New England Biolabs 2022. Protein Expression Using BL21(DE3) (C2527). **protocols.io** 

https://dx.doi.org/10.17504/protocols.io.bdd4i28w

**New England Biolabs** 

T7 Expression Strain, expression plasmid, transform expression plasmid into BL21(DE3), BL21, C2527

\_\_\_\_\_ protocol,

Mar 08, 2020

Feb 11, 2022

33948



1

## BL21(DE3) Genotype:

fhuA2 [lon] ompT gal ( $\lambda$  DE3) [dcm]  $\Delta$ hsdS  $\lambda$  DE3 =  $\lambda$  sBamHlo  $\Delta$ EcoRl-B int::(lacl::PlacUV5::T7 gene1) i21  $\Delta$ nin5

#### **NOTES:**

- 1. **Caution:** This product contains DMSO, a hazardous material. Review the MSDS before handling.
- 2. **Storage and Handling**: Competent cells should be stored at -80°C. Storage at -20°C will result in a significant decrease in transformation efficiency. Cells lose efficiency whenever they are warmed above -80°C, even if they do not thaw.

### **MATERIALS**

BL21(DE3) Competent E.coli - 20x0.05 ml New England

**Biolabs Catalog #C2527H** 

■BL21(DE3) Competent E.coli - 6x0.2 ml New England

Biolabs Catalog #C25271

Please refer to the Safety Data Sheets (SDS) for health and environmental hazards.

**Caution:** This product contains DMSO, a hazardous material. Review the MSDS before handling.

- 1 Transform expression plasmid into BL21(DE3).
- Plate on antibiotic selection plates.



Incubate @ Overnight at § 37 °C.

4

Resuspend a single colony in 10 mL liquid culture with antibiotic.

5

Incubate at § 37 °C until OD<sub>600</sub> reaches 0.4–0.8.



Induce with  $\Box 4 \mu L$  or  $\Box 40 \mu L$  of [M]100 Milimolar (mM) IPTG stpck (final concentration of 40 or 400  $\mu$ M).

7

Induce for 3-5 hours at § 37 °C.

8 ~

Check for expression either by Coomassie stained protein gel, Western Blot or activity assay. Check expression in both the total cell extract (soluble + insoluble) and the soluble fraction only.

If a fraction of the target protein is insoluble, repeat expression at a lower temperature (15 to 30°C) or test expression in Lemo21(DE3) (NEB #C2528).

- 9 For large scale, inoculate □1 L of liquid medium (with antibiotic) with a freshly grown colony or □10 mL of freshly grown culture.
- 10

Incubate at & 37 °C until OD<sub>600</sub> reaches 0.4–0.8.

11

Add [M]40 Micromolar ( $\mu$ M) or [M]400 Micromolar ( $\mu$ M) IPTG and express protein using optimal time/temperature determined in a small scale trial.