



AUG 11, 2023

# Reverse transcription, primer pools preparation and multiplex PCR steps for CHIKV serotype genomic sequencing

Gabriel Luz

Laís Ceschini<sup>1</sup>, Luisa Maria Inácio da Silva<sup>1</sup>, Wallau<sup>1</sup>

<sup>1</sup>Entomology department, Instituto Aggeu Magalhães (IAM), FIOCRUZ



Luisa Maria Inácio da Silva

OPEN ACCESS



**DOI:**  
[dx.doi.org/10.17504/protocols.io.5qpvo3rm7v4o/v1](https://dx.doi.org/10.17504/protocols.io.5qpvo3rm7v4o/v1)

**Protocol Citation:** Laís Ceschini, Luisa Maria Inácio da Silva, Gabriel Luz Wallau 2023. Reverse transcription, primer pools preparation and multiplex PCR steps for CHIKV serotype genomic sequencing . **protocols.io** <https://dx.doi.org/10.17504/protocols.io.5qpvo3rm7v4o/v1>

**MANUSCRIPT CITATION:** Machado, Lais Ceschini, et al. "Genome Sequencing Reveals Coinfection by Multiple Chikungunya Virus Genotypes in a Recent Outbreak in Brazil." PLOS Neglected Tropical Diseases, edited by Fabrice Simon, vol. 13, no. 5, May 2019, p. e0007332, doi:<https://doi.org/10.1371/journal.pntd.0007332>.

## ABSTRACT

This step-by-step protocol describes the cDNA synthesis, primer pools preparation and multiplex PCR conditions with the main goal to sequence the complete genome of CHIKV serotype strains.

## MATERIALS

Reverse transcription: SuperScript™ IV First-Strand Synthesis System. (200 reactions) Cat:18091200 Invitrogen

Multiplex PCR: Q5® High-Fidelity 2X Master Mix. Cat: M0492L NEB, H2O Ultrapure, primers described in table 1.

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
We use this protocol and it's working

**Created:** Aug 10, 2023

**Last Modified:** Aug 11, 2023

**PROTOCOL integer ID:**  
86340

## Reverse transcription

- 1 Using a 2mL tube prepare the **Mix 1** described below for 96 samples:

A	B	C
Mix 1 Reverse transcription	Vol. (1x)	96 samples (+2 = 98 to keep some extra due to pipetting issues)
Random Hexamers (50µM)	1µL	98µL
dNTPs mix (10mM each)	1µL	98µL
Total	2 µL	194µL

- 2 Using 0,2mL PCR tubes or 96 wells plates add 11-16µL of extracted RNA from RT-PCR positive samples. Add **2µL** of **Mix 1** to the tube/well and take it to the thermocycler with the following set up:

65°C ---- 5 minutes

**3** Take the tubes/wells to ice for 1 minute. (you can prepare a water bath with ice cubes to have a uniform temperature distribution).

**4** Using a 2mL tube prepare **Mix 2**:

A	B	C
Mix 2 Reverse Transcription	Vol. (1x)	96 samples (+2 = 98 to keep some extra due to pipetting issues)
5x SSIV Buffer	4µL	392µL
100mM DTT	1µL	98µL
RNaseOUT or RNase Inhibitor	1µL	98µL
SSIV Reverse Transcriptase	1µL	98µL
Total	7µL	686µL

**5** Add **7µL of Mix 2** to the tubes containing the **Mix 1** plus RNA and take it to the thermocycler following the set up below:

Step1:

42°C ---- 50 minutes

70°C ---- 10 minutes

4°C ---- Hold

**6** Store the cDNA at -20°C.

Observation:.. As a suggestion, to improve the final results only samples RT-PCR positive showing a Ct value of < 30 should be used for cDNA conversion and genomic amplification.

## Pools of primers

**7** Select two 0,6mL tubes for each pool.

- 8 Using the original 100uM primer solution eluted individually, put them together following the table below containing each primer volume.
- 9 Pool 1 will have a final volume of 469µl and pool 2 of 460µl.
- 10 In order to prepare the solution to use in the Multiplex PCR, dilute each pool 1:10. That is, 10µl of pool 1 and 90µl of ultrapure water.

**TABLE 1: Primers and pool order**

A	B	C	D	E
Primer	Sequence	Concentration inside of the pool *	Volume of primer within the pool	Pool
400_1_LEFT_1	TGACACACG TAGCCTACC AGTT	0,015uM	10ul	1
400_1_RIGHT_1	CGCATCGGG CAAACGCAG TGGTA	0,015uM	10ul	1
400_3_LEFT_3	GCAGACGTC GCGATATAC CAAG	0,015uM	10ul	1
400_3_RIGHT_3	CCAGCTCTT AAGTAGCAT GCGG	0,015uM	10ul	1
400_5_LEFT_0	GATGTGCAA GACTACCGA CACG	0,015uM	10ul	1
400_5_RIGHT_0	GACTGGGTA TCAGGCCTC TTGT	0,015uM	10ul	1
400_7_LEFT_4	CAAGAAGCC CAGGATGCT GAAA	0,015uM	10ul	1
400_7_RIGHT_4	GCTATGCGT ACACGTCTT CACT	0,015uM	10ul	1
400_9_LEFT_4	GCAGAGAGG ACAGAACAC GAGT	0,015uM	10ul	1
400_9_RIGHT_4	CTCTCTGTCT CATCACGTC GGT	0,015uM	10ul	1

A	B	C	D	E
400_11_LEFT _0	AGCAGTGCG GCTTCTTCA ATAT	0,015uM	10ul	1
400_11_RIGHT T_0	TGCCTAACT GCGTAACT CCTTT	0,015uM	10ul	1
400_13_LEFT _2	ATTAAGGAG TGGGAGGTG GAGC	0,015uM	10ul	1
400_13_RIGHT T_2	TCTAGAATG GACGCTGCC TCAG	0,015uM	10ul	1
400_15_LEFT _4	GGGGAAAGA ATGGAATGG CTGG	0,015uM	10ul	1
400_15_RIGHT T_4	CGTTCACTG GTTCTATCTG CGT	0,015uM	10ul	1
400_17_LEFT _0	AACTGAACG CAGCCTTTG TAGG	0,015uM	10ul	1
400_17_RIGHT T_0	ACACCTGTG GAGAGGAGA GGTA	0,015uM	10ul	1
400_19_LEFT _1	CATACAGAT GCGGACCCA AGTG	0,015uM	10ul	1
400_19_RIGHT T_1	GTTCAGGAG TCATGGCAT AACGG	0,015uM	10ul	1
400_21_LEFT _2	GCGCGTAAG TCCAAGGGA ATAC	0,015uM	10ul	1
400_21_RIGHT T_2	GTCTCCGCT GTTTCTTGTA CGG	0,015uM	10ul	1
400_23_LEFT _1	ACTTTCGGA GACTTCCTA CCCG	0,015uM	10ul	1
400_23_RIGHT T_1	ACAGCCTCT CTTTAGTCTC TGGA	0,015uM	10ul	1
400_25_LEFT _2	ACCAAATCA CCGATGAGT ATGATGC	0,015uM	10ul	1
400_25_RIGHT T_2	TCGTTATCCT GATAGGGCT GGC	0,015uM	10ul	1
400_27_LEFT _4	AGGCCTAAG GTGCAGGTT ATACA	0,015uM	10ul	1
400_27_RIGHT T_4	GCAGGTGAC AGCTGGAAA TCTC	0,015uM	10ul	1

A	B	C	D	E
400_29_LEFT_0	CGATGAATT GATGGCAGC CAGA	0,015uM	10ul	1
400_29_RIGH T_0	GCAAAGGTG GCCATGGAC ATTA	0,015uM	10ul	1
400_31_LEFT_1	TTCTACAATA GGAGGTACC AGCCT	0,015uM	10ul	1
400_31_RIGH T_1	TTCATGCAC ATTCTCTCTC TGCG	0,015uM	10ul	1
400_33_LEFT_3	GATACCCGT GCACATGAA GTCC	0,015uM	10ul	1
400_33_RIGH T_3	TTTTTCGTAG CAGCAGGGT GTG	0,015uM	10ul	1
400_35_LEFT_0	CCACAAGAC CGTACCTAG CTCA	0,015uM	10ul	1
400_35_RIGH T_0	TGGTGAAAT GGGTGCGTA CATG	0,015uM	10ul	1
400_37_LEFT_3	AATGTCACA ACAGTCCGG CAAT	0,015uM	10ul	1
400_37_RIGH T_3	TTGGGTGGT CAGGATACA GCAA	0,015uM	10ul	1
400_39_LEFT_1	GGCCACCCG CATGAGATA ATTC	0,015uM	10ul	1
400_39_RIGH T_1	ATAGGACAA TCAGGGCTG CCAG	0,015uM	10ul	1
400_41_LEFT_0	CTTGGAACC AACGCTATC GCTT	0,015uM	10ul	1
400_41_RIGH T_0	AGCAGCCAC AGTGATATT ATTCCT	0,015uM	10ul	1
400_43_LEFT_0	ACCAGGACA ATTTGGCGA CATC	0,015uM	10ul	1
400_43_RIGH T_0	ATACCTCAC ACGACATGT CCGT	0,015uM	10ul	1
400_45_LEFT_3	CTACACAAG TACACTGTG CAGCC	0,015uM	10ul	1
400_45_RIGH T_3	TGTTATTCAG GGGTTGTTC AGCC	0,015uM	10ul	1

A	B	C	D	E
400_2_LEFT_0	CCAGCAAGG AGGATGATG TCGGAC	0,015uM	10ul	2
400_2_RIGHT_0	TGTGTCGAA CCCTACCCA GTAC	0,015uM	10ul	2
400_4_LEFT_0	TGTTCTCAGT AGGGTCAAC GCT	0,015uM	10ul	2
400_4_RIGHT_0	GGATGCCGG TCATTTGATC ACA	0,015uM	10ul	2
400_6_LEFT_1	TGAGAAGCT TTTGGGGGT CAGA	0,015uM	10ul	2
400_6_RIGHT_1	ACATCTTCCT GTGCTGCCT GTA	0,015uM	10ul	2
400_8_LEFT_0	ACTTTCCCC GCAGACCGT ATTA	0,015uM	10ul	2
400_8_RIGHT_0	CAGCTTCTT CCTTCTTGC AGCA	0,015uM	10ul	2
400_10_LEFT_0	ACCTGGTGA CTAGCGGAA AGAA	0,015uM	10ul	2
400_10_RIGHT_0	GACGACACA ATGGCAGTC ACAG	0,015uM	10ul	2
400_12_LEFT_0	GAGGGTGGG TTAAACAAC TGCA	0,015uM	10ul	2
400_12_RIGHT_0	TTATCCCCG CTGTTTCGA GGAT	0,015uM	10ul	2
400_14_LEFT_1	ACGCGGATA ACCACTGGG ATAA	0,015uM	10ul	2
400_14_RIGHT_1	TTATAGCCG CTAACCAGG AGCA	0,015uM	10ul	2
400_16_LEFT_0	AGGTGACTC ACTGAGACT GCTC	0,015uM	10ul	2
400_16_RIGHT_0	ATCGTTCTTC GCGATGTCC ATG	0,015uM	10ul	2
400_18_LEFT_2	GGACCAAAC TTCTCAAATT ACACGGA	0,015uM	10ul	2
400_18_RIGHT_2	CCAAACTAC TGTCAGGGT GCAC	0,015uM	10ul	2

A	B	C	D	E
400_20_LEFT_4	CAGAAATGC CCGGTGGAT GATG	0,015uM	10ul	2
400_20_RIGH T_4	ATCGGCGCT TAGATCAAA CTGAC	0,015uM	10ul	2
400_22_LEFT_0	GAGGGAGAA ACCTGACCG TGAT	0,015uM	10ul	2
400_22_RIGH T_0	AGTCATAAC TCGTCTGCC GTGT	0,015uM	10ul	2
400_24_LEFT_4	CACGGCCAA TAGAAGCAG GTATC	0,015uM	10ul	2
400_24_RIGH T_4	TTGACGGAT TGAATGTCG CTCG	0,015uM	10ul	2
400_26_LEFT_0	ACCCACTTT GGACTCAGC AGTA	0,015uM	10ul	2
400_26_RIGH T_0	AGGACGGCG TTCAATCTCC TAA	0,015uM	10ul	2
400_28_LEFT_0	CCAGGATGA TTCAC TTGC GCTT	0,015uM	10ul	2
400_28_RIGH T_0	GGAGCTTTC TGGGATACA ACTGC	0,015uM	10ul	2
400_30_LEFT_1	GATGGCAAC GAACAGGGC TAAT	0,015uM	10ul	2
400_30_RIGH T_1	GGTCTGGGT CTGATGACT TGGA	0,015uM	10ul	2
400_32_LEFT_3	CCCCAAAA AGAAACCGG TTCA	0,015uM	10ul	2
400_32_RIGH T_3	GAGTACTGT ACTGCTCCG TGGT	0,015uM	10ul	2
400_34_LEFT_0	CGTCACGAA AATCACCCC TGAG	0,015uM	10ul	2
400_34_RIGH T_0	TCTGTCGCTT CGTTTCTGAT GC	0,015uM	10ul	2
400_36_LEFT_0	CCGTGCACG ATTACTGGA ACAA	0,015uM	10ul	2
400_36_RIGH T_0	CACAATTGC ACTTGTACC GCAC	0,015uM	10ul	2



A	B	C	D	E
400_38_LEFT _1	TCCTCTGGC AAATGTGAC ATGC	0,015uM	10ul	2
400_38_RIGH T_1	CACCCACCA TCGACAGGA GTAT	0,015uM	10ul	2
400_40_LEFT _1	TATACCTGT GGAACGAGC AGCA	0,015uM	10ul	2
400_40_RIGH T_1	TGTACCGCA GCATTTAC GTAC	0,015uM	10ul	2
400_42_LEFT _4	TCAGCATAC AGGGCTCAT ACCG	0,015uM	10ul	2
400_42_RIGH T_4	GACGGTCTC TGCAGTACC AGTT	0,015uM	10ul	2
400_44_LEFT _0	ATCTCCATC GACATACCG GACG	0,015uM	10ul	2
400_44_RIGH T_0	TGTGACGCC GGGTAATTG ACTA	0,015uM	10ul	2
400_46_LEFT _1	TCCCTAAAG AGACACACC GCAT	0,015uM	10ul	2
400_46_RIGH T_1	TCTTAGCTAT ATATGGTGT GTCTCTTAG GG	0,015uM	10ul	2

\*approximate concentration of each primer in the 25µl PCR reaction.

Note: The primers were designed using the <https://primalscheme.com> based on the KP164576, KP164571, KP164572, KP164568, KP164570 and KP164569 reference genomes (Machado, 2019).

## Multiplex PCR

- 11** Prepare the **Mix 1** for a Multiplex PCR for each **Pool 1 and Pool 2** using a Falcon tube of 15mL (~96 amostras) or a 2mL tube.

A	B	C	D
Mix 1 Multiplex PCR	Vol. Pool 1 (1x)	Vol. Pool 2 (1x)	96 samples (+2) (pool1 ou pool2)

A	B	C	D
Q5 Master Mix High fidelity 2X	12,5 µl	12,5 µl	1.225 µl
Pool primers (Pool1 or Pool2) /Use concentration	1,7 µl	1,7 µl	166,6 µl
Ultra Pure Water	8,3 µl	8,3 µl	813,3 µl
Total	22,5µl	22,5µl	2205µl

- 12** Add **2,5µl of cDNA** (totalling 5µl) in 22,5µl of the pool1 and pool2 reaction and take it to the thermocycler following the conditions bellow:

Step1:

98°C ---- 30 seconds

Step2: (45 cycles)

98°C ---- 15 seconds

58°C ---- 30 seconds

72°C ---- 5 minutes

Step3:

72°C ---- 2 minutes

Hold 4°C