



Oct 20, 2020

Swarming assay (high number)

Forked from a private protocol

Serena Ding¹¹Imperial College London**1** *Works for me* dx.doi.org/10.17504/protocols.io.uagesbw**Behavioural Genomics**

Serena Ding

ABSTRACT

For imaging swarming behaviour of a very high number (thousands+) of young adult *C. elegans* on agar using the Cyclops single worm tracker. Worms are synchronised by bleaching and refeeding for 72 hours, and then 300 μ L of young adult hermaphrodites are transferred by glass pipette onto a 90 mm regular NGM plate for imaging for 15 hours at 0.1Hz.

DOI

dx.doi.org/10.17504/protocols.io.uagesbw

PROTOCOL CITATION

Serena Ding 2020. Swarming assay (high number). **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.uagesbw>

MANUSCRIPT CITATION please remember to cite the following publication along with this protocol

FORK FROM

Forked from a private protocol

KEYWORDS

behaviour, *C. elegans*, swarming, high number

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Oct 03, 2018

LAST MODIFIED

Oct 20, 2020

PROTOCOL INTEGER ID

16424

ABSTRACT

For imaging swarming behaviour of a very high number (thousands+) of young adult *C. elegans* on agar using the Cyclops single worm tracker. Worms are synchronised by bleaching and refeeding for 72 hours, and then 300 μ L of young adult hermaphrodites are transferred by glass pipette onto a 90 mm regular NGM plate for imaging for 15 hours at 0.1Hz.

Bleach synchronising worms (-7 to -4 day)

- 1 Bleach synchronise 3 large plates (90mm) worth of animals, making sure lots of gravid hermaphrodites are present. Leave on rotator at 20 °C until use.

Re-feeding worms (-3 day PM)

- 2 Re-feed worms onto 3-4 large plates (90mm) using a glass pipette. Incubate at 20 °C.



90 mm plates are pre-seeded with undiluted OP50. Lawn does not have to be any specific shape - uniform, with gaps, all is ok, depending on experimental needs.

Imaging (Day 0 PM)

- 3 Set appropriate frame rate using CommonVision control box (0.1 Hz). Open Gecko and initialise experiments: fixed duration (60 minute recording, repeat 15 times), MJPG format.
- 4 Pool animals together from 3 large plates (90mm), wash in M9 twice, and aspire as much supernatant as possible after the last wash.
- 5 Transfer ~160 uL of worm onto an imaging plate (90 mm plate pre-seeded with undiluted OP50 at least one day prior) using a glass pipette, making sure to carry as little liquid as possible as this will disperse bacteria on the imaging plate.



Some liquid will inevitably spill onto the lawn, which is ok. This is where the worm aggregate/front initially forms.

- 6 Place the plate inside the tracker with an appropriate field of view (i.e. worm front at the edge of the illuminated circular field of field, so that they are expected to sweep over the field of view). Image overnight.

Imaging (Day 1 AM)

- 7 Close Gecko. Transfer the files to the server.