



Aug 17, 2020

A simple method for purification of immunoglobulin-Y (IgY).

Angel A Justiz-Vaillant¹¹University of the West Indies St. Augustine**1** Works for me dx.doi.org/10.17504/protocols.io.bjuwknxeUniversity of the West Indies angel.vaillant@sta.uwi.eduAngel Justiz-Vaillant
University of the West Indies St. Augustine

DOI

dx.doi.org/10.17504/protocols.io.bjuwknxe

PROTOCOL CITATION

Angel A Justiz-Vaillant 2020. A simple method for purification of immunoglobulin-Y (IgY).. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bjuwknxe>

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 17, 2020

LAST MODIFIED

Aug 17, 2020

PROTOCOL INTEGER ID

40566

- 1 Separate egg yolk from the egg white.
- 2 Add ten parts of distilled water to one part egg yolk.
- 3 Mix and stir slowly for 6.30 h at 4°C.
- 4 Centrifuge at 10 000 × g, at 4°C for 27 min to precipitate the lipids.
- 5 Collect the supernatant containing the IgY.

- 6 Slowly add potassium sulfate to the preparation or specimen, stirring constantly, to a final concentration of 500 mM.
- 7 Adjust pH to 7.6.
- 8 Pass the sample through a 0.45 µm filter just before applying it to the HiTrap™ Column column.