



Dec 13, 2021

DNA Extraction from Bacteriophages

Frej Larsen¹¹Copenhagen University

1



protocol .

FOOD Micro UCPH



Frej Larsen

DNA/RNA extraction is performed using the **QIAamp® Viral RNA Mini kit from Qiagen** without the addition of carrier RNA to the AVL buffer and an additional nuclease step prior to extraction. The nuclease step is introduced to degrade any DNA/RNA that may be in the sample after the virome extraction protocol ([dx.doi.org/10.17504/protocols.io.b2qaqdse](https://doi.org/10.17504/protocols.io.b2qaqdse)) has been performed. Please note that even though the kit is a viral RNA extraction kit, viral DNA will also be extracted.

Frej Larsen 2021. DNA Extraction from Bacteriophages. **protocols.io**
<https://protocols.io/view/dna-extraction-from-bacteriophages-b2tmqek6>



protocol ,

Dec 13, 2021

Dec 13, 2021







55885





Before beginning this protocol, ensure that wash buffers have been properly diluted with ethanol and that the centrifuge is available and not chilled as low temperatures may impede ethanol evaporation.

Note that, unlike when using the kit directly, carrier RNA should not be added to the AVL buffer.

The Pierce Universal nuclease should be diluted 100 times (ie. by mixing 1 uL nuclease with 99 uL SM buffer/sterile water)

Wear gloves when performing this protocol.

- 1 For each sample, pipet  **1 µL** 100x diluted Pierce Universal Nuclease (2.5units/uL) to an empty 1.5 mL tube.
- 2 Transfer  **140 µL** of each sample from the outer chamber of the CentrisArt filter tube to the tube and mix by pipetting
- 3 Incubate for at least  **00:03:00** at room temperature. 3m
- 4 Add  **540 µL** AVL buffer to inactivate nucleases and lyse phage heads. Mix by pipetting or pulse vortexing.
- 5 Incubate  **00:10:00** at room temperature 10m
- 6 Change gloves.
- 7 Briefly centrifuge samples with a microcentrifuge.
- 8 Add  **560 µL** absolute ethanol. Mix thoroughly by pulse vortexing.

- 9 Briefly centrifuge samples with a microcentrifuge
- 10 Transfer  **630 µL** of the sample to a spin column
- 11 Centrifuge  **6000 x g, 21°C, 00:01:00** or until all liquid has passed through the filter. 1m
- 12 Change the collection tube and  **go to step #10** until all of the sample has passed through the filter.
- 13 Add  **500 µL** AW1 wash buffer to each spin column.
- 14 Centrifuge  **6000 x g, 21°C, 00:01:00** or until all liquid has passed through the filter 1m
- 15 Add  **500 µL** AW2 wash buffer to each spin column.
- 16 Centrifuge  **20000 x g, 21°C, 00:03:00** or until all liquid has passed through the filter. 3m
- 17 Change the collection tube, then centrifuge  **20000 x g, 21°C, 00:01:00** to dry the membrane. 1m
- 18 Place the spin column in a new RNase-free 1.5 mL tube.
- 19 Add  **30 µL** AVE elution buffer directly onto the filter membrane in the spin column.

20 Incubate at room temperature for at least 🕒 **00:01:00** . 1m

21 Centrifuge 🌀 **6000 x g, 21°C, 00:01:00** 1m

22 Discard the spin column and store eluate at 🧊 **-80 °C**