



Oct 12, 2022

CTAB DNA Extraction Protocol for Mollusks

Lavanya M Vythalingam¹, Tim Regan¹, Tim Bean¹, Carolina Peñaloza¹¹The Roslin Institute and Royal (Dick) School of Veterinary Studies, University of Edinburgh

2 Works for me

 Sharedx.doi.org/10.17504/protocols.io.e6nvwjpr2lmk/v1 Lavanya M Vythalingam
University of Edinburgh

ABSTRACT

This is a comprehensive CTAB extraction protocol intended for DNA extraction of mollusks. It includes detailed methodology and several reagent preparation.

DOI

dx.doi.org/10.17504/protocols.io.e6nvwjpr2lmk/v1

PROTOCOL CITATION

Lavanya M Vythalingam, Tim Regan, Tim Bean, Carolina Peñaloza 2022. CTAB DNA Extraction Protocol for Mollusks. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.e6nvwjpr2lmk/v1>



LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Sep 06, 2022

LAST MODIFIED

Oct 12, 2022

PROTOCOL INTEGER ID

69614

Method

1

Slice tissue sample using a clean, sterile scalpel. Put sample into a 1.5 ml tube with

CTAB extraction

 **400 µL** solution Teknova Catalog #C2190

2

Add **10 µL** **Proteinase K** Contributed by users **20 mg/mL** and invert tube to mix.

3

Incubate at **50 °C** for 60 minutes until lysed (or longer depending on tissue thickness).
Invert tube occasionally to mix sample.

4

Add  400 μ L

Chloroform/Isoamyl Alcohol (24:1) **Acros**

Organics Catalog #327155000

Organics Catalog #327155000 . Gently shake tube to emulsify sample.





5

Leave at room temperature for 2 minutes and shake tube to mix.

6

Spin tube at 13 000 rpm for 2 minutes and transfer upper layer to a 1.5 ml tube.

7

Add  **5 μ L** of  **RNase A** **Contributed by users**  **10 mg/mL** . Leave for 30 minutes at  **37 $^{\circ}$ C** .

8

Add 1 equal volume of

Chloroform/Isoamyl Alcohol (24:1) **Acros**

Organics Catalog #327155000

9

Spin tube at 13 000 rpm for 2 minutes and transfer upper layer to a new 1.5 ml tube containing

▢ **900 µL CTAB dilution solution**. Gently invert tube to mix. Leave in fridge overnight (optional).

10

Spin tube at 13 000 rpm for 10 minutes.

11

Discard supernatant with a micropipette and add ▢ **1000 µL 0.4 M NaCl in TE**. Invert tube to wash.

12

Spin tube at 13 000 rpm for 5 minutes and remove supernatant.

13

Add ▢ **300 µL 1.42 M NaCl in TE** to dissociate the CTAB from the DNA. Invert tube until pellet becomes transparent.

May require gentle shaking at ⚡ **37 °C**

14

Add ▢ **600 µL** ⊗ **Ethanol Contributed by users** (previously left at ⚡ **-20 °C**). Invert manually until precipitation appears.

Leave precipitating overnight for a higher DNA concentration

15

Centrifuge at full speed for 8 minutes, then remove supernatant. Dry at room temperature or in heat block at ⚡ **37 °C**


Notes

Reagent preparation

Reagent preparation

Reagent preparation

Reagent preparation

 Sodium Chloride **Contributed by**
users Catalog #PubChem CID: 5234

Tris HCl Buffer 1M Solution, Sterile pH 7.5 Bio Basic
Inc. Catalog #SD8124.SIZE.450ml

5 mL Technologies Catalog #AM9260G

Water to 100 mL



Citation: Lavanya M Vythaligam, Tim Regan, Tim Bean, Carolina Pe             CTAB DNA Extraction Protocol for Mollusks
<https://dx.doi.org/10.17504/protocols.io.e6nvwipr2lmk/v1>

This is an open access protocol distributed under the terms of the **Creative Commons Attribution License** (<https://creativecommons.org/licenses/by/4.0/>), which permits unrestricted use, distribution, and reproduction in any medium,

Tris HCl Buffer 1M Solution, Sterile pH 7.5 **Bio Basic**
Inc. Catalog #SD8124.SIZE.450ml

⊗ EDTA (0.5 M), pH 8.0 Life
Technologies Catalog #AM9260G

Water to  200 mL

0.4M Nacl in TE

Tris HCl Buffer 1M Solution, Sterile pH 7.5 **Bio Basic**
Inc. Catalog #SD8124.SIZE.450ml

⊗ EDTA (0.5 M), pH 8.0 **Life**

Technologies Catalog #AM9260G

1.42 M NaCl in TE

Tris HCl Buffer 1M Solution, Sterile pH 7.5 Bio Basic
Inc. Catalog #SD8124.SIZE.450ml

200 µL

EDTA (0.5 M), pH 8.0 Life

Technologies Catalog #AM9260G

Water to 100 mL