

May 17, 2024

## Western blot - alpha-synuclein

 Forked from [Western blot - alpha-synuclein](#)

DOI

**[dx.doi.org/10.17504/protocols.io.rm7vzjb22lx1/v1](https://dx.doi.org/10.17504/protocols.io.rm7vzjb22lx1/v1)**

Pietro La Vitola<sup>1</sup>, Eva M. Szegö<sup>1</sup>

<sup>1</sup>German Center for Neurodegenerative Diseases (DZNE)

ASAP Collaborative Rese...



**Pietro La Vitola**

German Center for Neurodegenerative Diseases (DZNE), Alignin...

OPEN  ACCESS



DOI: **[dx.doi.org/10.17504/protocols.io.rm7vzjb22lx1/v1](https://dx.doi.org/10.17504/protocols.io.rm7vzjb22lx1/v1)**

**Protocol Citation:** Pietro La Vitola, Eva M. Szegö 2024. Western blot - alpha-synuclein. **protocols.io**  
**<https://dx.doi.org/10.17504/protocols.io.rm7vzjb22lx1/v1>**

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** February 27, 2023

**Last Modified:** May 17, 2024

**Protocol Integer ID:** 99998

**Keywords:** alpha-synuclein, western blot, dosal medulla oblungata

**Funders Acknowledgement:****Aligning Science Across  
Parkinson's****Grant ID: ASAP-000420**

## Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK


The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to **protocols.io** is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with **protocols.io**, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

## Abstract

This protocol describes how to detect alpha-synuclein protein in mouse dorsal medulla oblongata tissue by western blot

## Materials

### Laemmli SDS sample buffer, not reducing

 Laemmli SDS-Sample buffer, not reducing **Thermo Scientific Catalog #J60660.AC**

### Towbin transfer buffer

25 mM Tris

192 mM glycine

 8.3

20% methanol (vol/vol).

1 L of buffer:


- 800 mL distilled H<sub>2</sub>O
- 200 mL methanol
- 3.03 g Tris base
- 14.4 g Glycine



## Sample preparation

1d 4h 16m

### 1 *Tissue homogenization*

1.1 Homogenise samples in 150  $\mu$ L ice-cold lysis buffer (1% Triton X-100 in 0.1 M phosphate buffered saline solution,  7.6), supplemented with protease and phosphatase inhibitors.

1.2 Centrifuge samples  14000 x g, 4°C, 00:30:00

30m

1.3 Transfer supernatants into pre-cooled test tubes.


1.4 Measure protein concentration (e.g. BCA method).

1.5 Mix 4  $\mu$ g samples in Laemmli sample buffer supplemented with 5 % beta-mercaptoethanol.

1.6 Heat samples  95 °C  00:05:00

5m

### 2 *Electrophoresis and transfer*

2.1 Separate samples (4  $\mu$ g total protein) by polyacrylamide gel electrophoresis using precast Bolt™ 4-20% Bis-Tris, 1.0 mm, Mini Protein Gels at 120 V  01:20:00 or until dye front reaches the bottom of the gel. Run with pre-stained size markers

1h 20m

2.2 Soak nitrocellulose membrane (pore size: 0.2  $\mu$ m) with Towbin buffer (see materials)

2.3 Soak transfer sandwich components (2 sheets of filter paper and 2 blotting pads) in Towbin transfer buffer and assemble in the transfer cassette in the following order:  
cathode plate  
1 x blotting pad  
1 x filter paper  
gel



nitrocellulose membrane  
1 x filter paper  
1 x pad  
Use a roller to remove any air bubbles

- 2.4 Place cassette in transfer tank and transfer protein onto 0.2  $\mu$ m nitrocellulose membranes (300 mA, ⌚ 01:30:00 in cold Towbin buffer.

1h 30m

- 2.5 Rinse membranes with TBS ⌚ 00:00:30

30s

### 3 *Immunodetection of protein bands*

- 3.1 Block non-specific binding sites with blocking solution (TBS-T: TBS containing 2% BSA and 0.05 % Tween-20) for ⌚ 01:00:00

1h

- 3.2 Incubate membranes in blocking solution containing mouse anti- $\alpha$ -synuclein (1:1000; RRID:AB\_398108) and rabbit anti- $\beta$ -actin (1:1000; RRID:AB\_2305186) ⌚ 18:00:00

18h

🔥 Room temperature

- 3.3 Wash with TBS-T ⌚ 00:05:00 4x

5m

- 3.4 Incubate membrane with peroxidase-conjugated goat anti-mouse and goat anti-rabbit IgG in TBS-T (1:5000 each). ⌚ 00:00:00 RT

- 3.5 Wash with TBS-T ⌚ 00:05:00 4x

5m

- 3.6 Incubate membrane with ECL ⌚ 00:00:30

30s

- 3.7 Visualize ECL stained membrane using a BioRad ChemiDoc™.



## Protocol references

### CITATION

Szegö EM, Van den Haute C, Höfs L, Baekelandt V, Van der Perren A, Falkenburger BH (2022). Rab7 reduces  $\alpha$ -synuclein toxicity in rats and primary neurons..

LINK

<https://doi.org/10.1016/j.expneurol.2021.113900>

## Citations

Szegö EM, Van den Haute C, Höfs L, Baekelandt V, Van der Perren A, Falkenburger BH. Rab7 reduces  $\alpha$ -synuclein toxicity in rats and primary neurons.

**<https://doi.org/10.1016/j.expneurol.2021.113900>**