



VERSION 2

APR 12, 2023

OPEN  ACCESS

DOI:

[dx.doi.org/10.17504/protocols.io.36wgqjkkxvk5/v2](https://dx.doi.org/10.17504/protocols.io.36wgqjkkxvk5/v2)

**Protocol Citation:** Sarah K Buddenberg 2023. Infection of *Biomphalaria glabrata* snails with *Schistosoma mansoni* miracidia. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.36wgqjkkxvk5/v2> Version created by Sarah K Buddenberg

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
We use this protocol and it's working

**Created:** Apr 11, 2023

**Last Modified:** Apr 12, 2023

**PROTOCOL integer ID:**  
80330

## Infection of *Biomphalaria glabrata* snails with *Schistosoma mansoni* miracidia V.2

Sarah K Buddenberg<sup>1</sup><sup>1</sup>Wellcome Sanger Institute

Schistosoma mansoni



Sarah K Buddenberg

Wellcome Sanger Institute, University of Cambridge

### ABSTRACT

To infect *Biomphalaria glabrata* snails with miracidia hatched from *Schistosoma mansoni* eggs

**Keywords:** miracidia,  
schistosoma, biomphalaria

MATERIALS

- ⊗ MilliQ water **Contributed by users**
- ⊗ 24-well Clear TC-treated Multiple Well Plates Individually Wrapped Sterile **Costar Catalog #3524**
- ⊗ Blunt featherweight forceps wide tip **BioQuip Catalog #4750**
- ⊗ 1x DPBS **Gibco - Thermo Fischer Catalog #14190144**
- ⊗ Glass Pasteur pipettes with rubber bulbs **Contributed by users**

1x Aquarium Water (diluted from 10X; see recipes)  
1-2L flask  
Aluminum foil

*B. glabrata* snails 5-8mm  
*S. mansoni* infected livers in PBS  
28°C room or incubator

Incubator at 37°C and 5% CO2  
Bucket of ice

Equipment	
Zoom stereomicroscope with diascopic illumination stand	NAME
stereomicroscope	TYPE
Nikon	BRAND
SMZ800N	SKU

10X AQUARIUM WATER  
5.56g CaCl2  
12.28g MgSO4-7H2O  
0.43g K2SO4  
4.2g NaHCO3  
480µl FeCl3-6H2O (0.5g/100ml water)  
Fill to 10L and store in 1L bottles  
Dilute to 1x: 9L MQ water + 1L 10x

Liver preparation

- 1 Collect livers from patent mice into 50ml Falcon tubes containing pre-warmed 37°C 1x DPBS

- 2 Remove livers from PBS and place in large mortar or laboratory blender
- 3 Gently homogenise liver tissue
- 4 Place homogenised liver slurry into a flask (size of flask depends on how many livers you have) and fill with diH<sub>2</sub>O or aquarium water

## Hatching miracidia

- 5 Cover flask in aluminum foil and shine light horizontally across the opening of the flask for 1-2 hrs
- 6 Take a small aliquot of water from the top of the flask using a glass Pasteur pipette and place into a petri dish. Using a stereomicroscope, check for swimming miracidia

## Exposing snails to miracidia for regular life cycle mainten.

- 7 Under a microscope, collect 15 miracidia with a glass Pasteur pipette per well  
  
Alternatively (but less preferred) estimate the number of miracidia by counting 12-5µl aliquots in a petri dish. After collecting (shaking the tube in between) and placing the aliquots in the petri dish add 5µl of Lugol to each drop (this kills and stains the miracidia)
- 8 Place individual 5-8mm snails into 24-well culture plates and cover snails completely with 1x aquarium water

- 9 Expose snails to miracidia for at least 3 hours (up to overnight)
- 10 After exposure, remove the snails carefully using featherweight forceps and put them in a new tank with food

## Exposing snails to miracidia for monomiracidium infections

- 11 Dilute miracidia so that one miracidium can easily be collected in ~3-5µl of water
- 12 Collected single miracidium using fresh 10µl pipette each collection and place in 24-well plate
- 13 After a plate is filled, check each well under a microscope to verify there is a single miracidium in each well
- 14 Place individual 5-8mm snails in the wells containing confirmed single miracidium and cover snails completely with 1x aquarium water
- 15 Leave the plate overnight (inside incubator or room at 28°C)
- 16 The following day, transfer the snails to a new tank with food

