

Sep 25, 2024

## FLAM-seq (with Kapa mRNA enrichment)\_CMS\_edit-2024-09-25

Forked from a private protocol

DOI

**[dx.doi.org/10.17504/protocols.io.j8nlk8q46l5r/v1](https://dx.doi.org/10.17504/protocols.io.j8nlk8q46l5r/v1)**

Corneliu Sologon<sup>1</sup>

<sup>1</sup>University of Miami



Corneliu Sologon

University of Miami

---

OPEN  ACCESS



DOI: **[dx.doi.org/10.17504/protocols.io.j8nlk8q46l5r/v1](https://dx.doi.org/10.17504/protocols.io.j8nlk8q46l5r/v1)**

**Protocol Citation:** Corneliu Sologon 2024. FLAM-seq (with Kapa mRNA enrichment)\_CMS\_edit-2024-09-25. **protocols.io**  
**<https://dx.doi.org/10.17504/protocols.io.j8nlk8q46l5r/v1>**

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** In development

**We are still developing and  
optimizing this protocol**

**Created:** September 25, 2024

**Last Modified:** September 25, 2024

**Protocol Integer ID:** 108373

### Abstract

fork\_CMS\_edit-2024-09-25



## Materials

### ***Kapa mRNA HyperPre Kit Illumina (cat. KK8440, Roche; Replaces TruSeq mRNA)***

- mRNA Capture Beads
- mRNA Bead Binding Buffer (BBB)
- mRNA Bead Wash Buffer(BWB)

### ***Poly(A) Tail Length Assay Kit (cat. 764551KT, Thermo Fisher)***

- 5X tail buffer mix
- 10X tail enzyme mix
- Universal Reverse (Univ. RV Primer; used in Sect 4: cDNA Amplification)

### ***SMARTer PCR cDNA cDNA Synthesis Kit (cat. 634926, Takara)***

- 5X First strand buffer
- DTT 20mM
- dNTP mix 10 mM
- RNase Inhibitor
- SMARTScribe RT
- 5' PCR Primer II A (used in Sect 4: cDNA Amplification)

### ***Advantage 2 DNA polymerase mix (cat 639207, Takara)***

- 10X Advantage 2SA PCR buffer
- dNTP Mix (10 mM each)

### ***RNAClean XP Beads (cat. A63987, Beckman Coulter)***

### ***XP DNA beads (cat. A63881, Beckman Coulter)***

### **Custom Oligonucleotides**

***isoTSO***- Template switch oligo (where “i” indicates stereoisomers of dC and dG as in the associated publication):

iCiGiCAAGCAGTGGTATCAACGCAGAGTACATrGrGrG

### ***RT primer 1:***

GGTAATACGACTCACTATAGCGAGANNNNNNNNNNNNCCCCCCCCTTT

### ***PCR primer 1 (to use in combination with RT primer 1; not used; replaced by 5' PCR Primer II A):***

GGTAATACGACTCACTATAGCGAG

### ***RT primer 2 (to use in alternative to 1):***

TGAGTCGGCAGAGAACTGGCGAANNNNNNNNNNNNCCCCCCCCTTT,

### ***PCR primer 2 (to use in combination with RT primer 2; not used; replaced by 5' PCR Primer II A):***

TGAGTCGGCAGAGAACTGGCGAA



## Poly(A)+ RNA preparation (using Kapa mRNA Beads)

10m 10s

### 1 Prepare mRNA Beads

1.1 52.5 µL mRNA Beads can be scaled for multiple samples

500 rpm, 00:01:00

1.2 Room temperature 00:05:00 magnet

Discard supernatant

5m

1.3 Remove from magnet

52.5 µL BBB

500 rpm, 00:01:00

1.4 Room temperature 00:05:00 magnet

Discard supernatant

1.5 Remove from magnet

52.5 µL BBB

500 rpm, 00:01:00

2 2-10 µg ds DNA in 0.2 mL PCR tube

50 µL QS; NAF

3 50 µL re-suspended mRNA Beads (from 1.5)

1000 rpm, 20°C, 00:01:00

4 65 °C 00:02:00

20 °C 00:05:00

undetermined, 00:00:10 , pulse

7m 10s

5 Room temperature 00:05:00 magnet

Discard supernatant

5m



6	Remove from magnet 🧪 200 µL BWB 🌀 1000 rpm, 20°C, 00:01:00	
7	🌡 Room temperature ⌚ 00:05:00 magnet Discard supernatant	5m
8	🧪 50 µL NAF 🌀 1000 rpm, 20°C, 00:01:00	
9	🌡 70 °C ⌚ 00:02:00 🌡 20 °C ⌚ 00:05:00 🔴 undetermined, 00:00:10 , pulse	7m 10s
10	🧪 50 µL BBB 🌀 1000 rpm, 20°C, 00:01:00 🌡 20 °C ⌚ 00:05:00	5m
11	🌡 Room temperature ⌚ 00:05:00 magnet Discard supernatant	5m
12	Remove from magnet 🧪 200 µL BWB 🌀 1000 rpm, 20°C, 00:01:00	
13	🌡 Room temperature ⌚ 00:05:00 magnet Discard supernatant	5m
14	Remove from magnet 🧪 16 µL NAF 🌀 1000 rpm, 20°C, 00:01:00	
15	🌡 70 °C ⌚ 00:02:00	2m
16	🌡 Room temperature ⌚ 00:05:00 magnet	5m



17  16  $\mu$ L  ds DNA to new 0.2 mL PCR tube

## GI Tailing (using USB poly(A) tail length assay, Thermo Fisher).

18 Prepare on ice (20 uL total):

 14  $\mu$ L RNA - poly(A)+

 4  $\mu$ L 5X tail buffer mix

 2  $\mu$ L 10X tail enzyme mix

19  37 °C  01:00:00

1h



20  1.5  $\mu$ L Stop Solution



2m

 4 °C  00:02:00

21  38.7  $\mu$ L XP RNA Beads (1.8X)



5m

 Room temperature  00:05:00

22  Room temperature  00:03:00 magnet


3m


Discard supernatant

23  50  $\mu$ L 80% EtOH  00:00:30

30s

Discard supernatant


24  go to step #23 x1

25  undetermined, 00:00:10 , pulse


10s

Discard supernatant

26 Remove from magnet

 17  $\mu$ L NAF

 1000 rpm, 20°C, 00:01:00

27  16  $\mu$ L  ds DNA to new 0.2 mL PCR tube



## cDNA synthesis (using SMARTscribe reverse transcriptase kit, Clontech).

28 Prepare RT master mix at **room temperature** (22 uL total; 1.1X per sample):

🧴 8 µL 5X First STrand buffer

🧴 1.5 µL DTT (20 mM)

🧴 4 µL dNTP mix (10 mM)

🧴 2 µL RNase Inhibitor

🧴 2 µL isoTSO (12 uM)

🧴 2 µL SMARTScribe RT

🧴 2.5 µL NAF

29 In 0.2 mL PCR tube from **step 27**

🧴 16 µL 🧬 ds DNA

🧴 2 µL RT Primer 1 (10 uM); custom 🧬 ds DNA

30 🔄 1000 rpm, Room temperature , 00:01:00

10s

⚙️ undetermined, 00:00:10 , pulse

31 🌡️ 72 °C ⌚ 00:03:00

1h 13m

🌡️ 42 °C **HOLD** - add 🧴 22 µL RT master mix

🌡️ 42 °C ⌚ 01:00:00

🌡️ 70 °C ⌚ 00:10:00

🌡️ 4 °C **HOLD**

32 🧴 24 µL XP DNA Beads (0.6X)

5m

🌡️ Room temperature ⌚ 00:05:00

33 🌡️ Room temperature ⌚ 00:03:00 magnet






Discard supernatant

34 🧴 50 µL 80% EtOH ⌚ 00:00:30

Discard supernatant


























35 ➡️ go to step #23 x1



- 36  undetermined, 00:00:10 , pulse  
Discard supernatant
- 37 Remove from magnet  
 42  $\mu\text{L}$  NAF  
 1000 rpm, 20°C, 00:01:00
- 38  40  $\mu\text{L}$   ds DNA to new 0.2 mL PCR tube

## cDNA library amplification (using Advantage 2 PCR enzyme system, Clontech).

11m 25s

- 39 Prepare the cDNA library amplification rxn:  
 40  $\mu\text{L}$   ds DNA  
 10  $\mu\text{L}$  10X Advantage 2SA PCR buffer  
 2  $\mu\text{L}$  dNTP (10 mM)  
 2  $\mu\text{L}$  5' PCR Primer II A (SMARTer PCR kit)  
 2  $\mu\text{L}$  Univ. RV Primer (Poly(A) Tail Kit)  
 42  $\mu\text{L}$  NAF  
 2  $\mu\text{L}$  Advantage 2 Polymerase Mix
- 40 Perform PCR  
1.  98 °C **HOLD** - add cDNA rxn  
2.  98 °C  00:01:00  
3.  98 °C  00:00:10  
4.  63 °C  00:00:15  
5.  68 °C  00:03:00 **repeat 40.3 -40.5 24X**  
6.  68 °C  00:07:00  
7.  4 °C **HOLD**
- 41  60  $\mu\text{L}$  XP DNA Beads (0.6X)  
 Room temperature  00:05:00
- 42  Room temperature  00:03:00 magnet  
Discard supernatant

11m 25s



43  200  $\mu$ L 80% EtOH  00:00:30


Discard supernatant

44  go to step #23 x1

45  undetermined, 00:00:10 , pulse

Discard supernatant

46 Remove from magnet

 42  $\mu$ L NAF

 1000 rpm, 20°C, 00:01:00

47  40  $\mu$ L  ds DNA to new 0.2 mL PCR tube

## Sequencing library preparation (using the SMRTbell™ Template Prep Kit)

48 Performed in PacBio sequencing core