

Aug 28, 2024

Fluorescent Immunolabelling for endogenous mouse LRRK2

DOI

dx.doi.org/10.17504/protocols.io.261genqkog47/v1

Roberta Marongiu^{1,2}, andrew.west west^{3,2}

¹Department of Neurological Surgery, Weill Cornell Medical College, New York, NY 10065;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD 20815, USA;

³Duke University

ASAP Collaborative Rese...



Eileen Ruth Torres

Weill Cornell Medicine

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.261genqkog47/v1

Protocol Citation: Roberta Marongiu, andrew.west west 2024. Fluorescent Immunolabelling for endogenous mouse LRRK2. protocols.io <https://dx.doi.org/10.17504/protocols.io.261genqkog47/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 19, 2022

Last Modified: August 28, 2024

Protocol Integer ID: 68930

Keywords: ASAPCRN, LRRK2



Funders Acknowledgement:

**Aligning Science Across
Parkinson's**

Grant ID: 020608

Abstract

This protocol is used to label endogenous LRRK2 in mouse.

It has been edited from West et al., 2014

Protocol optimized for fresh fixed brain tissue sectioned at 40 um.

Materials

⊗ Primary Antibody anti-LRRK2 c41-2 **Abcam Catalog #ab133474**

⊗ DAPI **Thermo Fisher Scientific Catalog #D1306**

⊗ Donkey serum **Sigma Aldrich Catalog #S30-100ML**

Solutions:














Tris Buffer Solution TBS - 2L

- 50mM Tris pH 7.5 – 100 mL of 1M stock

- 150 mM NaCl – 17.53 g

- Distilled Water – 1.9 L



- 1 Select brain slices, from fresh, unfrozen, fixed tissue - slices are 40um thick
- 2 Wash 3 times for  00:10:00 in TBS 10m
- 3 Sections were “quenched” in 100% Methanol for  00:15:00 at  4 °C on shaker. 15m
- 4 Wash 3 times for  00:10:00 in TBS 10m
- 5 Incubate the slices with Sodium Citrate – 10 mM Sodium Citrate pH 6.0 w/0.05% Tween for  00:30:00 at  37 °C on shaker. 30m
- 6 Wash 3 times for  00:10:00 in TBS 10m
- 7 During the washes, prepare the blocking solution : 5% normal donkey serum (matched to secondary AB) in TBS with 0.3% Triton
- 8 Incubate the slides with blocking solution for  01:00:00 at  4 °C on shaker. 1h
- 9 Wash 3 times for  00:10:00 in TBS 10m
- 10 During the washes, prepare the primary antibody solution : 5% normal donkey serum in TBS with NO TRITON + Primary Antibody anti-LRRK2 (1:500 or 1:1000)
- 11 Incubate at  4 °C on shaker for  24:00:00  48:00:00 3d
- 12 Wash 3 times for  00:10:00 in TBS 10m
- 13 During the washes, prepare the secondary antibody solution : 5% normal donkey serum in TBS with NO TRITON + Alexa Fluor secondary AB (1:1000)



-
- 14 Incubate at 4 °C on shaker for 18:00:00 18h
- 15 Wash 2 times for 00:05:00 in TBS 5m
- 16 During the washes, prepare the DAPI solution : DAPI (1mg/ml): 1:10,000, in TBS
- 17 Incubate the slides with the DAPI solution for 00:10:00 10m
- 18 Wash 3 times for 00:10:00 in TBS 10m
- 19 Mount sections on superfrost + slides and coverslip with Invitrogen Prolong Antifade reagent.