



Sep 18, 2020

# BSCI:414 Lab 3 Analyze Agarose Gel Results and Create RNase P Primers

Harley King<sup>1</sup>

<sup>1</sup>NIST Center for Neutron Research National Institute of Standards and Technology Gaithersburg, MD 20899 Department of Materials Science and Engineering University of Maryland College Park, MD 20742-2115

In Development

This protocol is published without a DOI.

USG Fall 2020 BSCI:414



Harley King

NIST Center for Neutron Research National Institute of Stand...

## ABSTRACT

BSCI:414 Lab 3

## PROTOCOL CITATION

Harley King 2020. BSCI:414 Lab 3 Analyze Agarose Gel Results and Create RNase P Primers . **protocols.io**  
<https://protocols.io/view/bsci-414-lab-3-analyze-agarose-gel-results-and-cre-bmh8k39w>

## LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Sep 18, 2020

## LAST MODIFIED

Sep 18, 2020

## PROTOCOL INTEGER ID

42272

### Review Lab 2 Experiment; Analyze Data

- 1 Review Lab 2 hypothesis and experimental setup.
- 2 Watch experiment process video, pausing to discuss.  
<https://youtu.be/SJoLg-GCeOY>
- 3
  1. Analyze PCR fragment lengths on agarose gel. You can find a copy of the gel in our ELMS course: Files>Lab Results>Lab 2 Agarose Gel>"Lab2. agarose.gel.50.samples.jpg"
  2. Which lanes have bands? Which lanes do not? Why?
  3. Discuss cases in which F and R primers point in same direction or away from each other.
- 4 Copy results to lab notebook. Is your band size in agreement with your hypothesis?

### Find CDC Primers in RNase P Gene

- 5 In the "BSCI:414 Plasmids" under the class root folder, find gene "RNaseP\_gene\_NM\_006413". Copy this gene to a new folder you create under "BSCI:414 Lab 3" using your name e.g. "F20\_HarleyKing\_Lab 3". Also give the plasmid a new

name like "F20\_HarleyKing\_RNAseP\_gene\_NM\_006413".

- 6 Navigate to the "CDC Coronavirus Real-Time RT-PCR Primers and Probes" site:  
<https://www.cdc.gov/coronavirus/2019-ncov/lab/rt-pcr-panel-primer-probes.html>

- 7 Make primers in Benchling using the primer and probe sequences from "RNAse P". Make sure to save the primers underneath your Lab 3 folder. The primers should have a memorable name like, "F20\_HarleyKing\_Lab3\_RNAseP-FWD".

|      |                        |  |                   |       |
|------|------------------------|--|-------------------|-------|
| RP-F | RNAse P Forward Primer | AGA TTT GGA CCT GCG AGC G                        | None              | 500nM |
| RP-R | RNAse P Reverse Primer | GAG CGG CTG TCT CCA CAA GT                       | None              | 500nM |
| RP-P | RNAse P Probe          | FAM – TTC TGA CCT GAA GGC TCT GCG CG – BHQ-1     | FAM, BHQ-1        | 125nM |
| RP-P | RNAse P Probe          | FAM-TTC TGA CCT /ZEN/ GAA GGC TCT GCG CG-3IABkFQ | FAM, ZEN, 3IABkFQ | 125nM |

Find these sequences in Benchling and make primers. Use the "Manual" method to make primers.

- 8 Answer the following questions in your lab notebook with a table:
1. What is the T<sub>m</sub> of the forward primer? How many bp is it?
  2. What is the T<sub>m</sub> of the reverse primer? How many bp is it?
  3. What is the T<sub>m</sub> of the probe quenched with BHQ1? How many bp is it?
  4. What is the T<sub>m</sub> of the probe quenched with 3IABkFQ? How many bp is it?
- 9 Copy and paste the emission spectra of the FAM probe in your lab notebook. What is its excitation/emission wavelength?