

Sep 26, 2020

# EC isolation from human mesenteric artery for scRNA-seq

Xiaofang Tang<sup>1</sup>, Kiran KS Sriram<sup>1</sup>, Yingjun Luo<sup>1</sup>, Zhen Chen<sup>1</sup><sup>1</sup>City of Hope

1 Works for me This protocol is published without a DOI.

Xiaofang Tang

## ABSTRACT

This protocol outlines a method to efficiently isolate the human mesenteric arterial endothelial cells. Following the physical scraping, endothelial cells from the inner surface of the vessels are separately liberated using a digestion buffer solution. The procedure takes approximately one hour and has been validated for down-stream applications such as FACS, VE-cad immunostaining and single cell RNA profiling.

## PROTOCOL CITATION

Xiaofang Tang, Kiran KS Sriram, Yingjun Luo, Zhen Chen 2020. EC isolation from human mesenteric artery for scRNA-seq. **protocols.io**  
<https://protocols.io/view/ec-isolation-from-human-mesenteric-artery-for-scrn-bf37jqrn>

## KEYWORDS

endothelial cells, isolation, EC isolation, human, mesenteric artery, human mesenteric artery, scRNA-seq

## LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

May 06, 2020

## LAST MODIFIED

Sep 26, 2020

## PROTOCOL INTEGER ID

36703

## GUIDELINES

Cells should be washed and resuspended in LoBind 2.0 mL eppendorf tubes and wide-bore pipette tips for better cell viability.

Pipette the cell suspension very gently and slowly for all steps involved.

Each 4-5 cm vessel, use **1 mL TrypLE**.

Keep the tissues/cells **On ice**.

## MATERIALS

NAME	CATALOG #	VENDOR
DPBS (10X), calcium, magnesium	14080055	Thermo Fisher
Medium 199	M2520-10X	Sigma Aldrich
40 µm cell strainer	14100150	Fisher Scientific
Trypan Blue Solution	MT25900CI	Corning
Disposable Safety Scalpels	6008TR-10	Mycro Instrumentation
Fisher BioReagents™ Bovine Serum Albumin Heat Shock Treated	BP1600-100	Fisher Scientific

NAME	CATALOG #	VENDOR
Pipette Tips RT LTS 1000µL LW 768A/8	30389219	Rainin
DNA LoBind Tube 2.0 mL	022431048	Eppendorf
Pipet-Lite Pipette Unv. SL-1000XLS	17014407	Rainin
30 G X 1 in Needle	305128	BD Biosciences
Pipette Tips RT LTS 1000µL FLW 768A/8	30389218	Mettler Toledo
TrypLE™ Express Enzyme (1X) phenol red	12605010	Thermo Fisher Scientific

#### STEPS MATERIALS

NAME	CATALOG #	VENDOR
Pipette Tips RT LTS 1000µL FLW 768A/8	30389218	Rainin

#### MATERIALS TEXT

Additional materials required:

- 10-cm petri-dish
- 5 mL tubes
- Sterile surgical tools
- Candle platform
- 0,04% BSA-DPBS

#### SAFETY WARNINGS

For hazard information and safety warnings, please refer to the SDS (Safety Data Sheet).

#### BEFORE STARTING

Warm up  **1 mL TrypLE** at  **37 °C** for  **00:05:00** before use.

Prepare at least  **6 mL 0.04% BSA-DPBS** .

#### EC isolation from human mesenteric artery for scRNA-seq

1 


Wash freshly isolated human mesenteric artery (~5 cm in length) with DPBS in a 10-cm petri-dish. Use **sterile surgical tools** to clean the vessel (Removing the fat and outer connective tissues).



Freshly isolated human mesenteric artery.

2 Cut open the vessel and expose the inner lumen. Anchor the four corners of the vessel using 23G needles on a black

wax.



Anchored vessel.

3 

Add **1 mL TrypLE** to the lumen side of the artery.

4 Use a scalpel to scrape the lumen 10-15 times (the purpose is to dissociate the endothelium of the intima).

5 

Keep the remaining vessel in DPBS **On ice**. The rest of the vessel can be scraped again to collect the vascular smooth muscle cells in the intima layer.

6 

Transfer **1 mL cell suspension** to a 5 mL tube.

7 

Incubate with shaking at ~ **150 rpm, 37°C** for **00:05:00**.

8  

Add equal volume ( **1 mL** ) of M199 medium (or any EC growth media) to stop enzyme digestion.

9 

Centrifuge at **600 x g, 4°C** for **00:05:00**.

10 

Remove the supernatant and re-suspend the pellet gently and thoroughly with **1 mL 0.04% BSA-DPBS** using wide-orifice LiteTouch pipette tips.



Pipette Tips RT LTS 1000µL FLW  
768A/8  
by Rainin  
Catalog #: 30389218

11



Pass the cell suspension through a 40 µm strainer to remove any cell debris. Wash the strainer by adding another **4 mL 0.04% BSA-DPBS**.

12



Centrifuge at **600 x g, 4°C** for **00:05:00**.

13



Remove the supernatant and re-suspend the pellet with **0.5 mL 0.04% BSA-DPBS** by using wide-orifice LiteTouch pipette tips and gently pipetting up and down along the inner wall of the tubes.



Pipette Tips RT LTS 1000µL FLW  
768A/8  
by Rainin  
Catalog #: 30389218

14



Count the cell viability:

Mix **10 µl trypan blue** + **10 µl cell suspension**. Take **10 µl mixture** to count cell number.



Ensure there's no cell debris, significant number of doublets or cell clumps in the suspension.

**Always keep the cells On ice.**



Usually we could get at ~ **10<sup>5</sup> cells in total**.

15



The EC isolation results.



2.png

Morphology of isolated mesenteric arterial endothelial cells after two passages (Passages 2 and 3).

flow-1.png

CD31 flow cytometry of isolated mesenteric arterial ECs, HUVEC was used as a positive control

- 16 Use appropriate number of cells for scRNA library preparation following the *10X Genomics 3' Expression protocol*. The rest of the cells can be cultured or used for other downstream assays.