



APR 15, 2024

## Fiber Photometry Protocol

In 2 collections

Sasha Burwell<sup>1</sup>

<sup>1</sup>Duke University

ASAP Collaborative Research Network



Sasha Burwell  
Duke University

OPEN ACCESS



### ABSTRACT

This protocol details collection of fiber photometry from VTA dopamine neurons.

### GUIDELINES

Recordings collected with a Tucker-Davis Technologies RZ10X and TDT Synapse software.

#### DOI:

[dx.doi.org/10.17504/protocols.io.eq2lyw39qvx9/v1](https://dx.doi.org/10.17504/protocols.io.eq2lyw39qvx9/v1)

**Protocol Citation:** Sasha Burwell 2024. Fiber Photometry Protocol. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.eq2lyw39qvx9/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



**Protocol status:** Working

**Created:** Apr 08, 2024

**Last Modified:** Apr 15, 2024

## On every day of behavior (see Behavior protocol):

5h

- 1 Photobleach the mono fiberoptic patchcords (Doric Lenses Inc., MFP\_400/430/1100-0.57\_1mm\_FCM-MF1.25\_LAF) (Synapse -> Preview Mode -> Fiber Bleaching) for at least  01:00:00 , and up to  04:00:00 , in both the 415nm and 465nm channels.
  - This can be done the night before to save time during experimental hours
- 2 Use a BNC cable to connect the solenoid output circuitry (see Behavior Hardware Setup protocol) and a ground connection to one of the RZ10x BNC ports.
  - Use a User Input Gizmo to record input from the correct BNC port
  - This will mean that, whenever the solenoid gets a +5V signal to open from the behavior MATLAB code, the RZ10X will also record an analog +1 signal. This enables later alignment of the fiber photometry recordings with the behavior data for analysis.
- 3 Use the PM1 power meter (Preview Mode -> Power Meter) to check and adjust, as necessary, the LED parameters:
  - The 415nm LED level should be at an amperage that produces a 20 $\mu$ W reading
  - The 465nm LED level should be at an amperage that produces a 25 $\mu$ W reading
- 4
 

**Note**

The system is ready to record. Prior to starting a behavior session with any given mouse:

Headfix the mouse.
- 5 Gently clean the tips of the implanted mono fiberoptic cannula (Doric Lenses Inc., MFC\_400/430-0.66\_5mm\_MF1.25\_FLT) and the mono fiberoptic patchcord with a cotton tipped applicator soaked in 70% ethanol.

- 6 Dry with an air can, then connect the cannula and the patchcord using a zirconia sleeve (Doric Lenses Inc., SLEEVE\_ZR\_1-BK).
  - Make sure the cannula and patchcord are firmly touching – if there is a gap in between, the recording performance will be low.
- 7 Insert the correct Mouse ID into Synapse.
- 8 Switch the Synapse software into Record Mode, and turn on both the 415 and 465 LED drivers.
  - You will be able to clearly see light at the tip of the patchcord, even through the zirconia sleeve, when the LED drivers are on and correctly working.
- 9 Run the behavior session (see Behavior protocol).
  - Start the fiber photometry recording before starting behavior, and end it after behavior is done, so that the entire behavior session time is collected as fiber photometry time points
- 10 Once the behavior session is complete, end the recording by returning the software to Idle Mode.
- 11 Remove the patchcord and zirconia sleeve connector from the mouse's cannula, and return the mouse to its home cage.
- 12 Repeat steps 4-11 with every mouse undergoing behavior.